

# Content

## 内容

<b>1. ECG Monitoring Cables &amp; Lead wires</b>	<b>05</b>
<b>心电监护电缆与导联线</b>	
1. 1 ECG Cables & Lead wires Cross reference 心电电缆导联线对照表	07
1. 2 Major types or brands ECG lead wires 主要类型或品牌心电导联线	16
<b>2. Telemetry ECG Cables</b>	<b>19</b>
<b>心电遥测电缆</b>	
<b>3. Disposable ECG Monitoring lead wires</b>	<b>21</b>
<b>一次性监护心电导联线</b>	
3. 1 Disposable ECG lead wires 一次性监护分体心电导联线	22
3. 2 Disposable One piece ECG lead wires 一次性一体化心电导联线	25
<b>4. ECG Lead Adapter</b>	<b>26</b>
<b>心电导联转换头</b>	
<b>5. Diagnosis ECG Cables &amp; Accessories</b>	<b>30</b>
<b>诊断用心电电缆及附件</b>	
5. 1 One piece EKG fixed cable 一体化心电电缆	33
5. 2 EKG Trunk Cable 心电主电缆	35
5. 3 EKG Lead wire 心电导联线	37
5. 4 EKG Accessories 附件	39

# ECG Cables & Accessories

## 心电电缆及附件

### Certification

#### 产品认证



CE Mark  
CE 标志



FDA 510(k) clearance No: K082959  
(Patient Cable & Leadwire 病人电缆&导联线)



SFDA医疗器械注册证:  
证书号码 No. 1210080

### Compliance with

#### 符合标准

ANSI/AAMI EC53-1998(R)2001  
EN 60601-1:1990  
ISO10993-1:2009E  
ISO10993-5:2009  
ISO10993-10:2009



# ① ECG Monitoring Cables & Lead wires

## 心电监护电缆与导联线

### Specification 产品规格

Cable material 线材	TPU Jacket, Light gray TPU外被, 浅灰色
Connector & Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Main Cable length 主电缆长度	2.5m
Lead wire length 导联线长度	0.9m

### Compatibility 兼容性

Major brands ECG or multi-parameter monitoring system  
主要品牌的心电和多参数监护系统

Shielded Cable & lead wire 屏蔽型电缆和导联线	Photo 图片	REF# 型号	REF # Rule Description 型号命名规则说明
One piece Fixed Cable 一体化电缆 (Shielded 屏蔽)		2586P	2: Shielded( 1=unshielded) 5: 5-lead(3= 3 lead) 86: factory/model code N: with neonate clip end(Default is adult) S: snap end(P=grabber end) I: IEC code (Default is AHA) 2: 屏蔽 (1=非屏蔽) 5: 5-导联(3=3-导联) 86: 原厂/规格型号 N: 新生儿夹 (缺省: 成人) 见下面图片 S: 扣式末端(P=钳式) I: 欧标 (缺省: 美标)
Trunk Cable 主电缆 (Shielded 屏蔽)		AA-2585	AA: cable yoke type, up to lead wire 2: Shielded( 1=unshielded) 5: 5-lead(3= 3 lead) 86: factory/model code Default : AHA code AA: 分线盒类型, 取决于导联线 2: 屏蔽 (1=非屏蔽) 5: 5-导联(3=3-导联) 86: 原厂/规格型号 缺省: 美标
Lead Wire 导联线 (Shielded 屏蔽)		NKB6-90S-I	NK: lead wire type( Acc to yoke) 6: 6-lead (3= 3 lead,5=5lead) 90: 90cm(wire length) S: snap end(P=grabber end) I: IEC code(Default : AHA code) NK: 导联线类型 (按分线盒) 6: 6-导联(3=3-导联,5=5-导联) 90: 90cm(导联线长度) S: 扣式末端 (P=钳式) I: 欧标 (缺省: 美标)











Unshielded Cable & lead wire 非屏蔽型电缆和导联线	Photo 图片	REF# 型号	REF # Rule Description 型号命名规则说明
Din trunk Cable D型主电缆 (1.5mm, Unshielded 非屏蔽)  Used with Din lead wire especially for neonate ECG monitoring 与D型导联线配合使用, 常用在 新生儿ECG监护系统		D-1540-I	D: Din style cable 1: Unshielded 5: 5-lead(3= 3 lead) 40: factory/model code I: IEC code (Others see below) D: D型电缆 1: 非屏蔽 5: 5-导联(3=3-导联) 40: 兼容原厂代码 I: 欧标 (其它如下)
Din lead wire D型导联线 (1.5mm, Unshielded 非屏蔽)  Used with Din trunk cable only 仅与D型主电缆配合使用		D5-90S-I	D: Din lead wire 5: 5-lead (3= 3 lead) 90: 90cm(wire length) S: snap end(P=grabber end) I: IEC code (Others see below) D: D型导联线 5: 5-导联(3=3-导联) 90: 90cm(导联线长度) S: 扣式末端 (P=钳式) I: 欧标 (其它如下)
Neonate lead wire 新生儿导联线 (1.5mm, Unshielded 非屏蔽)		D3-90NP-I	D: Din lead wire 3: 3-lead (5= 5 lead) 90: 90cm(wire length) NP=mini grabber end only I: IEC code (Others see below) D: D型导联线 3: 3-导联(5=5-导联) 90: 90cm(导联线长度) NP=迷你钳式末端 I: 欧标 (其它如下)











ECG Color Code on electrodes (Regular grabber/Mini grabber/Regular snap) ECG电极颜色规则 (常规钳式/迷你钳式/常规扣式)						
Code 代码	AHA(A) 美标	Alternate AHA(H) 替代美标	Standard IEC(I) 欧标	Alternate IEC(IH) 替代欧标	AHA(A) 美标	Standard IEC(I) 欧标
Grabber 钳式						
Mini Grabber 新生儿钳式						
Snap 扣式						













# 1.1 ECG Cables & Lead wires Cross reference 心电电缆导联线对照表


OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
<b>AAMI</b>									
Without resistance (无电阻) BCI Datascope MDE/Invivo Mindray Spacelabs WelchAllyn Protocol		3	D-1340		D3-90S	D3-90P	2340S	2340P	
		Price	D-1340-I		D3-90S-I	D3-90P-I	2340S-I	2340P-I	
		5	D-1540		D5-90S	D5-90P	2540S	2540P	
		Price	D-1540-I		D5-90S-I	D5-90P-I	2540S-I	2540P-I	
With 1K resistance (1K电阻) Same as above		3	D-1340R		D3-90S	D3-90P	2340RS	2340RP	
		Price	D-1340R-I		D3-90S-I	D3-90P-I	2340RS-I	2340RP-I	
		5	D-1540R		D5-90S	D5-90P	2540RS	2540RP	
		Price	D-1540R-I		D5-90S-I	D5-90P-I	2540RS-I	2540RP-I	
LL style		3	LL-2340		LL3-90S	LL3-90P	2340S	2340P	D-1340
		Price	LL-2340-I		LL3-90S-I	LL3-90P-I	2340S-I	2340P-I	D-1340-I
		5	LL-2540		LL5-90S	LL5-90P	2540S	2540P	D-1540
		Price	LL-2540-I		LL5-90S-I	LL5-90P-I	2540S-I	2540P-I	D-1540-I





OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
<b>Artema/S&amp;W</b>									
8000,Diascope1,-2, Diascope NT, Athena		3	NA	NA	NA	NA	2397S	2397P	D-1397
		Price					2397S-IH	2397P-IH	D-1397-IH
<b>BCI See AAMI</b>									
<b>Biolight宝莱特(China)</b>									
M7000,M9500 AS Modular		5	NA	NA	NA	NA	2557S	2557P	D-1557
		Price					2557S-I	2557P-I	D-2557-I
AnyView A2/A3/A5/A6/A8 Q3/Q5/Q7/V6 Vital Sign Monitor		3	NA	NA	NA	NA	2367S	2367P	D-1367
		Price					2367S-I	2367P-I	D-1367-I
<b>Bionet(Korea)</b>									
		5	NA	NA	NA	NA	2514S	2514P	D-1514
		Price					2514S-I	2514P-I	D-2514-I
<b>Biosys(Korea)</b>									
		3	NA	NA	NA	NA	2335S	2335P	D-1335
		Price					2335S-I	2335P-I	D-1335-I
<b>CAS MEDICAL</b>									
		3	NA	NA	NA	NA	2339S	2339P	D-1339
		Price					2339S-IH	2339P-IH	D-1339-IH
<b>Colin/Omron(Japan)</b>									
BP 308/308E/408	See AAMI								
BP88S		3	CL-2520		CL3-90S	CL3-90P	2320S	2320P	D-1320
		Price	CL-2520-I		CL3-90S-I	CL3-90P-I	2320S-I	2320P-I	D-1320-I
		5	CL-2520		B7	B7	C4	C4	
		Price	CL-2520-I		CL5-90S	CL5-90P	2520S	2520P	D-1520
BP88 BP306		3	CL-2521		CL3-90S	CL3-90P	2321S	2321P	D-1321
		Price	CL-2521-I		CL3-90S-I	CL3-90P-I	2321S-I	2321P-I	D-1321-I
		5	CL-2521		B7	B7	C4	C4	
		Price	CL-2521-I		CL5-90S	CL5-90P	2521S	2521P	D-1521
<b>Criticare See AAMI</b>									
602A-A, Poet 2, 504US,Scholar,506, 507,508,1100,1508, 8100									
<b>Creative科瑞康(China)</b>									
		3	NA	NA	NA	NA	2330S	2330P	D-1330
		Price					2330S-I	2330P-I	D-1330-I
		5	NA	NA	NA	NA	2530S	2530P	D-1530
		Price				2530S-I	2530P-I	D-1530-I	

OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
Datascope /Mindray See AAMI									
Datascope -Mindray V12, V21		3	DT-2526 DT-2526-I		DT3-90S	DT3-90P	2326S	2326P	D-1326
		Price	C6		DT3-90S-I	DT3-90P-I	2326S-I	2326P-I	D-1326-I
		5	DT-2526 DT-2526-I		DT5-90S	DT5-90P	2526S	2526P	D-1526
		Price	C6		DT5-90S-I	DT5-90P-I	2526S-I	2526P-I	D-1526-I
Dixtal									
Ep12			NA	NA	NA	NA	21070S	21070P	NA
		Price					21070S-I	21070P-I	
Drager									
UM1/3/3.1, RM1, PM8010, PM8014, Babyguard8000, Cato, PB8800, Series 8000, Cicero, Parameterbox		3	NA	NA	NA	NA	2377S	2377P	D-1377
		Price					2377S-I	2377P-I	D-1377-I
							C4	C4	C4
ECAFIX									
		5	NA		NA	NA	2552S 2552S-I	2552P 2552P-I	NA
Fukuda Danshi									
		3	SM-2599 SM-2599-I		SM5-90S	SM5-90S	2599S	2599P	D-1599
		Price	C6		SM5-90S-I	SM5-90S-I	2599S-I	2599P-I	D-1599-I
		12	MQ-21099		MQ5-90S	MQ5-90P	21099S	21099P	NA
		Price	D3		MQ-21099-I	MQ5-90S-I	MQ5-90P-I	21099S-I	
-DS5100E, DS5300W, DS5400, DS7001 Mindray/Datascope /Expert		5	NA	NA	NA	NA	2559S 2559S-I	2559P 2559P-I	D-1559 D-1559-I
		Price					C6	C6	C6
GE/Critikon									
Dianmap Pro 1000 MPS Select Portable monitors		3	DX-2340R DX-2340R-I		DX3-90S	DX3-90P	2340RS	2340RP	D-1340R
		Price	C6		DX3-90S-I	DX3-90P-I	2340RS-I	2340RP-I	D-1340R-I
					B7	B7	C1	C1	C1

OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
GE-Datex-Ohmeda									
AS/3, CS/3, Cardiocap II, Satlite Plus S5 Without Resister for respiration		3	DX-2392 DX-2392-I		DX3-90S	DX3-90P	2392S	2392P	D-1392
		Price	C6		DX3-90S-I	DX3-90P-I	2392S-I	2392P-I	D-1392-I
		5	DX-2592 DX-2592-I		DX5-90S	DX5-90P	2592S	2592P	D-1592
		Price	C6		DX5-90S-I	DX5-90P-I	2592S-I	2592P-I	D-1592-I
AS/3, CS/3, Cardiocap II, Satlite Plus S5 With resisters		3	DX-2395 DX-2395-I		DX3-90S	DX3-90P	2395S	2395P	D-1395
		Price	C6		DX3-90S-I	DX3-90P-I	2395S-I	2395P-I	D-1395-I
		5	DX-2595 DX-2595-I		DX5-90S	DX5-90P	2595S	2595P	D-1595
		Price	C6		DX5-90S-I	DX5-90P-I	2595S-I	2595P-I	D-1595-I
		12	MQ-21095		MQ5-90S	MQ5-90P	21095S	21095P	NA
		Price			MQ5-90S-I	MQ5-90P-I	21095S-I	21095P-I	
					MQ-21095-I	MQ5-90S-I	MQ5-90P-I	21095S-I	21095P-I
					NA	2XC2	2XC2	NA	NA
GE/Hellige									
Eagle 1000 Def.. SCP 900/ Cardioserv SMC, SMS, SMK, SMU, SMV		3	DX-2355 DX-2355-I		DX3-90S	DX3-90P	2355S	2355P	D-1355
		Price	C6		DX3-90S-I	DX3-90P-I	2355S-I	2355P-I	D-1355-I
		5	DX-2555 DX-2555-I		DX5-90S	DX5-90P	2555S	2555P	D-1555
Tram 7200		5	MQ-2584 MQ-2584-I		MQ5-90S	MQ5-90P	2584S	2584P	D-1584
		Price	C6		MQ5-90S-I	MQ5-90P-I	2584S-I	2584P-I	D-1584-I
					C2	C2	C6	C6	C6
GE/Marquette									
Eagle, Solar, Dash monitor MAC-Lab Cath Lab Datex-Ohmeda M-Prestrn & "E" modules		3	MQ-2586 MQ-2586-I		MQ3-90S	MQ3-90P	2386S	2386P	D-1386
		Price	C6		MQ3-90S-I	MQ3-90P-I	2386S-I	2386P-I	D-1386-I
		5	MQ-2586 MQ-2586-I		MQ5-90S	MQ5-90P	2586S	2586P	D-1586
		Price	C6		MQ5-90S-I	MQ5-90P-I	2586S-I	2586P-I	D-1586-I
		6	MQ-2686 MQ-2686-I		MQ6-90S	MQ6-90P	2686S	2686P	NA
		Price	C8.2		MQ6-90S-I	MQ6-90P-I	2686S-I	2686P-I	
		12	MQ-21086 MQ-21086-I		MQ5-90S	MQ5-90P	21086S	21086P	NA
		Price	D3		MQ5-90S-I	MQ5-90P-I	21086S-I	21086P-I	
					2XC2	2XC2	D3	D3	
		3	MQA-2382 MQA-2382-I		MQA3-90S	MQA3-90P	2382S	2382P	NA
		Price	C6		MQA3-90S-I	MQA3-90P-I	2382S-I	2382P-I	
		3	MQB-2381 MQB-2381-I		MQB3-90S	MQB3-90P	2381S	2381P	NA
		Price	C6		MQB3-90S-I	MQB3-90P-I	2381S-I	2381P-I	

OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
<b>Nihon Kohden</b>									
Life Scope: OPV-1500		3	(JC-103TA) (JC-103T)		NKB3-90S	NKB3-90P	2316S	2316P	D-1316
		Price	NA		NKB3-90S-I	NKB3-90P-I	2316S-1	2316P-I	D-1316-I
		3	NK-2540 NK-2540-I		NKA3-90S	NKA3-90P	2340S	2340P	D-1340
		Price	C6		NKA3-90S-I	NKA3-90P-I	2340S-I	2340P-I	D-1340-I
		5	NK-2540 NK-2540-I		NKA5-90S	NKA5-90P	2540S	2540P	D-1540
		Price	C6		NKA5-90S-I	NKA5-90P-I	2540S-I	2540P-I	D-1540-I
<b>Philips/HP</b>									
43100,43200A, 77020A/AC/AR Serie 78300, serie 78500 Serie 78600, M2400A,M2406A, M2424A		3	AA-23029 AA-23029-I		AA3-90S	AA3-90P	23029S	23029P	D-13029
		Price	C6		AA3-90S-I	AA3-90P-I	23029S-I	23029P-I	D-13029-I
		5	AA-25029 AA-25029-I		AA5-90S	AA5-90P	25029S	25029P	D-15029
Merlin 78352C/78354C/ 78834CM1001A/ M1002A, Virdia, Omicare CMA, ML65A,66A Heartstream Codemaster XL ...		3	AA-2385 AA-2385-I		AA3-90S	AA3-90P	2385S	2385P	D-1385
		Price	C6		AA3-90S-I	AA3-90P-I	2385S-I	2385P-I	D-1385-I
		5	AA-2585 AA-2585-I		AA5-90S	AA5-90P	2585S	2585P	D-1585
Modules: M1001B, M1002B, Multi-Measurement Modules: M3000A, M3001A, M3002A, MP2: M8102A, MP5: M8105A, SureSigns Monitors: 863064, 863065, 863066, 863068, 863087, 863088		3	HP-2385 HP-2385-I		HP3-90S	HP3-90P	2385S	2385P	D-1385
		Price	C7		HP3-90S-I	HP3-90P-I	2385S-I	2385P-I	D-1385-I
		5	HP-2585 HP-2585-I		HPA5-90S	HPA5-90P	2585S	2585P	D-1585
Multi-Measurement Modules: M3000A, M3001A Opt C12 or Opt C18, M3002A Opt C12 or Opt C18; Patient Monitors: MP5: M8105A Opt C12, MP2: M8102A Opt C12, ECG Modules: M1001B, M1002B, 860335, 860338		10	HP-21085		HPA5-90S (blue)	HPA5-90P (Blue)	21085S	21085P	
		Price	D3		HPB5-90S (Red)	HPB5-90P (Red)	21085S-I	21085P-I	
					C4X2	C4X2	D3	D3	

OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
<b>Physio Control</b>									
Lifepak5,6,7,8,9,9P, VSM2,1440,Vingmed 5	See GE Critikon	3							
		Price							
LIFEPAK12, LIFEPAK20E, LIFEPAK15		3						2327S	2327P
		Price						2327S-I	2327P-I
		5						2527S	2527P
		Price						2527S-I	2527P-I
With resistors							C6	C6	
	NA	12	MD-21027		MDA4-120S MDB6-90S			21027S	21027P-I
		Price	D3		MDA4-120S-I MDB6-90S-I			21027S-I	21027P-I
		Price	D3		D3			D3	D3
<b>Schiller</b>									
<b>Siemens</b>									
Sirecust series 400 Sirecust 341,610/ 620/630/720/730/ 732/900/960/961 1260/1281/1801, 400-500 ECG		3	SM-2376 SM-2376-I		SM3-90S	SM3-90P	2376S	2376P	D-1376
		Price	C6		SM3-90S-I	SM3-90P-I	2376S-I	2376P-I	D-1376-I
		5	SM-2576 SM-2576-I		SM5-90S	SM5-90P	2576S	2576P	D-1576
SC5000,SC6000, SC7000 SC8000,SC9000 Multi-Parameter Cable SM-STE5		3	SM-STE5		SMB3-90S	SMB3-90P			
		Price	D8		SMB3-90S-I	SMB3-90P-I			
		5	SM-STE5		B7	B7			
	See the following page	5	SM-STE5		SM5-90S	SM5-90P			
		Price	D8		SM5-90S-I	SM5-90P-I			
		3	SM-2323 SM-2323-I		SM3-90S	SM3-90P	2323S	2323P	D-1323
Spacelabs	See AAMI	3	SL-2396 SL-2396-I		SL3-90S	SL3-90P	2396S	2396P	D-1396
		Price	C6		SL3-90S-I	SL3-90P-I	2396S-I	2396P-I	D-1396-I
		5	SL-2596 SL-2596-I		SL5-90S	SL5-90P	2596S	2596P	D-1596
Ultraview 90492, 91220, 90496, 91496		5	SL-2596 SL-2596-I		SL5-90S-I	SL5-90P-I	2596S-I	2596P-I	D-1596-I
		Price	C6		C2	C2	C6	C6	C6
		12	MQ-21096 MQ-21096-I		MQ5-90S MV5-90S	MQ5-90P MV5-90P	21096S	21096P	
	NA	12	MQ-21096 MQ-21096-I		MQ5-90S-I MV5-90S-I	MQ5-90P-I MV5-90P-I	21096S-I	21096P-I	
		Price	D3		2XC2	2XC2	D3	D3	






OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
Welch allyn									
Protocol	See AAMI								
Micropaq Propaq LT		3	D-1353 D-1353-I		D3-90S D3-90S-I	D3-90P D3-90P-I	2353S 2353S-I	2353P 2353P-I	
		Price	C7	A8	B3	C1	C1		
		5	D-1553 D-1553-I		D5-90S D5-90S-I	D5-90P D5-90P-I	2553S 2553S-I	2553P 2553P-I	
		Price	C7	B4	B6	C2	C2		
Zoll									
R SERIES E SERIES M SERIES M CCT SERIES		5	NA		NA	NA	21069S 21069S-I	21069P 21069P-I	
		Price					D3	D3	











One piece cables with neonate electrode end for most brands are available upon request.  
For neonate application, Din style neonate leads are also the alternatives.  
有各品牌新生儿一体化电缆及D型电缆+新生儿导联线，适用新生儿心电监护。


















## 1.2 Major types or brands ECG lead wires




### 主要类型或品牌心电导联线

Types of Major ECG lead wires 主要类型导联线	Lead No. 导联数	REF# Snap(AHA/IEC) 扣式 (美标/欧标)	OEM REF# OEM原厂型号	Price code 代格代码	REF# Grabber(AHA/IEC) 钳式 (美标/欧标)	OEM REF# OEM原厂型号	Price code 代格代码
Din unshielded Individual set leads D型分离式导联线(非屏蔽)							
	3	D3-90S D3-90S-I		A8	D3-90P D3-90P-I		B3
	5	D5-90S D5-90S-I		B4	D5-90P D5-90P-I		B6
	7	D7-90S/D7-90S-I		B6	D7-90P/D7-90P-I		B7
	10	D10-90S/D10-90S-I		B8	D10-90P/D10-90P-I		C2
Din unshielded individual set (neonate) D型分离式新生儿导联线 (非屏蔽)							
	3				D3-90NP D3-90NP-I (Mini Grabber)		B4
	5				D5-90NP D5-90NP-I (Mini Grabber)		B6
LL style (shielded) LL型分离式导联线(屏蔽)							
	3	LL3-90S LL3-90S-I		B7	LL3-90P LL3-90P-I		B7
	5	LL5-90S LL5-90S-I		C2	LL5-90P LL5-90P-I		C2
Mindray/Datascope lead 心电导联线							
	3	MR3-90S MR3-90S-I		C1	MR3-90P MR3-90P-I		C1
	5	MR5-90S MR5-90S-I	012-00-1503-02	C2	MR5-90P MR5-90P-I		C2
	3	DT3-90S DT3-90S-I	0012-00-1503-04/5/6 45cm,61cm,100cm	C2	DT3-90P DT3-90P-I	0012-00-1514-04/5/6 45cm,61cm,100cm	C2
	5	DT5-90S DT5-90S-I	0012-00-1503-01/2/3 45cm,61cm,100cm	C2	DT5-90P DT5-90P-I	0012-00-1514-01/2/3 45cm,61cm,100cm	C2
Datex-GE Integrated leads (DX300 series) 集成式导联线 (DX300系列)							
	3	DX3-90S DX3-90S-I	545327	B7	DX3-90P DX3-90P-I	545317 545315	B7
	5	DX5-90S DX5-90S-I	545328	C2	DX5-90P DX5-90P-I	545318 545316	C2

OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
<b>Innocare</b>									
IM200:36281		5	NA		NA	NA	2550S 2550S-I	2550P 2550P-I	NA
		Price					C2	C2	
<b>Kontron</b>									
Micromon, Minimon Vitalmon, Supermon Kolormon		3	NA		NA	NA	2398S 2398S-IH	2398P 2398P-IH	D-1398 D-1398-IH
		Price					C4	C4	C4
105, 107, 108, 125, 128, VSM510. Old generation		3	NA		NA	NA	2393S 2393S-I	2393P 2393P-IH	D-1393 D-1393-IH
		Price					C4	C4	C4
<b>Lohmerier</b>									
M010,M011/M211 Defi 501		3	LM-2358 LM-2358-I		LM3-90S LM3-90S-I	LM3-90P LM3-90P-I	2358S 2358S-I	2358P 2358P-I	D-1358 D-1358-I
		Price	C6		B7	B7	C4	C4	C4
		5	LM-2558 LM-2558-I		LM5-90S LM5-90S-I	LM5-90P LM5-90P-I	NA	NA	NA
		Price	NA		C2	C2			
<b>MEK</b>									
MP500, MP600, MP1000		3	NA		NA	NA	2310S 2310S-I	2310P 2310P-I	D-1310 D-1310-I
		Price					C2	C2	C2
		5	NA		NA	NA	2510S 2510S-I	2510P 2510P-I	D-1510 D-1510-I
		Price					C4	C4	C4
		3	NA		NA	NA	2311S 2311S-I	2311P 2311P-I	D-1311 D-1311-I
		Price					C1	C1	C1
		5	NA		NA	NA	2511S 2511S-I	2511P 2511P-I	D-1511 D-1511-I
		Price					C2	C2	C2
<b>Mennen</b>									
Series 9000, 1100,2000		3	NA		NA	NA	2374S 2374S-I	2374P 2374P-I	D-1374 D-1374-I
		Price					C4	C4	C4
		5	NA		NA	NA	2574S 2574S-I	2574P 2574P-I	D-1574 D-1574-I
		Price					C6	C6	C6
Envoy		3	AA-2375 AA-2375-I		MN3-90S MN3-90S-I	MN3-90P MN3-90P-I	2375S 2375S-I	2375P 2375P-I	D-1375 D-1375-I
		Price	C6		B7	B7	C4	C4	C4
		5	AA-2575 AA-2575-I		MN5-90S MN5-90S-I	MN5-90P MN5-90P-I	2575S 2575S-I	2575P 2575P-I	D-1575 D-1575-I
		Price	C6		C2	C2	C6	C6	C6

OEM & Models OEM原厂型号	Monitor Plug 仪器端接头	Lead No. 导联数	Trunk Cable 主电缆	Lead yoke Conn. 导联线接头	Lead wire/w Snap 导联线: 扣式	Lead wire/w Grabber 导联线: 钳式	One piece/w Snap 一体化扣式	One piece/w Grabber 一体化钳式	Din Trunk Cable D型主电缆
<b>Mindray</b>									
PM7000/8000/9000	See AAMI								
PM5000 PM6000 BeneView T5/T6/T8		3	MR-2512 MR-2512-I		MR3-90S MR3-90S-I	MR3-90P MR3-90P-I	2312S 2312S-I	2312P 2312P-I	D-1312 D-1312-I
		Price	C6		B7	B7	C4	C4	C4
		5	MR-2512 MR-2512-I		MR5-90S MR5-90S-I	MR5-90P MR5-90P-I	2512S 2512S-I	2512P 2512P-I	D-1512 D-1512-I
		Price	C6		C2	C2	C6	C6	C6
<b>NEC</b>									
YCE205		3	NA		NECA3-90S NECA3-90S-I	NECA3-90P NECA3-90P-I	NA	NA	NA
		Price			B7	B7			
47502		3	NA		NECB3-90S NECB3-90S-I	NECB3-90P NECB3-90P-I			
		Price			B7	B7			
47504		3	NA		NECC3-90S NECC3-90S-I	NECC3-90P NECC3-90P-I	NA	NA	NA
		Price			B7	B7			
<b>Nihon Kohden</b>									
Def. Cardiolife Lifescopes series 6-8-9-11 6102-6205-6120		3	NK-2508 NK-2508-I		NKA3-90S NKA3-90S-I	NKA3-90P NKA3-90P-I	2308S 2308S-I	2308P 2308P-I	D-1308 D-1308-I
		Price	C6		B7	B7	C4	C4	C4
		5	NK-2508 NK-2508-I		NKA5-90S NKA5-90S-I	NKA5-90P NKA5-90P-I	2508S 2508S-I	2508P 2508P-I	D-1508 D-1508-I
		Price	C6		C2	C2	C6	C6	C6
		3	NK-2507 NK-2507-I		NKA3-90S NKA3-90S-I	NKA3-90P NKA3-90P-I	2307S 2307S-I	2307P 2307P-I	D-1307 D-1307-I
		Price	C6		B7	B7	C4	C4	C4
		5	NK-2507 NK-2507-I		NKA5-90S NKA5-90S-I	NKA5-90P NKA5-90P-I	2507S 2507S-I	2507P 2507P-I	D-1507 D-1507-I
		Price	C6		C2	C2	C6	C6	C6
TEC-7621C/TEC-7631 TEC-7721C/TEC-7731 7700		3	NA		NA	NA	2313S 2313S-I	2313P 2313P-I	
		Price					C6	C6	
BSM-2300 lifescopes I BSM-4100 lifescopes P BSM-5100 lifescopes A BSM-9510 lifescopes M BSM-9800 lifescopes S		3	NK-2509 NK-2509-I		NKB3-90S NKB3-90S-I	NKB3-90P NKB3-90P-I	2309S 2309S-I	2309P 2309P-I	D-1309 D-1309-I
		Price	C8		B7	B7	C6	C6	C6
		6	NK-2509 NK-2509-I		NKB6-90S NKB6-90S-I	NKB6-90P NKB6-90P-I	2609S 2609S-I	2609P 2609P-I	
		Price	C8		C3	C3	C8	C8	

Nihon Kohden 下一页续接

Types of Major ECG lead wires 主要类型导联线	Lead No. 导联数	REF# Snap(AHA/IEC) 扣式(美标/欧标)	OEM REF# OEM原厂型号	Price code 代格代码	REF# Grabber(AHA/IEC) 钳式(美标/欧标)	OEM REF# OEM原厂型号	Price code 代格代码
<b>Datex Individual set leads (DX300 series) 分离式导联线 (DX300系列)</b>							
	3	DXB3-90S DXB3-90S-I	545327	B8	DXB3-90P DXB3-90P-I	545317 545315	B8
	5	DXB5-90S DXB5-90S-I	545328	C4	DXB5-90P DXB5-90P-I	545318 545316	C4
<b>GE Marquette/Multilink-integrated leads 集成式导联线</b>							
	3	MQ3-90S MQ3-90S-I	411203-001, 411203-002(130cm) 411203-003, 411203-004(130cm)	C2	MQ3-90P MQ3-90P-I	412682-001, 412682-002(130cm) 412682-003, 412682-004(130cm)	C2
	5	MQ5-90S MQ5-90S-I (Gray)	411202-001, 411202-002(130cm) 411200-003, 411200-004	C2	MQ5-90P MQ5-90P-I (Gray)	412681-001, 412681-002(130cm) 414556-003, 414556-004	C2
	5 V2-V6 C2-C6	MV5-90S MV5-90S-I (Brown)	411202-001, 411202-002(130cm) 411200-003, 411200-004	C2 C2	MV5-90P MV5-90P-I (Brown)	416467-001/ 416467-002 416467-003/ 416467-004	C2 C2
<b>GE Multilink-individual set leads 分离式导联线</b>							
	5	MQB5-90S MQB5-90S-I (Gray)	411200-001, 411200-002(130cm) 411200-003, 411200-004	C4	MQB5-90P MQB5-90P-I (Gray)	414556-001, 414556-002(130cm) 414556-003, 414556-004	C4
	5 V2-V6 C2-C6	MVB5-90S MVB5-90S-I (Brown)	416468-001(130cm), 416468-003, 416468-004, 416468-002	C4	MVB5-90P MVB5-90P-I (Brown)	416467-001(130cm), 416467-003 416467-004, 416467-002	C4
<b>Nihon Kohden Integrated leads 光电心电导联线</b>							
	3	NKA3-90S NKA3-90S-I	019P	B7	NKA3-90P NKA3-90P-I		B7
	5	NKA5-90S NKA5-90S-I		C2	NKA5-90P NKA5-90P-I		C2
	2	NKB2-90S NKB2-90S-I	BR-912PA BR-912P	NA	NKB2-90P NKB2-90P-I	BR-902PA BR-902P	NA
	3	NKB3-90S NKB3-90S-I	BR-913PA BR-913P	B7	NKB3-90P NKB3-90P-I	BR-903PA BR-903P	B7
	6	NKB6-90S NKB6-90S-I	BR-916PA BR-916P	C3	NKB6-90P NKB6-90P-I	BR-906PA BR-906P	C3
<b>Philips/HP Integrated leads 集成式导联线</b>							
	3	AA3-90S AA3-90S-I	M1605A M1615A	B7	AA3-90P AA3-90P-I	M1603A M1613A	B7
	5	AA5-90S AA5-90S-I	M1615A M1635A	C2	AA5-90P AA5-90P-I	M1623A M1633A	C2
<b>Philips/HP Individual set leads 分离式导联线</b>							
	3	AAB3-90S AAB3-90S-I	M1605A M1615A	B8	AAB3-90P AAB3-90P-I	M1603A M1613A	B8
	5	AAB5-90S AAB5-90S-I	M1615A M1635A	C4	AAB5-90P AAB5-90P-I	M1623A M1633A	C4

Types of Major ECG lead wires 主要类型导联线	Lead No. 导联数	REF# Snap(AHA/IEC) 扣式(美标/欧标)	OEM REF# OEM原厂型号	Price code 代格代码	REF# Grabber(AHA/IEC) 钳式(美标/欧标)	OEM REF# OEM原厂型号	Price code 代格代码
<b>Philips IntelliVue Integrated leads 心电导联线</b>							
	3	HP3-90S HP3-90S-I(Gray)	M1673A M1674A	C2	HP3-90P HP3-90P-I(Gray)	M1671A M1672A	C2
	5	HPA5-90S HPA5-90S-I(Blue)	M1644A M1645A	C4	HPA5-90P HPA5-90P-I(Blue)	M1968A M1971A	C4
	5 V2-V6 C2-C6	HPB5-90S HPB5-90S-I(red)	M1602A M1604A	C4	HPB5-90P HPB5-90P-I(red)	M1976A M1978A	C4
<b>Phyiso Control/Meditronic integrated leads 集成式导联线</b>							
	4	MDA4-120S MDA4-120S-I		C2			
	6	MDB6-90S MDB6-90S-I		C6			
<b>Siemens/VS,Drager Integrated leads 集成式导联线</b>							
	3	SM3-90S SM3-90S-I		B7	SM3-90P SM3-90P-I		B7
	5	SM5-90S SM5-90S-I		C2	SM5-90P SM5-90P-I		C2
	3 (5/3)	SMB3-90S SMB3-90S-I		C1	SMB3-90P SMB3-90P-I		C1
<b>Siemens/VS,Drager Individual set leads 分离式导联线</b>							
	3	SMD3-90S SMD3-90S-I		B8	SMD3-90P SMD3-90P-I		B8
	5	SMD5-90S SMD5-90S-I		C4	SMD5-90P SMD5-90P-I		C4
	3 (5/3)	SMF3-90S SMF3-90S-I		C1	SMF3-90P SMF3-90P-I		C1
<b>Drager Integrated leads 集成式导联线</b>							
	3	DG3-90S DG3-90S-I		C5	DG3-90P DG3-90P-I		C5
	5	DG5-90S DG5-90S-I		C5	DG5-90P DG5-90P-I		C5
	6	DG6-90S DG6-90S-I		C5	DG6-90P DG6-90P-I		C5
<b>Spacelabs Integrated leads 集成式导联线</b>							
	3	SL3-90S SL3-90S-I	700-0007-00,60cm 700-0007-02,100cm 700-0007-01,60cm 700-0007-03,100cm	B7	SL3-90P SL3-90P-I	700-0006-00,60cm 700-0006-02,100cm 700-0006-01,60cm 700-0006-03,100cm	B7
	5	SL5-90S SL5-90S-I	700-0007-08,60cm 700-0007-10,100cm 700-0007-09,60cm 700-0007-11,100cm	C2	SL5-90P SL5-90P-I	700-0006-08,60cm 700-0006-10,100cm 700-0006-09,60cm 700-0006-11,100cm	C2
<b>Spacelabs Individual set leads 分离式导联线</b>							
	3	SLB3-90S SLB3-90S-I	700-0007-00, 700-0007-02, 700-0007-01, 700-0007-03	B8	SLB3-90P SLB3-90P-I	700-0006-00, 700-0006-02, 700-0006-01, 700-0006-03	B8
	4	SLB4-90S SLB4-90S-I	012-0499-06, 60cm	C3	SLB4-90P SLB4-90P-I	700-0006-06 700-0006-07	C3
	5	SLB5-90S SLB5-90S-I	700-0007-08, 700-0007-10, 700-0007-09, 700-0007-11	C4	SLB5-90P SLB5-90P-I	700-0006-08, 700-0006-10, 700-0006-09, 700-0006-11	C4
	6	SLB6-90S SLB6-90S-I		C5	SLB6-90P SLB6-90P-I	700-0006-14 700-0006-15	C5





## ② Telemetry ECG Cables

### 心电遥测电缆


#### Specification 产品规格

Cable material 线材	TPU Jacket, Light gray TPU外壳, 灰白色
Connector & Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Lead wire length 导联线长度	0.9m

#### Compatibility 兼容性

Major brands ECG transmitter & Holter system  
主要品牌遥测及动态心电图系统

OEM & Models OEM原厂设备	Photo 图片	Lead# 导联数	Unimed REF# 尤迈型号	OEM PN 原厂型号	Price code 价格代码
GE Healthcare Apexpro FH Telemetry Patient net wireless telemetry systems		3	AP3-90S	394111-002	C6
			AP3-90S-I		C6
		5	AP5-90S	394111-006	C7
			AP5-90S-I		C7
		6	AP6-90S	394111-008	C8
			AP6-90S-I		C8
GE Healthcare ApexPro Telemetry		6	MQB6-90S	421930-001	C6
			MQB6-90S-I		C6
		6	MQB6-120S	421930-002	C6
			MQB6-120S-I		C6
	MQB6-120P	421930-002	C6		
Mindray TMS-6012, TEL-100 Mennen		5	TM5-90S		C7
			TM5-90S-I		C7
Philips/HP Medical Legacy M2600A/ M2601A, M2636A Transmitters		3	AT3-90S	M2590A	C7
			AT3-90P	M2591A	C7
			AT3-90S-I	M2594A	C7
			AT3-90P-I	M2595A	C7
		5	AT5-90S	M2592A	C8
			AT5-90P	M2593A	C8
	AT5-90S-I	M2596A	C8		
	AT5-90P-I	M2597A	C8		
Telemetry 3/5-lead combiner (U.S. only 仅美国区域)	1 pcs/package(线框)	3/5	Y05-T202	M2598A / M2599A	B8
Philips Medical IntelliVue Telemetry		3	HT3-90S	989803151991	C2
			HT3-90P	989803173141	C2
			HT3-90S-I	989803152001	C2
			HT3-90S-I	989803151981	C2
		5	HT5-90S	989803173151	C4
			HT5-90P	989803152071	C4
	HT5-90S-I	989803152081	C4		
	HT5-90P-I	989803152061	C4		
Shower Shield (防水罩)	10 pcs/package	3	Y05-T200	989803153031	C3
		5	Y05-T087	989803153041	C3
		6	Y05-T201	989803153051	C4
Single/double ECG alignment guide protective insert, untethered	10 pcs/package	1	Y05-T088	989803140401	A2
		2	Y05-T089	98980314042	A3
MGY美高仪/MGY-M1 蓝港/Visual Stethoscope 3.5mm Jack plug		3	H353-90S		B8
			H353-90P		B8
			H353-90S-I		B8
	H353-90P-I		B8		
Creative 科瑞康 PC-80 Easy ECG Monitor 2.5mm jack plug		3	H253-90S		B8
			H253-90P		B8
			H253-90S-I		B8
	H253-90P-I		B8		



- Full compatibility of the lead wires enable to be used with original ECG trunk Cables or devices, no extra adapter required.
- Built-in ribbon cable for patient comfort and ease of use
- Easy to read labeling facilities easy connection and lead placement
- Unique light blue color indicated the single use
- 全面兼容原装主电缆或仪器，无需其它转接线
- 可撕式带状电缆，舒适易用
- 标签清晰易读、连接放置简单
- 独特的浅蓝色插头设计

### 3.1 Disposable ECG lead wires 一次性监护分体心电导联线

## ③ Disposable ECG Monitoring lead wires 一次性心电导联线

Patient Safety  
安全

Patient Comfort  
舒适

Electrical Performance  
电气性能

Ease of Use  
易使用

Single patient use (disposable) from patient to monitor, Reduces the risk of cross contamination from reusable cables and leads, improving patient care and cost control

连接病人与监护仪的一次性心电导联线（供单个病人使用），减少了重复性电缆及导联线在病人间交叉感染的风险，改善了护理环境并有节省成本。

Descriptions 描述	Probe Connector 导联线接头端	Lead # 导联	REF# Snap(AHA/IEC) 扣式 (美标/欧标)	REF# Grabber(AHA/IEC) 钳式 (美标/欧标)	Packaging 包装	Unit code 单根价格
<b>Datascope</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联心电导联线，可分离排线，屏蔽，90厘米		3	DT3-90DS DT3-90DS-I	DT3-90DP DT3-90DP-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联心电导联线，可分离排线，屏蔽，90厘米		5	DT5-90DS DT5-90DS-I	DT5-90DP DT5-90DP-I	10	B2
<b>Drager /Siemens</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联心电导联线，可分离排线，屏蔽，90厘米		3	DG3-90DS DG3-90DS-I	DG3-90DP DG3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联心电导联线，可分离排线，屏蔽，90厘米		5	DG5-90DS DG5-90DS-I	DG5-90DP DG5-90DP-I	10	B3
Disposable 6-lead shielded Ribbon cable, 90cm 一次性6导联心电导联线，可分离排线，屏蔽，90厘米		6	DG6-90DS DG6-90DS-I	DG6-90DP DG6-90DP-I	10	B5
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联心电导联线，可分离排线，屏蔽，90厘米		3	SM3-90DS SM3-90DS-I	SM3-90DP SM3-90DP-I	10	B2
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联心电导联线，可分离排线，屏蔽，90厘米		5	SM5-90DS SM5-90DS-I	SM5-90DP SM5-90DP-I	10	B2

Descriptions 描述	Probe Connector 导联线接头端	Lead # 导联	REF# Snap(AHA/IEC) 扣式 (美标/欧标)	REF# Grabber(AHA/IEC) 钳式 (美标/欧标)	Packaging 包装	Unit code 单根价格
<b>GE Multi-link</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联, 可分离排线, 屏蔽, 90厘米		3	MQ3-90DS MQ3-90DS-I	MQ3-90DP MQ3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联导联线, 可分离排线, 屏蔽, 90厘米		5	MQ5-90DS MQ5-90DS-I	MQ5-90DP MQ5-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, V or C lead ,90cm 一次性5导联, 可分离排线, V 或C 导联, 屏蔽, 90厘米		5	MV5-90DS MV5-90DS-I	MV5-90DP MV5-90DP-I	10	B5
Disposable 6-lead shielded Ribbon cable for Apex, 90cm 一次性6导联导联线, 可分离排线, 用于APEX, 屏蔽, 90厘米		6	MQB6-90DS MQB6-90DS-I	MQB6-90DP MQB6-90DP-I	10	B6
<b>Mindray</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联, 可分离排线, 屏蔽, 90厘米		3	MR3-90DS MR3-90DS-I	MR3-90DP MR3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联导联线, 可分离排线, 屏蔽, 90厘米		5	MR5-90DS MR5-90DS-I	MR5-90DP MR5-90DP-I	10	B5
<b>Nihon Kohden</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联心电图导联线, 可分离排线, 屏蔽, 90厘米		3	NKB3-90DS NKB3-90DS-I	NKB3-90DP NKB3-90DP-I	10	B3
Disposable 6-lead shielded Ribbon cable, 90cm 一次性6导联, 可分离排线, 屏蔽, 90厘米		6	NKB6-90DS NKB6-90DS-I	NKB6-90DP NKB6-90DP-I	10	B5
<b>Philips IntelliVue</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联, 可分离排线, 屏蔽, 90厘米		3	HP3-90DS HP3-90DS-I	HP3-90DP HP3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联导联线, 可分离排线, 屏蔽, 90厘米		5	HPA5-90DS HPA5-90DS-I	HPA5-90DP HPA5-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, V or C lead ,90cm 一次性5导联, 可分离排线, V 或C 导联, 屏蔽, 90厘米		5	HPB5-90DS HPB5-90DS-I	HPB5-90DP HPB5-90DP-I	10	B5
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联心电图导联线, 可分离排线, 屏蔽, 90厘米		3	HT3-90DS HT3-90DS-I	HT3-90DP HT3-90DP-I	10	B4
Disposable 5-lead shielded Ribbon V or C-lead cable, 90cm 一次性5导联心电图导联线, 可分离排线, 屏蔽, V或C-导, 90厘米		5	HT5-90DS HT5-90DS-I	HT5-90DP HT5-90DP-I	10	B6



Descriptions 描述	Probe Connector 导联线接头端	Lead # 导联	REF# Snap(AHA/IEC) 扣式 (美标/欧标)	REF# Grabber(AHA/IEC) 钳式 (美标/欧标)	Packaging 包装	Unit code 单根价格
<b>Philips Twin Pin(AA style)</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联, 可分离排线, 屏蔽, 90厘米		3	AAB3-90DS AAB3-90DS-I	AAB3-90DP AAB3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联, 可分离排线, 屏蔽, 90厘米		5	AAB5-90DS AAB5-90DS-I	AAB5-90DP AAB5-90DP-I	10	B6
<b>Philips Veridia Telemetry</b>						
Disposable 3-lead shielded Ribbon cable, 90cm with combiner 一次性3导联导联线, 带线框, 可分离排线, 屏蔽, 90厘米		3	AT3-90DS AT3-90DS-I	AT3-90DP AT3-90DP-I	10	B5
Disposable 5-lead shielded Ribbon cable, 90cm with combiner 一次性5导联导联线, 带线框, 可分离排线, 屏蔽, 90厘米		5	AT5-90DS AT5-90DS-I	AT5-90DP AT5-90DP-I	10	B6
<b>Spacelabs</b>						
Disposable 3-lead shielded Ribbon cable, 90cm 一次性3导联, 可分离排线, 屏蔽, 90厘米	NA	3	SLB3-90DS SLB3-90DS-I	SLB3-90DP SLB3-90DP-I	10	B3
Disposable 5-lead shielded Ribbon cable, 90cm 一次性5导联导联线, 可分离排线, 屏蔽, 90厘米		5	SL5-90DS SL5-90DS-I	SL5-90DP SL5-90DP-I	10	B5

## 3.2 Disposable ECG One piece lead wires 一次性诊断心电导联线



OEM & Models 原厂型号	EKG machine Plug 机器接头端	Lead # 导联	REF# Snap(AHA/IEC) 扣式 (美标/欧标)	REF# Grabber(AHA/IEC) 钳式 (美标/欧标)	Packaging 包装	Unit code 单根价格
<b>AAMI</b>						
Merlin 78352C/78354C/ 78834CM1001A		3	2340DS 2340DS-I	2340DP 2340DP-I	8	B6
		5	2540DS 2540DS-I	2540DP 2540DP-I	8	B7
<b>GE Multi-link</b>						
Eagle,Solar,Dash monitor MAC-Lab Cath Lab Datex-Ohmeda M-Prestn & "E" modules		3	2386DS 2386DS-I	2386DP 2386DP-I	8	B5
		5	2586DS 2586DS-I	2586DP 2586DP-I	8	B6
<b>Mindray / Datascope</b>						
PM5000 PM6000 BeneView T5/T6/T8		3	2312DS 2312DS-I	2312DP 2312DP-I	8	B6
		5	2512DS 2512DS-I	2512DP 2512DP-I	8	B7
Datascope-Mindray V12, V21		3	2326DS 2326DS-I	2326DP 2326DP-I	8	B6
		5	2526DS 2526DS-I	2526DP 2526DP-I	8	B7
<b>Nihon Kohden</b>						
BSM-2300 iflescope I BSM-4100 iflescope P BSM-5100		3	2309DS 2309DS-I	2309DP 2309DP-I	8	B6
		6	2609DS 2609DS-I	2609DP 2609DP-I	8	B7
<b>Philips</b>						
Merlin 78352C/78354C/ 78834CM1001A/		3	2385DS 2385DS-I	2385DP 2385DP-I	8	B6
		5	2585DS 2585DS-I	2585DP 2585DP-I	8	B7
<b>Spacelabs</b>						
Ultraview 90492, 91220, 90496, 91496		3	2396DS 2396DS-I	2396DP 2396DP-I	8	B6
		5	2596DS 2596DS-I	2596DP 2596DP-I	8	B7



## ④ ECG Lead Adapter 心电导联转换头

The unique design of the ECG lead wires adapters enable the hospitals to use the one type lead wires to fit major brands patient monitoring system running in one hospital by attaching the proper ECG adapter to the trunk cable. No confusion on lead wires any more!


心电导联线转换器的独特设计，让医院实现了用一种导联线适配各种主要品牌医用监护系统的可能。只须让合适的转换器连上主电缆，一切关于导联线的选用困扰都将不复存在！

### Specification 产品规格

Cable material 线材	ABS Jacket, Light gray housing ABS, 浅灰色外壳
Connector & Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Link Tie 绑带	12cm length with trunk cable clip 长12cm, 主电缆固定带

### Compatibility 兼容性

Major brands ECG monitoring leads system  
主要品牌心电监护导联线



OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
<b>Convert To Din style lead wires D型导联线转换头</b>								
Colin	AD-CL3-D3			B5	AD-CL5-D5			B5
GE Marquette	AD-MQ5-D3			B5	AD-MQ5-D5			B5
GE Datex-Ohmeda	AD-DX3-D3			B5	AD-DX5-D5			B5
Nihon Kohden	AD-NKA3-D3			B5	AD-NKA5-D5			B5
Nihon Kohden	AD-NKB3-D3			B5				
Philips medical	AD-AA3-D3			B5				
Philips medical	AD-HP3-D3			B5				
Siemens/Drager	AD-SM3-D3			B5	AD-SM5-D5			B5
Spacelabs	AD-SL3-D3			B5	AD-SL5-D5			B5
<b>Convert To VS/Siemens lead VS/Siemens导联线转换头</b>								
Colin	AD-CL3-SM3			B5	AD-CL5-SM5			B5
GE Marquette	AD-MQ5-SM3			B5	AD-MQ5-SM5			B5
GE Datex-Ohmeda	AD-DX3-SM3			B5	AD-DX5-SM5			B5
Nihon Kohden	AD-NKA3-SM3			B5	AD-NKA5-SM5			B5
Nihon Kohden	AD-NKB3-SM3			B5				
Philips medical	AD-AA3-SM3			B5				
Philips medical	AD-HP3-SM3			B8				
Spacelabs	AD-SL3-SM3			B5	AD-SL5-SM5			B5

OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
<b>Convert To Philips/HP style lead Philips/HP型转换头</b>								
Colin	AD-CL3-HP3			B6	AD-CL5-HP5			B7
GE Marquette	AD-MQ5-HP3			B6	AD-MQ5-HP5			B7
GE Datex-Ohmeda	AD-DX3-HP3			B6	AD-DX5-HP5			B7
Nihon Kohden	AD-NKA3-HP3			B6	AD-NKA5-HP5			B7
Nihon Kohden	AD-NKB3-HP3			B6				
Philips medical	AD-AA3-HP3			B6				
Siemens/Drager	AD-SM3-HP3			B6	AD-SM5-HP5			B7
Spacelabs	AD-SL3-HP3			B6	AD-SL5-HP5			B7
<b>Convert To GE Multi-Link lead GE Multi-Link转换头</b>								
Colin					AD-CL5-MQ5			B5
GE Datex-Ohmeda					AD-DX5-MQ5			B5
Nihon Kohden					AD-NKA5-MQ5			B5
Siemens/Drager					AD-SM5-MQ5			B5
Spacelabs					AD-SL5-MQ5			B5

Application Figure 应用图示





Datex 300 style lead



Spacelabs style lead



OEM Manufacturer OEM 制造商	3-lead REF# 3-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price 价格	5-lead REF# 5-导联型号	Trunk Cable end 主电缆端接头	Lead wire end 导联线端接头	Price code 价格代码
Convert To Datex 300 style lead								
Colin	AD-CL3-DX3			B5	AD-CL5-DX5			B5
GE Marquette				B5	AD-MQ5-DX5			B5
Nihon Kohden	AD-NKA3-DX3			B5	AD-NKA5-DX5			B5
Nihon Kohden	AD-NKB3-DX3			B5				
Philips medical	AD-AA3-DX3			B5				
Philips medical	AD-HP3-DX3			B6				
Siemens/Drager	AD-SM3-DX3			B5	AD-SM5-DX5			B5
Spacelabs	AD-SL3-DX3			B5	AD-SL5-DX5			B5
Convert To Spacelabs style lead								
Colin	AD-CL3-SL3			B5	AD-CL5-SL5			B5
GE Marquette	AD-MQ5-SL3			B5	AD-MQ5-SL5			B5
GE Datex-Ohmeda	AD-DX3-SL3			B5	AD-DX5-SL5			B5
Nihon Kohden	AD-NKA3-SL3			B5	AD-NKA5-SL5			B5
Nihon Kohden	AD-NKB3-SL3			B5				
Philips medical	AD-AA3-SL3			B5				
Philips medical	AD-HP3-SL3			B6				
Spacelabs	AD-SM3-SL3			B5	AD-SM5-SL5			B5

## ⑤ Diagnosis ECG Cables & Accessories

### 诊断用心电电缆及附件


#### Specification 产品规格

Cable material 线材	TPU Jacket, Light gray TPU外被, 浅灰
Connector & Yoke 接头和分线器	Nylon plug & Golden plated pins 尼龙插头和镀金插针
Trunk Cable length 主电缆长度	2m
Lead wire length 导联线长度	1.2m(limb lead);0.8m(Chest lead) 1.2m(四肢); 0.8m(胸导)

#### Compatibility 兼容性

Major brands ECG diagnosis system  
主要品牌的心电图机, 运动平板等心电图诊断设备

Application Diagram  
应用示意图



Product name 产品类型	Photo 图片	REF# Example 型号#例举	Reference to REF# Descriptions 型号命名参考说明
One piece Fixed Cable 一体化电缆		E10R-HP-B/I	E: EKG diagnostic cable or lead wire 10: 10-lead(12=12 lead) R: built in resistors for protection(default: no resistors) HP: Philips/HP, factory/model code 1: Connector configure; B: banana end(N= needle, P=grabber, S=snap) I: IEC code (Default is AHA code)
Trunk Cable 主电缆		E10R-MQ-I	E: EKG diagnostic cable or lead wire 10: 10-lead(12=12 lead) R: built in resistors for protection(default: no resistors) MQ: GE Marquette, factory/model code 1: Connector configure (default: only one connector available) I: IEC code (Default is AHA code)
Lead Wire 导联线		MQ10-LB-I	MQ: Marquette, factory/model code 10: 10-lead(12=12 lead) L: lead wires B: banana end(N= needle, P=grabber, S=snap) I: IEC code (Default is AHA code)

All one piece cable or lead wires are available with banana end, needle end, grabber end and snap end  
供应各式一体化电缆和导联线 (香蕉头/针式/钳式/扣式)

Code 代码	AHA(A) 美标	Standard IEC(I) 欧标	Alternate IEC(IH) 替代欧标
4mm Banana end 香蕉头 (B)			
3mm Needle end 针式 (N)			
Grabber End 钳式 (P)			
Snap End 扣式 (S)			

# 5.1 One piece EKG fixed cable

## 一体化心电电缆



OEM & Models 原厂及型号	EKG machine Plug 机器端插头	Lead# 导联数	With 4mm banana 一体化香蕉头4mm	With 3mm needle 一体化针头3mm	With grabber end 一体化钳式	With Snap end 一体化扣式	Price code 价格代码
<b>Burdick/Cardiac Science</b>							
EK-10 and Elite, E350, E350i, E550, E560, Eclipse™ LE/LE.		10	E10-BK1-B E10R-BK1-B	E10-BK1-N E10R-BK1-N	E10-BK1-P E10R-BK1-P	E10-BK1-S E10R-BK1-S	D3
		10	E10-BK1-B/I E10R-BK1-B/I	E10-BK1-N/I E10R-BK1-N/I	E10-BK1-P/I E10R-BK1-P/I	E10-BK1-S/I E10R-BK1-S/I	D3
Burdick Quinton Eclipse 4, 8, 400, 800, 850, 850i, plus, Eclipse LE II and ATRIA PLUS Quinton		10	E10-BK2-B E10R-BK2-B	E10-BK2-N E10R-BK2-N	E10-BK2-P E10R-BK2-P	E10-BK2-S E10R-BK2-S	D3
		10	E10-BK2-B/I E10R-BK2-B/I	E10-BK2-N/I E10R-BK2-N/I	E10-BK2-P/I E10R-BK2-P/I	E10-BK2-S/I E10R-BK2-S/I	D3
<b>Fukuda Denshi</b>							
FX101, FX102F, FX328U, FCP2201		10	E10-FD1-B E10R-FD1-B	E10-FD1-N E10R-FD1-N	E10-FD1-P E10R-FD1-P	E10-FD1-S E10R-FD1-S	D3
		10	E10-FD1-B/I E10R-FD1-B/I	E10-FD1-N/I E10R-FD1-N/I	E10-FD1-P/I E10R-FD1-P/I	E10-FD1-S/I E10R-FD1-S/I	D3
<b>Fukuda ME</b>							
KP-500, KP-500D		10	E10-FD2-B E10R-FD2-B	E10-FD2-N E10R-FD2-N	E10-FD2-P E10R-FD2-P	E10-FD2-S E10R-FD2-S	D3
		10	E10-FD2-B/I E10R-FD2-B/I	E10-FD2-N/I E10R-FD2-N/I	E10-FD2-P/I E10R-FD2-P/I	E10-FD2-S/I E10R-FD2-S/I	D3

OEM & Models 原厂及型号	EKG machine Plug 机器端插头	Lead# 导联数	With 4mm banana 一体化香蕉头4mm	With 3mm needle 一体化针头3mm	With grabber end 一体化钳式	With Snap end 一体化扣式	Price code 价格代码
<b>GE-Marquette</b>							
MAC 500, MAC 1100, MAC 1200		10	E10-MQ-B E10R-MQ-B	E10-MQ-N E10R-MQ-N	E10-MQ-P E10R-MQ-P	E10-MQ-S E10R-MQ-S	D3
		10	E10-MQ-B/I E10R-MQ-B/I	E10-MQ-N/I E10R-MQ-N/I	E10-MQ-P/I E10R-MQ-P/I	E10-MQ-S/I E10R-MQ-S/I	D3
<b>Hellige / Siemens Hormann / Bosch</b>							
		10	E10-HL-B E10R-HL-B	E10-HL-N E10R-HL-N	E10-HL-P E10R-HL-P	E10-HL-S E10R-HL-S	D3
		10	E10-HL-B/I E10R-HL-B/I	E10-HL-N/I E10R-HL-N/I	E10-HL-P/I E10R-HL-P/I	E10-HL-S/I E10R-HL-S/I	D3
<b>Kenz</b>							
103,106 CARDIOLINE (Italy)		10	E10-KZ1-B E10R-KZ1-B	E10-KZ1-N E10R-KZ1-N	E10-KZ1-P E10R-KZ1-P	E10-KZ1-S E10R-KZ1-S	D3
		10	E10-KZ1-B/I E10R-KZ1-B/I	E10-KZ1-N/I E10R-KZ1-N/I	E10-KZ1-P/I E10R-KZ1-P/I	E10-KZ1-S/I E10R-KZ1-S/I	D3
ECG 108/110 Cardico 302/601/ 1210/1211		10	E10R-KZ2-B	E10R-KZ-N	E10R-KZ2-P	E10R-KZ2-S	D3
		10	E10R-KZ2-B/I	E10R-KZ2-N/I	E10R-KZ2-P/I	E10R-KZ2-S/I	D3
<b>Mortara</b>							
ELI-200, ELI-50		10	E10-MT-B E10R-MT-B	E10-MT-N E10R-MT-N	E10-MT-P E10R-MT-P	E10-MT-S E10R-MT-S	D3
		10	E10-MT-B/I E10R-MT-B/I	E10-MT-N/I E10R-MT-N/I	E10-MT-P/I E10R-MT-P/I	E10-MT-S/I E10R-MT-S/I	D3
<b>Nihon Kohden</b>							
CadioFax 6151, 6851,6353 CadioFax 8110		10	E10-NK1-B E10R-NK1-B	E10-NK1-N E10R-NK1-N	E10-NK1-P E10R-NK1-P	E10-NK1-S E10R-NK1-S	D3
		10	E10-NK1-B/I E10R-NK1-B/I	E10-NK1-N/I E10R-NK1-N/I	E10-NK1-P/I E10R-NK1-P/I	E10-NK1-S/I E10R-NK1-S/I	D3
CadioFax Q ECG-9110K ECG-9130K Doing Juang Cardioline (New)		10	E10-NK2-B E10R-NK2-B	E10-NK2-N E10R-NK2-N	E10-NK2-P E10R-NK2-P	E10-NK2-S E10R-NK2-S	D3
		10	E10-NK2-B/I E10R-NK2-B/I	E10-NK2-N/I E10R-NK2-N/I	E10-NK2-P/I E10R-NK2-P/I	E10-NK2-S/I E10R-NK2-S/I	D3
<b>Philips/HP</b>							
M1770A, M1771A, M1772A Pagewrite100/200/ 200i		10	E10-HP-B E10R-HP-B	E10-HP-N E10R-HP-N	E10-HP-P E10R-HP-P	E10-HP-S E10R-HP-S	D3
		10	E10-HP-B/I E10R-HP-B/I	E10-HP-N/I E10R-HP-N/I	E10-HP-P/I E10R-HP-P/I	E10-HP-S/I E10R-HP-S/I	D3
<b>Schiller</b>							
AT1, AT2, AT2 Plus, AT4 Resistor 10K Long Screw		10	E10-SH1-B E10R-SH1-B	E10-SH1-N E10R-SH1-N	E10-SH1-P E10R-SH1-P	E10-SH1-S E10R-SH1-S	D3
		10	E10-SH1-B/I E10R-SH1-B/I	E10-SH1-N/I E10R-SH1-N/I	E10-SH1-P/I E10R-SH1-P/I	E10-SH1-S/I E10R-SH1-S/I	D3
AT3, AT6, CS6, AT5, AT10, AT60 Short Screw		10	E10-SH2-B E10R-SH2-B	E10-SH2-N E10R-SH2-N	E10-SH2-P E10R-SH2-P	E10-SH2-S E10R-SH2-S	D3
		10	E10-SH2-B/I E10R-SH2-B/I	E10-SH2-N/I E10R-SH2-N/I	E10-SH2-P/I E10R-SH2-P/I	E10-SH2-S/I E10R-SH2-S/I	D3



## 5.2 EKG Trunk Cable

### 心电主电缆









OEM& Models 原厂及型号	EKG machine Plug 机器端插头	Lead# 导联数	REF# 型号	Used with lead wire set/GE Multi-link set 与导联线/GE Multi-link导联线一起使用	Price code 价格代码
<b>Burdick</b>					
EK10-E350		10	<b>E10-BK1</b> <b>E10R-BK1</b>	10-lead, chest /0.9m;Limb/1.2m MQ10-LB,MQ10-LN,MQ10-LP,MQ10-LS or MQB10-LB, ...	D3
		10	<b>E10-BK1-I</b> <b>E10R-BK1-I</b>	10-lead, chest /0.9m;Limb/1.2m MQ10-LB-I,MQ10-LN-I,MQ10-LP-I,MQ10-LS-I or MQB10-LB, ... (See the next page of 10-lead wire sets for detail) (详细细节见下页10-导联线)	D3
Burdick Quinton Eclipse 4,8,400,800, 850,850i,plus, Eclipse LE II and ATRIA PLUS		10	<b>E10-BK2</b> <b>E10R-BK2</b>	Same as above 同上	D3
		10	<b>E10-BK2-I</b> <b>E10R-BK2-I</b>	Same as above 同上	D3
<b>Fukuda Denshi</b>					
FX101, FX102F, FX328U, FCP2201		10	<b>E10-FD1</b> <b>E10R-FD1</b>	Same as above 同上	D3
		10	<b>E10-FD1-I</b> <b>E10R-FD1-I</b>	Same as above 同上	D3
<b>Fukuda ME</b>					
KP-500、 KP-500D		10	<b>E10-FD2</b> <b>E10R-FD2</b>	Same as above 同上	D3
		10	<b>E10-FD2-I</b> <b>E10R-FD2-I</b>	Same as above 同上	D3

OEM& Models 原厂及型号	EKG machine Plug 机器端插头	Lead# 导联数	REF# 型号	Used with lead wire set/GE Multi-link set 与导联线/GE Multi-link导联线一起使用	Price code 价格代码
<b>GE-Marquette</b>					
MAC 500, MAC 1100, MAC 1200		10	<b>E10-MQ</b> <b>E10R-MQ</b>	Same as above 同上	D3
		10	<b>E10-MQ-I</b> <b>E10R-MQ-I</b>	Same as above 同上	D3
ELI-200、 ELI-50		10	<b>E10-MT</b> <b>E10R-MT</b>	Same as above 同上	D3
		10	<b>E10-MT-I</b> <b>E10R-MT-I</b>	Same as above 同上	D3
<b>Hellige / Siemens Hormann / Bosch</b>					
		10	<b>E10-HL</b> <b>E10R-HL</b>	Same as above 同上	D3
		10	<b>E10-HL-I</b> <b>E10R-HL-I</b>	Same as above 同上	D3
<b>Kenz</b>					
103,106 CARDIOLINE (Italy)		10	<b>E10-KZ1</b> <b>E10R-KZ1</b>	Same as above 同上	D3
		10	<b>E10-KZ1-I</b> <b>E10R-KZ1-I</b>	Same as above 同上	D3
ECG 108/110 Cardico 302/601/ 1210/1211		10	<b>E10R-KZ2</b>	Same as above 同上	D3
		10	<b>E10R-KZ2-I</b>	Same as above 同上	D3
<b>Nihon Kohden</b>					
CaduoFax 6151, 6851,6353 CaduoFax 8110		10	<b>E10-NK1</b> <b>E10R-NK1</b>	Same as above 同上	D3
		10	<b>E10-NK1-I</b> <b>E10R-NK1-I</b>	Same as above 同上	D3
CaduoFax Q ECG-9110K ECG-9130K Doing Juang Cardioline (New)		10	<b>E10-NK2</b> <b>E10R-NK2</b>	Same as above 同上	D3
		10	<b>E10-NK2-I</b> <b>E10R-NK2-I</b>	Same as above 同上	D3
<b>Philips/HP</b>					
M1770A,M1771A, M1772A Pagewrite100/200/ 200i		10	<b>E10-HP</b> <b>E10R-HP</b>	Same as above 同上	D3
		10	<b>E10-HP-I</b> <b>E10R-HP-I</b>	Same as above 同上	D3
<b>Schiller</b>					
AT1, AT2, AT2 Plus,AT4 Resistor 10K Long Screw		10	<b>E10-SH1</b> <b>E10R-SH1</b>	Same as above 同上	D3
		10	<b>E10-SH1-I</b> <b>E10R-SH1-I</b>	Same as above 同上	D3
AT3,AT6,CS6,AT5, AT10,AT60 Short Screw		10	<b>E10-SH2</b> <b>E10R-SH2</b>	Same as above 同上	D3
		10	<b>E10-SH2-I</b> <b>E10R-SH2-I</b>	Same as above 同上	D3
<b>GE Healthcare-Resting sys.</b>					
CAM14 coiled patient Cable Used with GE MAC 5000 resting system with CAM14 Acquisition Module		10	<b>ECAM-GE14-S</b> 2016560-001 1.4m+0.3m <b>ECAM-GE14</b> 2016560-003 4.5m+0.65m	Used with CAM 14 patient lead wires on next page. 与下页的CAM 14导联线配合使用	D8
		10	<b>EAM-GE4-S</b> 700044-203 1.3m+0.25m <b>EAM-GE4</b> 700044-202 4.6m+0.25m	Used with CAM 4 patient lead wires on next page. 与下页的CAM 4导联线配合使用	D8


## 5.3 EKG Lead wire 心电导联线

Pictures 产品图片	OEM & Models 原厂型号	Electrode end 病人电极端	REF# (AHA) 型号 (美标)	REF# (IEC) 型号 (欧标)	Pieces/Set 根/套	Price code 价格代码
<b>GE Multi-Link 10-Lead-integrated leads 集成式导联线</b>						
	Banana End 900112-001(AHA) 900112-002(IEC)	4mm banana	MQ10-LB 900112-001	MQ10-LB-I 900112-002	10	D1
		3mm needle	MQ10-LN	MQ10-LN-I	10	D1
		Grabber	MQ10-LP	MQ10-LP-I	10	D1
		Snap	MQ10-LS	MQ10-LS-I	10	D1
<b>GE CAM 14 patient lead wire 导联线</b>						
		4mm banana	MQB10-LB	MQB10-LB-I	10	D1
		3mm needle	MQB10-LN	MQB10-LN-I	10	D1
		Grabber	MQB10-LP	MQB10-LP-I	10	D1
		Snap	MQB10-LS	MQB10-LS-I	10	D1
	420101-002	CAM/Universal	MQB10-AD	MQB10-AD-I	10	D1
	420101-001 Used with 900178-001/ 900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	MQB14-AD	MQB14-AD-I	14	D3
	2001925-005(66cm)	CAM/Universal	X-MQB-66DF	66cm,no color coded	10	C8
	2001925-002(74cm)	CAM/Universal	X-MQB-80DF	80cm,no color coded	10	C8
	2001925-006(90cm)	CAM/Universal	X-MQB-90DF	90cm,no color coded	10	C8
	2001925-003(102cm)	CAM/Universal	X-MQB-100DF	100cm,no color coded	10	C8
	2001925-004(130cm)	CAM/Universal	X-MQB-120DF	120cm,no color coded	10	C8
	2001925-003(102cm)	CAM/Universal	X-MQB-80DF	80cm, no color coded	10	C8
	2001925-004(130cm)	CAM/Universal	X-MQB-120DF	120cm, no color coded	10	C8
<b>GE AM4/AM5 patient lead wires 导联线</b>						
	Banana End 900112-001(AHA) 900112-002(IEC)	4mm banana	MQD10-LB	MQD10-LB-I	10	D1
		3mm needle	MQD10-LN	MQD10-LN-I	10	D1
		Grabber	MQD10-LP	MQD10-LP-I	10	D1
		Snap	MQD10-LS	MQD10-LS-I	10	D1
	900177-001 Used with 900178-001/ 900178-002 or Unimed Adapter as electrodes end (see below)	CAM/Universal	MQD12-AD	MQD12-AD-I	12	D3
	406581-004(90cm)	CAM/Universal	X-MQD-80DF	80cm, no color coded	10	C8
	406581-006(130cm)	CAM/Universal	X-MQD-120DF	120cm, no color coded	10	C8
<b>GE CAM14/AM4/AM5 Electrode Adapter 电极转接头</b>						
	900178-001,AHA 900178-101,IEC 900179-201,AHA 900178-201,IEC For the above leads set or GE original CAM14/AM4/AM5	4mm banana, 10 pcs/set	AD10-DB	AD10-DB-I	10	B8
		4mm banana, 14 pcs/set	AD14-DB	AD14-DB-I	14	C3
		4mm banana, 12 pcs/set	AD12-DB	AD12-DB-I	12	C2

Pictures 产品图片	OEM & Models 原厂型号	Electrode end 病人电极端	REF# (AHA) 型号 (美标)	REF# (IEC) 型号 (欧标)	Pieces/Set 根/套	Price code 价格代码
<b>Philips Lead wire 导联线</b>						
	PageWriter Touch: 860284, PageWriter Trim Cardiographs, 860286, 860287, 860288, 860290, 860297	4mm banana	PW10-LB 989803129161	PW10-LB-I 989803129191	10	D1
		3mm needle	PW10-LN	PW10-LN-I	10	D1
		Grabber	PW10-LP	PW10-LP-I	10	D1
		Snap	PW10-LS	PW10-LS-I	10	D1
	989803151631/989803151641	4mm banana	PWC10-LB	PWC10-LB-I	10	D1
		3mm needle	PWC10-LN	PWC10-LN-I	10	D1
		Grabber	PWC10-LP	PWC10-LP-I	10	D1
		Snap	PWC10-LS	PWC10-LS-I	10	D1
<b>QUINTON Lead wire 导联线</b>						
	Q-710 Stress Test System 031257-002	4mm banana	QT10-LB	QT10-LB-I	10	D1
		3mm needle	QT10-LN	QT10-LN-I	10	D1
		Grabber	QT10-LP	QT10-LP-I	10	D1
		Snap	QT10-LS	QT10-LS-I	10	D1
<b>Mortara Lead wire 导联线</b>						
	9293-046-60 ELI150C ELI230 ELI250C ELI350	4mm banana	MT10-LB	MT10-LB-I	10	D1
		3mm needle	MT10-LN	MT10-LN-I	10	D1
		Grabber	MT10-LP	MT10-LP-I	10	D1
		Snap	MT10-LS	MT10-LS-I	10	D1
<b>Norav (AA style) Lead wire 导联线</b>						
		4mm banana	AA10-LB	AA10-LB-I	10	D1
		3mm needle	AA10-LN	AA10-LN-I	10	D1
		Grabber	AA10-LP	AA10-LP-I	10	D1
		Snap	AA10-LS	AA10-LS-I	10	D1
<b>Welch Allyn Lead wire 导联线</b>						
	CP100/CP200 401129	4mm banana	WA10-LB	WA10-LB-I	10	D1
		3mm needle	WA10-LN	WA10-LN-I	10	D1
		Grabber	WA10-LP	WA10-LP-I	10	D1
		Snap	WA10-LS	WA10-LS-I	10	D1

# 5.4 EKG Accessories

## 附件



ECG Electrodes 心电极	REF# 型号	Package 包装	Price code 价格代码	REF to Figures 图例
Adult Limb Clamps/IEC, universal 成人四肢夹/欧标, 通用	UU-LC-A	4 pcs/box	B6	A
Adult Limb Clamps/AHA, all clamps in green 成人四肢夹/美标, 绿色	UU-LC-AA	4 pcs/box	B6	E
Pediatric Limb Clamps/IEC, universal 儿童四肢夹/欧标, 通用	UU-LC-P	4 pcs/box	B6	B
Pediatric Limb Clamps/AHA, all clamps in green 儿童四肢夹/美标, 绿色	UU-LC-PA	4 pcs/box	B6	未显示
Adult up Limb Clamps, for snap electrode, 3-lead 成人上肢夹, 接扣式电极, 3极	UU-DLC-A	2 pcs/box	B6	F
Adult suction electrodes, Universal 成人胸吸球, 通用	UU-SE-A	6 pcs/box	B8	C
Pediatric suction electrodes, Universal 儿童胸吸球, 通用	UU-SE-P	6 pcs/box	B8	D
Adult blue bulb 成人蓝色吸球	Y05-F027	6 pcs/box	A8	G
Pediatric green bulb 儿童绿色吸球	Y05-F022	6 pcs/box	A8	H

ECG Electrodes 心电极	REF# 型号	Package 包装	Price code 价格代码	REF to Figures 图例
* Flat clip/red Adapter/w 4mm banana 平式转换头/红色/ 4mm 香蕉头	AD-BF-10	10 pcs/bag	B7	A
* Flat clip/blue adapter/ w 3mm needle 平式转换头/蓝色/ 3mm 针头	AD-NF-10	10 pcs/bag	B7	B
Alligator clip/red Adapter/w 4mm banana 鳄鱼夹/ 红色/ 4mm 香蕉头	AD-BA-10	10 pcs/bag	B7	C
Alligator clip/blue adapter/w 3mm needle 鳄鱼夹/ 蓝色/ 3mm 针头	AD-NA-10	10 pcs/bag	B7	D
Adapter/4mm with grabber end 转接头/ 4mm 钳式	AD-BP-10	10 pcs/bag	B8	E
Adapter/4mm with snap end 转接头/ 4mm 扣式	AD-BS-10	10 pcs/bag	B8	F
Wide mouth Multi-function adapter/4mm with banana end 多功能心电转换钳/ 4mm 香蕉头	AD-BW-10	10 pcs/bag	B8	I
Wide mouth Multi-function adapter/3mm with banana end 多功能心电转换钳/ 3mm 香蕉头	AD-NW-10	10 pcs/bag	B8	J
Adapter wire/4mm with grabber end, 0.3m 转接线/ 4mm/香蕉头>>钳式	AD-B-P	10 pcs/bag	B8	G
Adapter wire/3mm with grabber end, 0.3m 转接线/ 3mm/针式>>钳式	AD-N-P	10 pcs/bag	B8	G
Adapter wire/4mm with snap end, 0.3m 转接线/ 4mm/香蕉头>扣式	AD-B-S	10 pcs/bag	B8	H
Adapter wire/3mm with snap end, 0.3m 转接线/ 3mm/针式>>扣式	AD-N-S	10 pcs/bag	B8	H

Items indicated “\*” are not sold in the USA.  
带\*的产品均不在美国区域销售

# Content



## 内容

- 1. Reusable Spo2 Sensor 重复用血氧探头 ..... 43
- 2. Disposable Spo2 Sensor 一次性血氧探头 ..... 73
- 3. Interface Cable 血氧转接线 ..... 77
- 4. Accessories 血氧配件 ..... 84

# SpO<sub>2</sub> Sensor & Interface Cable

## 血氧探头与连接电缆

### Certification 产品认证



### Compliance with 符合标准

ISO EN 9919-2005  
ISO10993-1,5,10:2009E  
EN 60601-1:1990



# ① Reusable SpO<sub>2</sub> Sensor


## 重复性血氧探头

### Specification 产品规格

Sensor Connector 探头插头端	OEM compatible, Nylon adapter with gold plated pins 尼龙插头、镀金插针、与原装兼容
LED & Photo Detector 发光管&接收管	660nm/905nm/940nm plastic or ceramic frame LED; 4mm or 8mm size plastic flame PD 660nm/905nm/940nm塑料或陶瓷发光管; 4mm 或 8mm 尺寸接收管
Cable material & color 线材材质&颜色	TPU jacket; light gray or bright blue TPU 外皮; 灰白色或深蓝色
Cable length & diameter 线长&直径	1.1m/1.5m/3.0m; φ 3.20mm/4.00mm or 5.00x2.50mm two line cable 线长1.1m/1.5m/3.0m; 直径3.20mm/4.00mm 或 5.00x2.50mm 双并线
Sensor house 探头外壳	ABS shell(finger clip); Silicon rubber boots( soft and multi-site) ABS 外壳(手指夹); 硅胶软指套(软指套和多功能探头)

### Compatible with major patient monitor brands 兼容主要品牌

Nellcor (耐普特) Criticare  
Ohmeda (欧美达) BCI  
Philips (飞利浦) Datascope  
Marquette (马奎) Datex德恩  
Spacelabs (太空) Nihon Kohden (光电)  
Novamatrix (诺医) Palco (派克)  
.....



Series系列	U400	U100	U400S	U100S	U200S	U300
REF to Fig	A	B	C	D	E	F
Name 名称	Adult Finger 成人手指夹	Pediatric Finger 儿童手指夹	Adult soft tip 成人软指套	Pediatric Soft tip 儿童软指套	Infant Soft tip 婴儿软指套	Neonate Wrap 新生儿捆绑
Patient size 适用范围	+40kg	15-40kg	+40kg	15-40kg	3-15kg	1-3kg
Cable color 线材颜色	Gray Cable 灰白色线	Gray Cable 灰白色线	Blue Cable 蓝色线	Gray two line 灰白双并线	Gray two line 灰白双并线	Gray two line 灰白双并线
Series系列	U800	U900	U600	U600B	U800E	U800V
REF to Fig	G	H	I	J	K	L
Name 名称	Multi-site Y型多功能	Adult Ear clip 耳夹式	Adult/Neo. Wrap 成人/新生儿包裹式	Pedia/Inf. Wrap 儿童/婴儿包裹式	Multi-site Y型多功能	Forehead/Rectal 额贴/直肠
Patient size 适用范围	All patients 所有病人	+40kg	+40kg or -3kg	15-40kg	All patients 所有病人	Adult 成人
Cable color 线材颜色	Gray two line 灰白双并线		Gray two line 灰白双并线	Gray two line 灰白双并线	Sensormate Used with DPT tapes 仅与贴片式绑带配合使用	Gray Cable 灰白色线

## How to identify the REF# of Unimed Spo2 sensors 如何识别尤迈血氧探头编号

### Reference Number Rule Description

#### 产品型号说明
















U403S-01XE	
U	Unimed brand 尤迈品牌
4	adult sensor (1: pediatric; 2: infant; 3: neonate; 5: disposable; 6 wrap; 7: ext-cable; 8: multi-site; 9: ear clip) 成人探头 (1: 儿童; 2: 婴儿; 3: 新生儿; 5: 一次性; 6: 包裹式; 7: 转接线; 8: 多功能; 9: 耳夹式)
03	1.1m cable length (03: 1.1m; 05: 1.5m; 08: 2.2m; 10: 3.0m) 1.1米线长 (03: 1.1米; 05: 1.5米; 08: 2.2米; 10: 3.0米)
S	Soft sensor (B: U style shell; Default: shell sensor or neonate wrap, alternative wrap, multi-site) 软指套探头 (B: U型外壳; 默认(缺省): 指夹探头, 新生儿捆绑, 包裹式, Y型多功能)
01	Nellcor code (02: Ohmeda; 03: Novametrix; 04: Spacelabs; 05: Criticare; 06: BCI...) Nellcor代码 (02: Ohmeda; 03: Novametrix; 04: Spacelabs; 05: Criticare; 06: BCI...)
X	Oximax technology (M: M technology, D: digital technology; Default: standard) 技术 (M: M技术, D: 数字技术; 默认(缺省): 标准型)
E	Sensormate sensor (Default: non-Sensormate sensor) Sensormate探头 (默认: 非Sensormate探头)



#### Cross Reference 接头对照表

Ending "X" can be replaced by "XP" or "P" for fully compatible with all oximax oximetry  
以 "X" 结尾的型号, 与对应 "XP"、"P" 结尾的型号通用

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>3F Medical (China)</b>								
IRIS 300/IRIS 400 Phoebe PMS8300/8310 PhoebeProPMS8300/ 8310PMS8110A			U410-189	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-189	Adult soft 成人软指套/3.0m	C8			
			U810-189	Multi-site Y型多功能/3.0m	C7			
<b>Adecon (China)</b>								
MCF-21L, DK-8000D,DK-8000L			U410-120	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-120	Adult soft 成人软指套/3.0m	C8			
			U810-120	Multi-site Y型多功能/3.0m	C7			
<b>Anke 安科(China)</b>								
ASC-553			U410-37	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-37	Adult soft 成人软指套/3.0m	C8			
			U810-37	Multi-site Y型多功能/3.0m	C7			
ASC-553 (*Oximax tech)			U410-78P	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-78P	Adult soft 成人软指套/3.0m	D1			
			U810-78P	Multi-site Y型多功能/3.0m	C8			
ASC-553A			U410-82	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-82	Adult soft 成人软指套/3.0m	C8			
			U810-82	Multi-site Y型多功能/3.0m	C7			
ASC-553A (*Ohmeda Tech)			U410-107	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-107	Adult soft 成人软指套/3.0m	C8			
			U810-107	Multi-site Y型多功能/3.0m	C7			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Artema/S&amp;W</b>								
Diascope,8000, Athena (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-80		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-80	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-80	Adult soft 成人软指套/3.0m	D1			
		U810-80	Multi-site Y型多功能/3.0m	C8				
S&W			U410-19	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-19	Adult soft 成人软指套/3.0m	D1			
			U810-19	Multi-site Y型多功能/3.0m	C8			
<b>Baxter/Simed</b>								
ASAT, S-50/100/100a			U410-12	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-12	Adult soft 成人软指套/3.0m	D1			
			U810-12	Multi-site Y型多功能/3.0m	C8			
Monarch/Explorer (*Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
<b>BCI</b>								
3300,3301,3302, 3303,3304,3401 (*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	Depends U708-06 可选U708-06		
			U103-06	Pediatric finger 儿童指夹/1.1m	C5			
			U403S-06	Adult soft 成人软指套/1.1m	C5			
			U103S-06	Pediatric soft 儿童软指套/1.1m	C5			
			U203S-06	Infant soft 婴儿软指套/1.1m	C5			
			U303-06	Neonate wrap 新生儿捆绑/1.1m	C5			
		3043	U803-06	Multi-site Y型多功能/1.1m	C3			
			U603-06	Adult/Neo Wrap 成人/新生儿捆绑/1.1m	C3			
			U603B-06	Pedia/Infant wrap 儿童/婴儿捆绑/1.1m	C3			
		3078	U903-06	Adult ear clip 成人耳夹/1.1m	C5			
3300,3301 6100,9100 Advisor, Autocorr, MiniTorr, Mini Torr plus (*BCI*)		See above*	See the above sensor*	1.1m	See above sensors*	U708-76		
			U410-76	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-76	Adult soft 成人软指套/3.0m	C8			
			U810-76	Multi-site Y型多功能/3.0m	C7			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM	REF# UNIMED	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Biosys</b>								
BPM200 BPM700 Guardian		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-27	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-27	Adult soft 成人软指套/3.0m	D1			
			U810-27	Multi-site Y型多功能/3.0m	C8			
<b>Biolight 宝莱特 (China)</b>								
M6, M12, M69			U410-47	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-47	Adult soft 成人软指套/3.0m	C8			
			U810-47	Multi-site Y型多功能/3.0m	C7			
M66, M69			U410-57	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-57	Adult soft 成人软指套/3.0m	D1			
			U810-57	Multi-site Y型多功能/3.0m	C8			
(*Oximax)			U410-51P	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-51P	Adult soft 成人软指套/3.0m	D1			
			U810-51P	Multi-site Y型多功能/3.0m	C8			
M700 (*Digital Tech)			U403-62D	Adult finger 成人指夹/1.1m	C8	No need 无需转接线		
			U403S-62D	Adult soft 成人软指套/1.1m	D1			
			U803-62D	Multi-site Y型多功能/1.1m	C8			
M8000, M9000, M8000A M9000A, BTD-352A, new M69 (*Digital Tech)			U410-62D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U410S-62D	Adult soft 成人软指套/3.0m	D2			
			U810-62D	Multi-site Y型多功能/3.0m	D1			
M900, M7000, M9500 A8 Modular (*Digital Tech)			U410-68D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U410S-68D	Adult soft 成人软指套/3.0m	D2			
			U810-68D	Multi-site Y型多功能/3.0m	D1			
<b>Bionet</b>								
		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-58		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
BM3/BM5			U410-58	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-58	Adult soft 成人软指套/3.0m	C8			
			U810-58	Multi-site Y型多功能/3.0m	C7			
BM3/BM5(New)			U410-108	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-108	Adult soft 成人软指套/3.0m	C8			
			U810-108	Multi-site Y型多功能/3.0m	C7			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM	REF# UNIMED	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
<b>Biocare (China)</b>									
BP-12B, BP-12C (*Digital Tech)			U410-61D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			U410S-61D	Adult soft 成人软指套/3.0m	D2				
			U810-61D	Multi-site Y型多功能/3.0m	D1				
BM-9000, BM-9000S FM-801, FM-801V (*Digital Tech)			U410-62D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			U410S-62D	Adult soft 成人软指套/3.0m	D2				
			U810-62D	Multi-site Y型多功能/3.0m	D1				
BM-9000			U410-29	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U410S-29	Adult soft 成人软指套/3.0m	C8				
			U810-29	Multi-site Y型多功能/3.0m	C7				
<b>CAS Medical</b>									
Oscillomate 9002 9303, 9302S (*Nonin)		8000AA	U403-08	Adult finger 成人指夹/1.1m	C3	U708-01			
			U403S-08	Adult soft 成人软指套/1.1m	C5				
			U803-08	Multi-site Y型多功能/1.1m	C3				
740-750(M) (*M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15			
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4				
			LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1			
			LNOP DC-12	*U810-15	Multi-site Y型多功能/3.0m	C8			
<b>Charmcare</b>									
CX100 (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
<b>ChoiceMed 超思 (China)</b>									
MMED3001		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3				
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
MMED100H, MMED300A, MMED6000DP (*Digital tech)			U410-84D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			U410S-84D	Adult soft 成人软指套/3.0m	D2				
			U810-84D	Multi-site Y型多功能/3.0m	D1				
MMED6000DP (*BCI Tech)			U410-96	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U410S-96	Adult soft 成人软指套/3.0m	C8				
			U810-96	Multi-site Y型多功能/3.0m	C7				
<b>Contec 康泰 (China)</b>									
			U410-29	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U410S-29	Adult soft 成人软指套/3.0m	C8				
			U810-29	Multi-site Y型多功能/3.0m	C7				
CMS60C/CW/D/DW, CMS60, PM60, PM60D, CSM5000, CMS5000C, CMS5100, CMS6000, CMS7000/ 8000/9000(*Digital Tech)			U410-46D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			U410S-46D	Adult soft 成人软指套/3.0m	D2				
			U810-46D	Multi-site Y型多功能/3.0m	D1				

Contec 康泰 (China) 下一页续接

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Contec康泰(China)</b>								
CMS50F/CMS50FW			U401-149	Adult finger 成人指夹/0.3m	C7	No need 无需转接线		
			U401S-149	Adult soft 成人软指套/0.3m	C8			
<b>Colin/Omron科林 (欧姆龙)</b>								
BP306,BP88,BP508 (Nellcor Non-Oximax) M570		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5	U708-01		
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
BP88S (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-69		
			U410-69	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-69	Adult soft 成人软指套/3.0m	D1			
			U810-69	Multi-site Y型多功能/3.0m	C8			
BX 10NE (* Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
BX 10MA (* M tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS Y1	*U803-49	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
		LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1			
			*U810-15	Multi-site Y型多功能/3.0m	C8			
<b>Corpuls</b>								
08/16 (*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	NA		
			U403S-06	Adult soft 成人软指套/1.1m	C5			
		3043	U803-06	Multi-site Y型多功能/1.1m	C3			
08/16(M) (* Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	NA		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
08/16(M) M Tech			U410-152M	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-152M	Adult soft 成人软指套/3.0m	D1			
			U810-152M	Multi-site Y型多功能/3.0m	C8			
<b>Conmen科曼 (China)</b>								
Start-8000A/8000B/8000C Star-8000D,8000H Star-5000B/5000C,C50 (*Digital Tech)			U410-66D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U410S-66D	Adult soft 成人软指套/3.0m	D2			
			U810-66D	Multi-site Y型多功能/3.0m	D1			
Start-8000E			U410-183	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-183	Adult soft 成人软指套/3.0m	C8			
			U810-183	Multi-site Y型多功能/3.0m	C7			
C60,C80 (*Oximax Tech)			U410-187P	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			U410S-187P	Adult soft 成人软指套/3.0m	D2			
			U810-187P	Multi-site Y型多功能/3.0m	D1			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Creative科瑞康(China)</b>								
UP-6000,PC-9000 (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5	U708-29		
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-29	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-29	Adult soft 成人软指套/3.0m	C8			
			U810-29	Multi-site Y型多功能/3.0m	C7			
UP-7000,UP-9000,PC-900, PC-900A,PC-900S PC-900SNPC-3000, PC-5000A/B,PC-9000B, PC-66A,PC-66B, PC-66V(*Digital tech)			U410-61D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-61D	Adult soft 成人软指套/3.0m	D1			
			U810-61D	Multi-site Y型多功能/3.0m	C8			
<b>Criticare (CSI)</b>								
510,503,504, 507,508,581, 1100,8100,Poet, Scholar (CSI*)		934-10LN	U410-05	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-05	Pediatric finger 儿童指夹/3.0m	D1			
		975AD-10	U410S-05	Adult soft 成人软指套/3.0m	D1			
			U110S-05	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-05	Infant soft 婴儿软指套/3.0m	D1			
			U310-05	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-05	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
			U610B-05	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8			
			U810-05	Multi-site Y型多功能/3.0m	C8			
			U910-05	Adult ear clip 成人耳夹/3.0m	D1			
			U403-05	Adult finger 成人指夹/1.1m	C3			
			U103-05	Pediatric finger 儿童指夹/1.1m	C4			
			U403S-05	Adult soft 成人软指套/1.1m	C4			
			U103S-05	Pediatric soft 儿童软指套/1.1m	C4			
			U203S-05	Infant soft 婴儿软指套/1.1m	C4	U708-05		
			U303-75	Neonate wrap 新生儿捆绑/1.1m	C4			
			U603-05	Multi-site Y型多功能/1.1m	C3			
			U603B-05	Adult/Neo成人/新生儿捆绑/1.1m	C3			
			U803-05	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C4			
			U903-05	Adult ear clip 成人耳夹/1.1m	C4			
503SPOT,504Plus, 504DX,506DX (*CSI)		934-10DN	U410-75	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U110-75	Pediatric finger 儿童指夹/3.0m	C8			
			U410S-75	Adult soft 成人软指套/3.0m	C8			
			U110S-75	Pediatric soft 儿童软指套/3.0m	C8			
			U210S-75	Infant soft 婴儿软指套/3.0m	C8			
			U310-75	Neonate wrap 新生儿捆绑/3.0m	C8			
			U610-75	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C7			
			U610B-75	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C7			
			U810-75	Multi-site Y型多功能/3.0m	C7			
			U910-75	Adult ear clip 成人耳夹/3.0m	C8			



OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
<b>Datascope</b>									
Accusat, Accutor 3/4 sat Passport, Accutor plus (*Datascope)			<b>U403-07</b>	Adult finger 成人指夹/1.1m	C4	U708-07			
			<b>U103-07</b>	Pediatric finger 儿童指夹/1.1m	C6				
			<b>U403S-07</b>	Adult soft 成人软指套/1.1m	C6				
			<b>U103S-07</b>	Pediatric soft 儿童软指套/1.1m	C6				
			<b>U203S-07</b>	Infant soft 婴儿软指套/1.1m	C6				
			<b>U303-07</b>	Neonate wrap 新生儿捆绑/1.1m	C6				
			<b>U603-07</b>	Adult/Neo wrap 成人/新生儿捆绑/1.1m	C5				
			<b>U603B-07</b>	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C5				
			<b>U803-07</b>	Multi-site Y型多功能/1.1m	C5				
			<b>U903-07</b>	Adult ear clip 成人耳夹/1.1m	C4				
Accutor plus (*M Tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-15			
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4				
		LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8				No need 无需转接线
		LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1				
			<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8				
PassportII, Exper,t Pssport XG, Accutor plus (*Non-Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-69			
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5				
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3				
Accutor-Trio (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-01			
Trio Nouvelle Generation (* Oximax)		DS100A	<b>U403-01X</b>	Adult finger 成人指夹/1.1m	C5	U710-70			
		DS100A	<b>U403S-01X</b>	Adult soft 成人软指套/1.1m	C7				
		D-YS	<b>U803-01X</b>	Multi-site Y型多功能/1.1m	C5				
<b>Digicare</b>									
Digimax5500 (*BCI)		3044	<b>U403-06</b>	Adult finger 成人指夹/1.1m	C3	U708-76			
			<b>U403S-06</b>	Adult soft 成人软指套/1.1m	C5				
		3043	<b>U803-06</b>	Multi-site Y型多功能/1.1m	C3				
			<b>U410-76</b>	Adult finger 成人指夹/3.0m	C7				No need 无需转接线
			<b>U410S-76</b>	Adult soft 成人软指套/3.0m	C8				
			<b>U810-76</b>	Multi-site Y型多功能/3.0m	C7				
		<b>U410-52</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线				
		<b>U410S-52</b>	Adult soft 成人软指套/3.0m	C8					
		<b>U810-52</b>	Multi-site Y型多功能/3.0m	C7					

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
<b>Dima (Italia)</b>									
(*Non-Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-01			
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5				
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3				
<b>Dixtal(Brazil)-Philips</b>									
Full compatible with Novamatrix Sensas		DX-1NC01-0	<b>U410-113</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			<b>U410S-113</b>	Adult soft 成人软指套/3.0m	C8				
			<b>U810-113</b>	Multi-site Y型多功能/3.0m	C7				
<b>Dolphin</b>									
(*Digital Tech)			<b>U403-50D</b>	Adult finger 成人指夹/1.1m	C8	U708-50			
			<b>U403S-50D</b>	Adult soft 成人软指套/1.1m	D1				
			<b>U803-50D</b>	Multi-site Y型多功能/1.1m	C8				
(*Digital Tech)			<b>U410-50D</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			<b>U410S-50D</b>	Adult soft 成人软指套/3.0m	D2				
			<b>U810-50D</b>	Multi-site Y型多功能/3.0m	D1				
<b>Drager/Siemens德格/西门子</b>									
Oxipac 2000 (*BCI)		3044	<b>U403-06</b>	Adult finger 成人指夹/1.1m	C3	No need 无需转接线			
			<b>U403S-06</b>	Adult soft 成人软指套/1.1m	C5				
		3043	<b>U803-06</b>	Multi-site Y型多功能/1.1m	C3				
Vitalert 1000,icero (*Non-Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-28			
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5				
Narkomed-Oxysat 2 (*Non-Oximax)		See the above	See the above	See the above	See the above	U710-85			
Sirecust, Sirem Spo2 (*Non-Oximax)		See the above	See the above	See the above	See the above	U708-33			
SC6000, 7000,9000 (*Non-Oximax)		See the above	See the above	See the above	See the above	U704-23 U708-23			
Oxytrend (*M tech))		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-15			
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS Y1	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4				
		LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8				No need 无需转接线
		LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1				
			<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8				

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Edan理邦 (China)</b>								
M9,M99A,M8,M8A, M8B, M3,M3B (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U708X-191		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
EDAN Oximetry M3,M3A,M3B, M8,M8A,M8B M9,M9A,M50,M80 iM8,iM8A,iM8B iM9,iM9A,iM50,iM80 (*Non-Oximax)		MS3-109069	U410-29	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			U410S-29	Adult soft 成人软指套/3.0m	C8			
			U810-29	Multi-site Y型多功能/3.0m	C7			
		MS-109079	U403-01	Adult finger 成人指夹/1.1m	C3			
M3,M3B,M8,M8A,M8B M9,M9A,M50,M80 iM9,iM9A,iM50,iM80 (*Oximax)		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5	U708-29		
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
			U410-191P	Adult finger 成人指夹/3.0m	C8			
(*Digital tech)			U410S-191P	Adult soft 成人软指套/3.0m	D1	No need 无需转接线		
			U810-191P	Multi-site Y型多功能/3.0m	C8			
			U410-46D	Adult finger 成人指夹/3.0m	C8			
H100B (*Non Oximax)		MS-109079	U403-01	Adult finger 成人指夹/1.1m	C3	No need 无需转接线		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
			U410S-46D	Adult soft 成人软指套/3.0m	D1			
H100N (*Oximax)			U410S-46D	Adult soft 成人软指套/3.0m	D1	No need 无需转接线		
			U810-46D	Multi-site Y型多功能/3.0m	C8			
			U410-46D	Adult finger 成人指夹/3.0m	C8			
<b>Fukuda Denshi福田</b>								
DS7100-DS7300 (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
<b>GE Corometrics</b>								
188,129,129F, 128,128F (*Non Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3			
120,510,511, 118,556 (*Non Oximax)		See the above	See the above	See the above	See the above	U708-25 (4033CAX)		
128, 129 (*Non Oximax)		See the above	See the above	See the above	See the above	U708-21		
SERIE 250 (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710X-21		
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7			
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5			
		See Oximax Marquette	See Oximax Marquette	See Oximax Marquette	See GE	No need 无需转接线		

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
<b>GE Critikon</b>									
Dianmap 9700/9701 9702 PLUS, MPS ORTABLE,/COMPACT, (* Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-01 (EC-8)			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
Dianmap,ProSeries V1 Pro 1000 V2 (Nellcor Optional)		See the above	See the above	See the above	See the above	U708-69 (SCP-10)			
Dianmap,ProCare, ProSeries V2Pro1000 V3 (*Oximax)		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	U710-70 (DOC-10)			
		DS100A	U403S-01X	Adult soft 成人软指套/1.1m	C7				
		D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5				
Dianmap 8700/8710/8720Pulse, Oxyshuttle11 et 2		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	Not Available 不能提供			
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5				
		D-YS	U803-01	Multi-site Y型多功能/1.1m	C3				
				U410-10	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
				U410S-10	Adult soft 成人软指套/3.0m	D1			
				U810-10	Multi-site Y型多功能/3.0m	C8			
Dianmap ProCare Series V2 Pro 1000 V3 (*M Tech)		1863/ LNCS DCI	*U403-49	Adult finger 成人指夹/1.1m	C4	U708M-15 (2027263-001)			
		1863/ LNCS DCI	*U403S-49	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS YI	*U803-49	Multi-site Y型多功能/1.1m	C4				
	LNOP DC-12	*U410-15	Adult finger 成人指夹/3.0m	C8					
	LNOP DC-12	*U410S-15	Adult soft 成人软指套/3.0m	D1					
			*U810-15	Multi-site Y型多功能/3.0m	C8	No need 无需转接线			
<b>GE Datex-Ohmeda</b>									
AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPI-II Capnomac Ultima Light monitor Satlite, Oscar (Datex tech*)		OXY-F-DB	U403-09	Adult finger 成人指夹/1.1m	C3	U708-09			
			U103-09	Pediatric finger 儿童指夹/1.1m	C5				
			U403S-09	Adult soft 成人软指套/1.1m	C5				
			U103S-09	Pediatric soft 儿童软指套/1.1m	C5				
			U203S-09	Infant soft 婴儿软指套/1.1m	C5				
			U303-09	Neonate wrap 新生儿捆绑/1.1m	C5				
		OXY-W-DB	U603-09	Adult/Neo wrap成人/新生儿捆绑/1.1m	C3				
			U603B-09	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C3				
			U803-09	Multi-site Y型多功能/1.1m	C3				
		OXY-E-DB	U903-09	Adult ear clip 成人耳夹/1.1m	C5				

GE Datex-Ohmeda 下一页续接

GE Datex-Ohmeda 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>GE Datex-Ohmeda</b>								
AS/3,CS/3, AS/5 (Built in M module) CH/S,CH/RSCG/2S, CARDIOCAPI-II Capnomac Ultima Light monitor Satlite, Oscar (Datex tech*)		OXY-F4-N	<b>U410-09</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U110-09</b>	Pediatric finger 儿童指夹/3.0m	D1			
			<b>U410S-09</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U110S-09</b>	Pediatric soft 儿童软指套/3.0m	D1			
			<b>U210S-09</b>	Infant soft 婴儿软指套/3.0m	D1			
			<b>U310-09</b>	Neonate wrap 新生儿捆绑/3.0m	D1			
		OXY-W4-N	<b>U610-09</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
			<b>U610B-09</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8			
		OXY-W4-N	<b>U810-09</b>	Multi-site Y型多功能/3.0m	C8			
		OXY-E4-N	<b>U910-09</b>	Adult ear clip 成人耳夹/3.0m	D1			
M-NSAT Module (AS/3,CS/3,S/5) (* Non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-69 (MC-10)		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
M-NSAT Module (i4) (*Oximax)		DS100A	<b>U403-01X</b>	Adult finger 成人指夹/1.1m	C5	U710-70 (DOC-10)		
		DS100A	<b>U403S-01X</b>	Adult soft 成人软指套/1.1m	C7			
		D-YS	<b>U803-01X</b>	Multi-site Y型多功能/1.1m	C5			
3700,3710,3770, 3775,Biox3740,3760, 3800,3900,3900P, 4700,5250RGB, Cardiocap/5, S/5, AS/3 (built in Module N-XOSAT) CS/3 (Built in module M-OSAT) (*Ohmeda)		OXY-F4-H	<b>U410-02</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U110-02</b>	Pediatric finger 儿童指夹/3.0m	D1			
			<b>U410S-02</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U110S-02</b>	Pediatric soft 儿童软指套/3.0m	D1			
			<b>U210S-02</b>	Infant soft 婴儿软指套/3.0m	D1			
			<b>U310-02</b>	Neonate wrap 新生儿捆绑/3.0m	D1			
		OXY-W4-H	<b>U610-02</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
			<b>U610B-02</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8			
		OXY-W4-H	<b>U810-02</b>	Multi-site Y型多功能/3.0m	C8			
		OXY-E4-H	<b>U910-02</b>	Adult ear clip 成人耳夹/3.0m	D1			
Ohmeda Tuffsat (*Ohmeda)		OXY-F1-H	<b>U403-02</b>	Adult finger 成人指夹/1.1m	C8	No need 无需转接线		
		OXY-F1-H	<b>U403S-02</b>	Adult soft 成人软指套/1.1m	D1			
Ohmeda TruSat Ohmeda TruSat TD (*+Oxytip+ Technology)		OXY-F4-MC	<b>U410-60</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			<b>U110-60</b>	Pediatric finger 儿童指夹/3.0m	D2			
		OXY-F4-MC	<b>U410S-60</b>	Adult soft 成人软指套/3.0m	D2			
			<b>U110S-60</b>	Pediatric soft 儿童软指套/3.0m	D2			
			<b>U210S-60</b>	Infant soft 婴儿软指套/3.0m	D2			
			<b>U310-60</b>	Neonate wrap 新生儿捆绑/3.0m	D2			
		OXY-W4-MC	<b>U610-60</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	D1			
			<b>U610B-60</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	D1			
			<b>U810-60</b>	Multi-site Y型多功能/3.0m	D1			
		OXY-E4-MC	<b>U910-60</b>	Adult ear clip 成人耳夹/3.0m	D2			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>GE Hellige海力格</b>								
Cardioserv Eagle1000-3000 (*Non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-55 (30344358)		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
Cardioserv (* Oximax)		DS100A	<b>U403-01X</b>	Adult finger 成人指夹/1.1m	C5	NA (2026751-001)		
		DS100A	<b>U403S-01X</b>	Adult soft 成人软指套/1.1m	C7			
		D-YS	<b>U803-01X</b>	Multi-site Y型多功能/1.1m	C5			
SMU,SMK,SMV (*Ohmeda)		See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	No need 无需转接线		
<b>GE Marquette</b>								
Eagle 1000/3000/4000 Dash 2000/3000/4000 Solar 8000/8000M i9500(Tram 451/451N /851/851N) CardioLab/Mac-Lab, Corometrics 120/2120 (Tran 451/451N) Responder3000 (nellcor optional) (*Non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-21 (2006644-001)		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
			<b>U410-21</b>	Adult finger 成人指夹/3.0m	C8			
			<b>U110-21</b>	Pediatric finger 儿童指夹/3.0m	D1			
			<b>U410S-21</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U110S-21</b>	Pediatric soft 儿童软指套/3.0m	D1			
			<b>U210S-21</b>	Infant soft 婴儿软指套/3.0m	D1			
			<b>U310-21</b>	Neonate wrap 新生儿捆绑/3.0m	D1			
			<b>U610-21</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
	<b>U610B-21</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8					
	<b>U810-21</b>	Multi-site Y型多功能/3.0m	C8					
	<b>U910-21</b>	Adult ear clip 成人耳夹/3.0m	D1					
Dash2500/Dash 3000/ 4000/5000 Solar 8000/8000M i9500 (Tram 451N5/851N5) Corometrics 250 (GE-Oximax*) Ending " X" can be replaced by " P" for fully compatible with all Oximax oximetry		DS100A	<b>U403-01X</b>	Adult finger 成人指夹/1.1m	C5	U710X-21 (2021406-001)		
		DS100A	<b>U403S-01X</b>	Adult soft 成人软指套/1.1m	C7			
		D-YS	<b>U803-01X</b>	Multi-site Y型多功能/1.1m	C5			
			<b>U410-21P</b>	Adult finger 成人指夹/3.0m	D1			
			<b>U110-21P</b>	Pediatric finger 儿童指夹/3.0m	D2			
			<b>U410S-21P</b>	Adult soft 成人软指套/3.0m	D2			
			<b>U110S-21P</b>	Pediatric soft 儿童软指套/3.0m	D2			
			<b>U210S-21P</b>	Infant soft 婴儿软指套/3.0m	D2			
			<b>U310-21P</b>	Neonate wrap 新生儿捆绑/3.0m	D2			
			<b>U610-21P</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	D1			
	<b>U610B-21P</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	D1					
	<b>U810-21P</b>	Multi-site Y型多功能/3.0m	D1					
	<b>U910-21P</b>	Adult ear clip 成人耳夹/3.0m	D2					
7000/7010,Tram 200/400/600/800 (*Ohmeda Tech.)		See above *Ohmeda	See above *Ohmeda	See above *Ohmeda	See above *Ohmeda			

GE Marquette 下页续接

GE Marquette 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>GE Marquette</b>								
Dash 3000, 4000,5000 Solar 8000/8000M i/9500 (Tram 451M,851M) Cardiolab/Mac-Lab (Tram 451M) Corometrics 259 (*M tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-21 (2027263-002)		
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4			
			<b>*U410-21M</b>	Adult finger 成人指夹/3.0m	C8			
			<b>*U110-21M</b>	Pediatric finger 儿童指夹/3.0m	D1			
			<b>*U410S-21M</b>	Adult soft 成人软指套/3.0m	D1			
			<b>*U110S-21M</b>	Pediatric soft 儿童软指套/3.0m	D1			
			<b>*U210S-21M</b>	Infant soft 婴儿软指套/3.0m	D1			
			<b>*U310-21M</b>	Neonate wrap 新生儿捆绑/3.0m	D1			
			<b>*U610-21M</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
	<b>*U610B-21M</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8					
	<b>*U810-21M</b>	Multi-site Y型多功能/3.0m	C8					
	<b>*U910-21M</b>	Adult ear clip 成人耳夹/3.0m	D1					
<b>GE Ohmeda OxyTip+</b>								
GE S/5 (i4 E series module built in)		OXY-F4-GE	<b>U410-121</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			<b>U110-121</b>	Pediatric finger 儿童指夹/3.0m	D2			
			<b>U410S-121</b>	Adult soft 成人软指套/3.0m	D2			
			<b>U110S-121</b>	Pediatric soft 儿童软指套/3.0m	D2			
			<b>U210S-121</b>	Infant soft 婴儿软指套/3.0m	D2			
			<b>U310-121</b>	Neonate wrap 新生儿捆绑/3.0m	D2			
			<b>U610-121</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	D1			
			<b>U610B-121</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	D1			
			<b>U810-121</b>	Multi-site Y型多功能/3.0m	D1			
	<b>U910-121</b>	Adult ear clip 成人耳夹/3.0m	D2					
Ohmeda TruSat 3755,37XX series 3800,3900,3900P 4700, 5250 Cardiocap/5, S/5,AS/3 (N-XOSAT or M built in) CS/3 (M-OSAT built in) S/5 light, Satlite Puls II Capnomac UltmaS/5 ,Coro 259 (*TruSignal tech)		OXY-F-UN	<b>U403-127</b>	Adult finger 成人指夹/1.1m	C8	Cables required depending on the monitor models 必须配用转接线 取决于机器端接头	Depends on the monitor models See the page of Spo2 adapter cable 取决于机器端接头 参考后面血氧转接线	
			<b>U103-127</b>	Pediatric finger 儿童指夹/1.1m	D1			
		OXY-F-UN	<b>U403S-127</b>	Adult soft 成人软指套/1.1m	D1			
			<b>U103S-127</b>	Pediatric soft 儿童软指套/1.1m	D1			
			<b>U203S-127</b>	Infant soft 婴儿软指套/1.1m	D1			
			<b>U303-127</b>	Neonate wrap 新生儿捆绑/1.1m	D1			
		OXY-W-UN	<b>U603-127</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C8			
		OXY-W-UN	<b>U603B-127</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C8			
		OXY-W-UN	<b>U803-127</b>	Multi-site Y型多功能/1.1m	C8			
		OXY-E-UN	<b>U903-127</b>	Adult ear clip 成人耳夹/1.1m	D1			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>GE Trusignal</b>								
		TS-F-D	<b>U403-117</b>	Adult finger 成人指夹/1.1m	C8	Required 必须配用转接线	Depends on the monitor models See the page of Spo2 adapter cable 取决于机器端接头 参考后面血氧转接线	
			<b>U103-117</b>	Pediatric finger 儿童指夹/1.1m	D1			
		TS-SA-D	<b>U403S-117</b>	Adult soft 成人软指套/1.1m	D1			
			<b>U103S-117</b>	Pediatric soft 儿童软指套/1.1m	D1			
			<b>U203S-117</b>	Infant soft 婴儿软指套/1.1m	D1			
			<b>U303-117</b>	Neonate wrap 新生儿捆绑/1.1m	D1			
		TS-W-D	<b>U603-117</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C8			
		TS-W-D	<b>U603B-117</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C8			
		TS-W-D	<b>U803-117</b>	Multi-site Y型多功能/1.1m	C8			
		TS-E-D	<b>U903-117</b>	Adult ear clip 成人耳夹/1.1m	D1			
<b>General 杰纳瑞(China)</b>								
G1B (*Digital Tech)			<b>U403-42D</b>	Adult finger 成人指夹/1.1m	C5	No need 无需转接线		
			<b>U403S-42D</b>	Adult soft 成人软指套/1.1m	C7			
G2A,G3A,G3C,G3D, G3F,G3G,G3H,G3L, G6A,G6A plus, G9L (*Digital Tech)			<b>U410-65D</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-65D</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-65D</b>	Multi-site Y型多功能/3.0m	C8			
<b>Goldway金科威(China)</b>								
UT4000A,UT4000B, UT4000F UT4000ACS			<b>U410-37</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			<b>U410S-37</b>	Adult soft 成人软指套/3.0m	C8			
			<b>U810-37</b>	Multi-site Y型多功能/3.0m	C7			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-37		
		See Pace tech	See Pace tech	See Pace tech	See Pace tech	No need 无需转接线		
UT4000A,UT4000C, UT4000EUT4000F, UT4000APro, UT4000FPro, UT6000A			<b>U410-11</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-11</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-11</b>	Multi-site Y型多功能/3.0m	C8			
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-11		
		See Pace tech	See Pace tech	See Pace tech	See Pace tech	No need 无需转接线		
3100		Digital sensor 数字探头	<b>U410-61D</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-61D</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-67D</b>	Multi-site Y型多功能/3.0m	C8			
UT4000B,UT4000F			<b>U410-67D</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-67D</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-67D</b>	Multi-site Y型多功能/3.0m	C8			

Goldway金科威(China) 下页续接

Goldway金科威(China) 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头						
<b>Goldway金科威(China)</b>														
(*Oximax)			<b>U410-38P</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线								
			<b>U410S-38P</b>	Adult soft 成人软指套/3.0m	D2									
			<b>U810-38P</b>	Multi-site Y型多功能/3.0m	D1									
UT6000A G3 (*Philips tech)			<b>U410-91</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线								
			<b>U410S-91</b>	Adult soft 成人软指套/3.0m	D1									
			<b>U810-91</b>	Multi-site Y型多功能/3.0m	C8									
<b>Guoteng 国腾(China)</b>														
PM300		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	No need 无需转接线								
			<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5									
			<b>U410-171</b>	Adult finger 成人指夹/3.0m	C8									
GT3000,GT5000, GT6000GT7000, GT9000,GT9002			<b>U410S-171</b>	Adult soft 成人软指套/3.0m	D1	No need 无需转接线								
			<b>U410-179</b>	Adult finger 成人指夹/3.0m	C7									
			<b>U410S-179</b>	Adult soft 成人软指套/3.0m	C8									
			<b>U810-179</b>	Multi-site Y型多功能/3.0m	C7									
<b>IVY</b>														
405A,405D (* Non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U710-85								
			<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5									
			<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3									
405P,405T (*Non Oximax)		See the above	See the above	See the above	See the above	U708-01								
Vital-Guard 450CN(M) (*M tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-15								
			<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6									
			<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4									
		2258/ LNCS Y1	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8									
			<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1									
			<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8									
<b>Kontron康强</b>														
7840,7845,Minimon 7337Micromon, Supermon,Kolormon, Trakmon (*Non Oximax)		DS100A	<b>U410-13</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线								
			<b>U410S-13</b>	Adult soft 成人软指套/3.0m	D1									
			<b>U810-13</b>	Multi-site Y型多功能/3.0m	C8									
		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3				U708-13					
			<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5									
			<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3									
7150,7850, Ivy 2000(M)		See the above	See the above	See the above	See the above	U708-01								
			See the above	See the above	See the above	U708-53								










OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头				
<b>Lohmeier隆美尔</b>												
			<b>U410-64</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线						
			<b>U410S-64</b>	Adult soft 成人软指套/3.0m	D2							
			<b>U810-64</b>	Multi-site Y型多功能/3.0m	D1							
<b>*M Technology</b>												
Rad-8 Rad-5/Rad-5v SatShare			1269/ LNOP DCI	<b>*U403-125</b>	Adult finger 成人指夹/1.1m	C7	Required 必须配用转接线					
			1276/ LNOP DCIP	<b>*U103-125</b>	Pediatric finger 儿童指夹/1.1m	C8						
			1269/ LNOP DCI	<b>*U403S-125</b>	Adult soft 成人软指套/1.1m	C8						
			1276/ LNOP DCIP	<b>*U103S-125</b>	Pediatric soft 儿童软指套/1.1m	C8						
				<b>*U203S-125</b>	Infant soft 婴儿软指套/1.1m	C8						
				<b>*U303-125</b>	Neonate wrap 新生儿捆绑/1.1m	C8						
				<b>*U603-125</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C7						
				<b>*U603B-125</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C7						
			1504/ LNOP Y1	<b>*U803-125</b>	Multi-site Y型多功能/1.1m	C7						
			LNOP TC-I	<b>*U903-125</b>	Adult ear clip 成人耳夹/1.1m	C8						
			1863/ LNCS DCI 1864/ LNCS DCIP 1863/ LNCS DCI			<b>*U403-49</b>			Adult finger 成人指夹/1.1m	C4	Required 必须配用转接线	
						<b>*U103-49</b>			Pediatric finger 儿童指夹/1.1m	C6		
						<b>*U403S-49</b>			Adult soft 成人软指套/1.1m	C6		
						<b>*U103S-49</b>			Pediatric soft 儿童软指套/1.1m	C6		
						<b>*U203S-49</b>			Infant soft 婴儿软指套/1.1m	C6		
<b>*U303-49</b>	Neonate wrap 新生儿捆绑/1.1m	C6										
<b>*U603-49</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C4										
<b>*U603B-49</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C4										
2258/ LNCS Y1	<b>*U803-49</b>	Multi-site Y型多功能/1.1m				C4						
LNCS TC-I	<b>*U903-49</b>	Adult ear clip 成人耳夹/1.1m	C6									
M Tech *Refurbished Only			LNOP DCSC	<b>*U402-15</b>	Adult finger 成人指夹/0.9m	C8	No need 无需转接线	Handheld				
			LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8						
				<b>*U110-15</b>	Pediatric finger 儿童指夹/3.0m	D1						
			LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1						
				<b>*U110S-15</b>	Pediatric soft 儿童软指套/3.0m	D1						
				<b>*U210S-15</b>	Infant soft 婴儿软指套/3.0m	D1						
				<b>*U310-15</b>	Neonate wrap 新生儿捆绑/3.0m	D1						
				<b>*U610-15</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8						
				<b>*U610B-15</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8						
				<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8						
	<b>*U910-15</b>	Adult ear clip 成人耳夹/3.0m	D1									







OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>MDE/Invivo</b>								
Escort 1 Escort 300 (*Non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U710-85		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
20300,20401,20403, Prism (*Non Oximax)		See the above	See the above	See the above	See the above	U708-01		
Millenia, Escort,(M) M10TM,M12TM, M8TM (*M Tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4			
	LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8				
		LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1	No need 无需转接线		
			<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8			
4500, ,Omni Track Omeg 1445Millennia			<b>U410-14</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-14</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-14</b>	Multi-site Y型多功能/3.0m	C8			
4500, ,Omni Track Omeg 1445Millennia M Tech			<b>U410-14M</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-14M</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-14M</b>	Multi-site Y型多功能/3.0m	C8			
<b>Meditronic - Physio Control</b>								
Lifepack 12 (*non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-01		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
Lifepack 12, 20(M) (*M Tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
		LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1			
			<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8			
<b>Medlab</b>								
NANOX 10 POX 10 Perarl 10 (*BCI)		3044	<b>U403-06</b>	Adult finger 成人指夹/1.1m	C3	U708-76		
			<b>U403S-06</b>	Adult soft 成人软指套/1.1m	C5			
		3043	<b>U803-06</b>	Multi-site Y型多功能/1.1m	C3			
<b>MEK</b>								
MP 100,110, 400,1000			<b>U410-26</b>	Adult finger 成人指夹/3.0m	C6	No need 无需转接线		
			<b>U410S-26</b>	Adult soft 成人软指套/3.0m	C7			
			<b>U810-26</b>	Multi-site Y型多功能/3.0m	C6			
(*non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-56		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>MEK</b>								
			<b>U410-56</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-56</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-56</b>	Multi-site Y型多功能/3.0m	C8			
<b>Mennen 美能</b>								
Envoy (*non Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-54		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
Horizon 1100,2000XL Mercury,MR1300 (*non Oximax)		See the above	See the above	See the above	See the above	U710-85		
Horizon XL,Mercury 1200,MR1300 (*non Oximax)		See the above	See the above	See the above	See the above	U708-24		
Horizon XL (*Ohmeda)		OXY-F4-H	<b>U410-02</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-02</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-02</b>	Multi-site Y型多功能/3.0m	C8			
Horizon,Envoy, Vitalogik (M) (*M Tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-59		
			<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C6			
<b>Mindray迈瑞 ( China)</b>								
509B, MEC-1000, MEC-1200 New PM6000 PM7000/8000/9000, VS 800,PM600 PM9000 Express (*Mindray Tech)			<b>U410-29</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			<b>U110-29</b>	Pediatric finger 儿童指夹/3.0m	C8			
			<b>U410S-29</b>	Adult soft 成人软指套/3.0m	C8			
			<b>U110S-29</b>	Pediatric soft 儿童软指套/3.0m	C8			
			<b>U210S-29</b>	Infant soft 婴儿软指套/3.0m	C8			
			<b>U310-29</b>	Neonate wrap 新生儿捆绑/3.0m	C8			
			<b>U610-29</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	C7			
			<b>U610B-29</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C7			
			<b>U810-29</b>	Multi-site Y型多功能/3.0m	C7			
			<b>U910-29</b>	Adult ear clip成人耳夹/3.0m	C8			
		512D	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-29 (0010-20-42594)		
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5			
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3			
PM-5000,PM 6000 PM-8000 (*Mindray Tech)		See the above	See the above	See the above	See the above	U708-22 (0515-30-11221)		
			<b>U410-22</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-22</b>	Adult soft 成人软指套/3.0m	D1			
	<b>U810-22</b>	Multi-site Y型多功能/3.0m	C8					

Mindray迈瑞 ( China) 下一页续接














Mindray迈瑞 (China) 续接上页










OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头		
Mindray迈瑞 (China)										
PM-6000 (*M tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-22 (0010-30-12452)				
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6					
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4					
		<b>*U410-22M</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
		<b>*U410S-22M</b>	Adult soft 成人软指套/3.0m	D1						
New PM-6000, PM-8000 (*M tech)		<b>*U810-22M</b>	Multi-site Y型多功能/3.0m	C8	No need 无需转接线					
		<b>*U410-29M</b>	Adult finger 成人指夹/3.0m	C7						
		<b>*U410S-29M</b>	Adult soft 成人软指套/3.0m	C8						
New PM-6000 PM7000,8000,9000 PM-9000 Express VS 800 New PM6000 (*M tech)		See the above M sensor	See the above M sensor	See the above M sensor	See the above	U708M-30				
		<b>*U410-30M</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线					
		<b>*U410S-30M</b>	Adult soft 成人软指套/3.0m	C8						
New PM-6000 PM7000,8000,9000 VS 800 PM-9000 Express (*Oximax)		DS100A	<b>U403-01X</b>	Adult finger 成人指夹/1.1m	C5	U710X-180				
		DS100A	<b>U403S-01X</b>	Adult soft 成人软指套/1.1m	C7					
		D-YS	<b>U803-01X</b>	Multi-site Y型多功能/1.1m	C5					
		<b>U410-180P</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
<b>U410S-180P</b>	Adult soft 成人软指套/3.0m	D1								
Beneview T5/T8 (*Mindray Tech)		See the above Non Oximax Sensor	See the above Non Oximax Sensor	1.1m	See the above	U708-48				
		<b>U410-48</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
		<b>U410S-48</b>	Adult soft 成人软指套/3.0m	D1						
	<b>U810-48</b>	Multi-site Y型多功能/3.0m	C8							
	<b>U810-48</b>	Multi-site Y型多功能/3.0m	C8							
Beneview T5/T8 (*Oximax*)		See Nellcor Oximax	See Nellcor Oximax	1.1m	See Oximax	U710X-48				
Beneview T5/T8 (*M tech)		See M	See M	1.1m	See M	U708M-48				
503,507,508, 509,PM9000 (*Criticare tech)		<b>U410-05</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
		<b>U410S-05</b>	Adult soft 成人软指套/3.0m	D1						
		<b>U810-05</b>	Multi-site Y型多功能/3.0m	C8						
500A,509B, PM9000 (*BCI)		<b>U410-31</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
		<b>U410S-31</b>	Adult soft 成人软指套/3.0m	D1						
		<b>U810-31</b>	Multi-site Y型多功能/3.0m	C8						
PM-50,PM-60 (*Non-Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	No need 无需转接线				
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5					

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头			
Minolta美能达											
PULSOX-300 PULSOX-300i		SR-5C	<b>U401-105</b>	Adult finger 成人指夹/0.3m	C8	No need 无需转接线					
		SR-5C	<b>U401S-105</b>	Adult soft 成人软指套/0.3m	D1						
M&B麦邦(China)											
MB526,MB526+,MB526T8 MB526T12,MB526+tw MB800,CD2000,CD2000+			<b>U410-92</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线					
			<b>U410S-92</b>	Adult soft 成人软指套/3.0m	C8						
			<b>U810-92</b>	Multi-site Y型多功能/3.0m	C7						
CD2000,CD2000+ (* Oximax tech)			<b>U410-92P</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线					
			<b>U410S-92P</b>	Adult soft 成人软指套/3.0m	D1						
			<b>U810-92P</b>	Multi-site Y型多功能/3.0m	C8						
Nellcor-Tyco											
NPB 40, NPB20, NPB75 NPB 180, NPB 185, NPB-190, NPB-195, NPB 3900, NPB 4000 (*Non-Oximax)+		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U704-01 U708-01					
		OXI-P/I	<b>U103-01</b>	Pediatric finger 儿童指夹/1.1m	C5						
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5						
		OXI-P/I	<b>U103S-01</b>	Pediatric soft 儿童软指套/1.1m	C5						
		OXI-P/I	<b>U203S-01</b>	Infant soft 婴儿软指套/1.1m	C5						
		OXI-A/N	<b>U303-01</b>	Neonate wrap 新生儿捆绑/1.1m	C5						
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3						
		OXI-A/N	<b>U603-01</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C3						
		OXI-P/I	<b>U603B-01</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C3						
		D-YS+YSE	<b>U903-01</b>	Adult ear clip 成人耳夹/1.1m	C5						
			<b>U410-01</b>	Adult finger 成人指夹/3.0m	C5				No need 无需转接线		
			<b>U110-01</b>	Pediatric finger 儿童指夹/3.0m	C7						
			<b>U410S-01</b>	Adult soft 成人软指套/3.0m	C7						
	<b>U110S-01</b>	Pediatric soft 儿童软指套/3.0m	C7								
	<b>U210S-01</b>	Infant soft 婴儿软指套/3.0m	C7								
	<b>U310-01</b>	Neonate wrap 新生儿捆绑/3.0m	C7								
	<b>U610-01</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	C5								
	<b>U610B-01</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C5								
	<b>U810-01</b>	Multi-site Y型多功能/3.0m	C5								
	<b>U910-01</b>	Adult ear clip 成人耳夹/3.0m	C7								
N-100,N-200		See the above	See the above	1.1m	See the above	U710-85					

Nellcor-Tyco 下页续接

↳ Nellcor-Tyco 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
<b>Nellcor-Tyco</b>									
NPB3000-3100, Symphony, NPB-290, NPB-295, NPB-395		See the above	See the above	1.1m	See the above	U708-69			
			U410-69		Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-69		Pediatric finger 儿童指夹/3.0m	D1			
			U410S-69		Adult soft 成人软指套/3.0m	D1			
			U110S-69		Pediatric soft 儿童软指套/3.0m	D1			
			U210S-69		Infant soft 婴儿软指套/3.0m	D1			
			U310-69		Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-69		Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
			U610B-69		Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8			
	U810-69		Multi-site Y型多功能/3.0m	C8					
	U910-69		Adult ear clip 成人耳夹/3.0m	D1					
N-65 (*Oximax) Ending "X" can be replaced by "XP" or "P" for fully compatible with all Oximax oximetry		DS100A	U403-01X	Adult finger 成人指夹/1.1m	C5	No need 无需转接线			
		OXI-P/I	U103-01X	Pediatric finger 儿童指夹/1.1m	C7				
		OXI-A/N	U403S-01X	Adult soft 成人软指套/1.1m	C7				
		OXI-P/I	U103S-01X	Pediatric soft 儿童软指套/1.1m	C7				
		OXI-P/I	U203S-01X	Infant soft 婴儿软指套/1.1m	C7				
		OXI-A/N	U303-01X	Neonate wrap 新生儿捆绑/1.1m	C7				
		OXI-A/N	U603-01X	Adult/Neo wrap 成人/新生儿捆绑/1.1m	C5				
		OXI-P/I	U603B-01X	Pedia/infant wrap 儿童/婴儿捆绑/1.1m	C5				
	D-YS	U803-01X	Multi-site Y型多功能/1.1m	C5					
	D-YS+YSE	U903-01X	Adult ear clip 成人耳夹/1.1m	C7					
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U701-01X			
N550,N560,N595, N600,N751, N5500		See Nellcor Oximax	See Nellcor Oximax	1.1m	See Oximax	U710-70			
<b>Neusoft 东软(China)</b>									
NSC series (*Oximax Tech)			U410-78P	Adult finger 成人指夹/3.0m	D1	No need 无需转接线			
			U410S-78P	Adult soft 成人软指套/3.0m	D2				
			U810-78P	Multi-site Y型多功能/3.0m	D1				
			U410-107	Adult finger 成人指夹/3.0m	C7	No need 无需转接线			
			U410S-107	Adult soft 成人软指套/3.0m	C8				
			U810-107	Multi-site Y型多功能/3.0m	C7				

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头		
<b>New tech 纽泰克 (China)</b>										
Pulse Oximer NT1A,NT1B,NT1D (Digital Sensor)			U403-42D	Adult finger 成人指夹/1.1m	C5	No need 无需转接线				
			U403S-42D	Adult soft 成人软指套/1.1m	C7					
Patient monitor Digital sensor			U410-65D	Adult finger 成人指夹/3.0m	C8	No need 无需转接线				
			U410S-65D	Adult soft 成人软指套/3.0m	D1					
NT2A,NT2C			U810-65D	Multi-site Y型多功能/3.0m	C8	No need 无需转接线				
			U410-165D	Adult finger 成人指夹/3.0m	D1					
			U410S-165D	Adult soft 成人软指套/3.0m	D2					
			U810-165D	Multi-site Y型多功能/3.0m	D1					
NT3A,NT3B, NT3E,NT3F, NT3B60			U410-167	Adult finger 成人指夹/3.0m	C8	No need 无需转接线				
			U410S-167	Adult soft 成人软指套/3.0m	D1					
			U810-167	Multi-site Y型多功能/3.0m	C8					
NT3A,NT3B, NT3E,NT3F,NT3B60 (*Digital Tech)			U410-173D	Adult finger 成人指夹/3.0m	D1	No need 无需转接线				
			U410S-173D	Adult soft 成人软指套/3.0m	D2					
			U810-173D	Multi-site Y型多功能/3.0m	D1					
<b>Nihon Kohden 光电</b>										
BSM-1101/2101 /23101 7100,7200, OLV-1100 ALife scope			U410-16	Adult finger 成人指夹/3.0m	C8	No need 无需转接线				
			U110-16	Pediatric finger 儿童指夹/3.0m	D1					
			U410S-16	Adult soft 成人软指套/3.0m	D1					
			U110S-16	Pediatric soft 儿童软指套/3.0m	D1					
			U210S-16	Infant soft 婴儿软指套/3.0m	D1					
			U310-16	Neonate wrap 新生儿捆绑/3.0m	D1					
			U610-16	Adult/Neo wrap 成人/新生儿捆绑/3.0m	D1					
			U610B-16	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8					
			U810-16	Multi-site Y型多功能/3.0m	C8					
			U910-16	Adult ear clip 成人耳夹/3.0m	D1					
Lifescope S, BSS-9800 Lifescope M, BSM 9510 Lifescope P BSM 4100 Life scope A, BSM 5100 Lifescope L, BSM 2350 Lifescope I BSM 2301K, 2303K		TL-201T	U405-16	Adult finger 成人指夹/1.5m	C7	U710-36				
			U105-16	Pediatric finger 儿童指夹/1.5m	C8					
			U405S-16	Adult soft 成人软指套/1.5m	C8					
			U105S-16	Pediatric soft 儿童软指套/1.5m	C8					
			U205S-16	Infant soft 婴儿软指套/1.5m	C8					
			U305-16	Neonate wrap 新生儿捆绑/1.5m	C8					
			U605-16	Adult/Neo wrap 成人/新生儿捆绑/1.5m	C7					
			U605B-15	Pedia/infant wrap 儿童/婴儿捆绑/1.5m	C7					
		TL-220T	U805-16	Multi-site Y型多功能/1.5m	C7					
			U905-16	Adult ear clip 成人耳夹/1.5m	C8					
<b>Nihon Kohden 光电 下页续接</b>										















↳ Nihon Kohden光电 续接上页

























OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Nihon Kohden光电</b>								
Lifescop S, BSS-9800 Lifescop M, BSM 9510 Lifescop P BSM 4100 Life scope A, BSM 5100 Lifescop L, BSM 2350 Lifescop I BSM 2301K, 2303K			<b>U410-36</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			<b>U110-36</b>	Pediatric finger 儿童指夹/3.0m	D2			
			<b>U410S-36</b>	Adult soft 成人软指套/3.0m	D2			
			<b>U110S-36</b>	Pediatric soft 儿童软指套/3.0m	D2			
			<b>U210S-36</b>	Infant soft 婴儿软指套/3.0m	D2			
			<b>U310-36</b>	Neonate wrap 新生儿捆绑/3.0m	D2			
			<b>U610-36</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	D1			
			<b>U610B-36</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	D1			
			<b>U810-36</b>	Multi-site Y型多功能/3.0m	D1			
	<b>U910-36</b>	Adult ear clip 成人耳夹/3.0m	D2					
OPV-1500		The same as above	The same as above	1.1m	The same as above	U710-106		
<b>Nonin燕牌</b>								
8500 2500 9840 Avant 4000,9600, 9700 Avant 2120 PulseSense Model 7500,7500FO (Nonin*)		8000AA	<b>U403-08</b>	Adult finger 成人指夹/1.1m	C3	No need 无需转接线		
		8000AP	<b>U103-08</b>	Pediatric finger 儿童指夹/1.1m	C5			
		8000SL	<b>U403S-08</b>	Adult soft 成人软指套/1.1m	C5			
		8000SM	<b>U103S-08</b>	Pediatric soft 儿童软指套/1.1m	C5			
			<b>U203S-08</b>	Infant soft 婴儿软指套/1.1m	C5			
			<b>U303-08</b>	Neonate wrap 新生儿捆绑/1.1m	C5			
			<b>U603-08</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C5			
			<b>U603B-08</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C3			
			<b>U803-08</b>	Multi-site Y型多功能/1.1m	C3			
		8000Q2	<b>U903-08</b>	Adult ear clip 成人耳夹/1.1m	C5			
	<b>U410-08L</b>	Adult finger 成人指夹/3.0m	C5					
8600		The same as above	The same as above	1.1m	The same as above	NA		
8604		The same as above	The same as above	1.1m	The same as above	U708-08		
			<b>U410-08</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
Avant 4000 Wearable Digital Pulse Oximetry			<b>U401-146</b>	Adult finger 成人指夹/0.3m	D1	No need 无需转接线		
			<b>U101-146</b>	Pediatric finger 儿童指夹/0.3m	D1			
			<b>U401S-146</b>	Adult soft 成人软指套/0.3m	D1			
			<b>U101S-146</b>	Pediatric soft 儿童软指套/0.3m	D1			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Novametrix诺医</b>								
512/513/710/715, MARS			<b>U403-73</b>	Adult finger 成人指夹/1.1m	C5	No need 无需转接线		
			<b>U103-73</b>	Pediatric finger 儿童指夹/1.1m	C7			
			<b>U403S-73</b>	Adult soft 成人软指套/1.1m	C7			
			<b>U103S-73</b>	Pediatric soft 儿童软指套/1.1m	C7			
			<b>U203S-73</b>	Infant soft 婴儿软指套/1.1m	C7			
			<b>U303-73</b>	Neonate wrap 新生儿捆绑/1.1m	C7			
			<b>U603-73</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C5			
			<b>U603B-73</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C5			
			<b>U803-73</b>	Multi-site Y型多功能/1.1m	C5			
			<b>U903-73</b>	Adult ear clip 成人耳夹/1.1m	C7			
			<b>U403-103</b>	Adult finger 成人指夹/1.1m	C8	NA		
			<b>U103-103</b>	Pediatric finger 儿童指夹/1.1m	D1			
			<b>U403S-103</b>	Adult soft 成人软指套/1.1m	D1			
			<b>U103S-103</b>	Pediatric soft 儿童软指套/1.1m	D1			
			<b>U203S-103</b>	Infant soft 婴儿软指套/1.1m	D1			
			<b>U303-103</b>	Neonate wrap 新生儿捆绑/1.1m	D1			
			<b>U603-103</b>	Adult/Neo wrap成人/新生儿捆绑/1.1m	C8			
			<b>U603B-103</b>	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C8			
			<b>U803-103</b>	Multi-site Y型多功能/1.1m	C8			
			<b>U903-103</b>	Adult ear clip 成人耳夹/1.1m	D1			
505/510/511 /515/B/C/520 7000/7100,Oxyleth (Novametrix*)			<b>U410-03</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U110-03</b>	Pediatric finger 儿童指夹/3.0m	D1			
			<b>U410S-03</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U110S-03</b>	Pediatric soft 儿童软指套/3.0m	D1			
			<b>U210S-03</b>	Infant soft 婴儿软指套/3.0m	D1			
			<b>U310-03</b>	Neonate wrap 新生儿捆绑/3.0m	D1			
			<b>U610-03</b>	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			<b>U610B-03</b>	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			<b>U810-03</b>	Multi-site Y型多功能/3.0m	C8			
			<b>U910-03</b>	Adult ear clip 成人耳夹/3.0m	D1			
<b>PaceTech/Generra</b>								
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-38		
			<b>U410-38</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
		<b>U410S-38</b>	Adult soft 成人软指套/3.0m	C8				
		<b>U810-38</b>	Multi-site Y型多功能/3.0m	C7				

↳ PaceTech/Generra 下页续接

PaceTech/Generra 续接上页

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
PaceTech/Generra								
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-39		
			U410-39	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-39	Adult soft 成人软指套/3.0m	D1			
			U810-39	Multi-site Y型多功能/3.0m	C8			
Mediaid/Palco								
M300/M305/340, M5305/M5340			U410-18	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-18	Pediatric finger 儿童指夹/3.0m	D1			
			U410S-18	Adult soft 成人软指套/3.0m	D1			
			U110S-18	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-18	Infant soft 婴儿软指套/3.0m	D1			
			U310-18	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-18	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			U610B-18	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			U810-18	Multi-site Y型多功能/3.0m	D1			
M30/40,M100,M130, M900,M960P, M6 Aster,M4 Nova, M1000			U410-148	Adult finger 成人指夹/3.0m	D1.3	No need 无需转接线		
			U110-148	Pediatric finger 儿童指夹/3.0m	D1.6			
			U410S-148	Adult soft 成人软指套/3.0m	D1.6			
			U110S-148	Pediatric soft 儿童软指套/3.0m	D1.6			
			U210S-148	Infant soft 婴儿软指套/3.0m	D1.6			
			U310-148	Neonate wrap 新生儿捆绑/3.0m	D1.6			
			U610-148	Adult/Neo wrap成人/新生儿捆绑/3.0m	D1.3			
			U610B-148	Pedia/infant wrap儿童/婴儿捆绑/3.0m	D1.3			
			U810-148	Multi-site Y型多功能/3.0m	D1.6			
Philips/HP飞利浦								
78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax)		DS100A	U403-01	Adult finger 成人指夹/1.1m	C3	U708-20		
		DS100A	U403S-01	Adult soft 成人软指套/1.1m	C5			
		D-Y-S	U803-01	Multi-site Y型多功能/1.1m	C3			
		M1196A	U403-91 U410-91	Adult finger 成人指夹/1.1m 3.0m	C8	U708-40		
		M1192A	U103-91	Pediatric finger 儿童指夹/1.1m	D1			
		M1191A	U403S-91 U410S-91	Adult soft 成人软指套/1.1m 3.0m	D1			
		M1192A	U103S-91	Pediatric soft 儿童软指套/1.1m	D1			
		M1195A	U203S-91	Infant soft 婴儿软指套/1.1m	D1			
		M1193A	U303-91	Neonate wrap 新生儿捆绑/1.1m	D1			
			U603-91	Adult/Neo wrap成人/新生儿捆绑/1.1m	C8			
			U603B-91	Pedia/infant wrap儿童/婴儿捆绑/1.1m	C8			
			U803-91	Multi-site Y型多功能/1.1m	D1			
		M1194A	U903-91	Adult ear clip 成人耳夹/1.1m	D1			

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
Philips/HP飞利浦								
78352C,78354C, 78834C,M1020A, M1025B,M1025A, M1350B,M1772A/B, M1732A/B,M2475B (*Non-Oximax)		M1190A	U410-20	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U110-20	Pediatric finger 儿童指夹/3.0m	D1			
		M1190A	U410S-20	Adult soft 成人软指套/3.0m	D1			
			U110S-20	Pediatric soft 儿童软指套/3.0m	D1			
			U210S-20	Infant soft 婴儿软指套/3.0m	D1			
			U310-20	Neonate wrap 新生儿捆绑/3.0m	D1			
			U610-20	Adult/Neo wrap成人/新生儿捆绑/3.0m	C8			
			U610B-20	Pedia/infant wrap儿童/婴儿捆绑/3.0m	C8			
			U810-20	Multi-site Y型多功能/3.0m	C8			
			U910-20	Adult ear clip 成人耳夹/3.0m	D1			
C1,C3,M3,M4, VirDia, MP20,MP30, MP40/50/60/70/90. VM4,VM6,VM8 (Philips*)		See Philips' sensor	See Philips' sensor	See Philips' sensor	See Philips' sensor	U708-41		
		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-43		
A1, A3 (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-01		
MP20/30/40/50 PM60/70/90 VN4,VM6,VM8 (*Oximax)		DS100A	U403-01P	Adult finger 成人指夹/1.1m	C5	U710X-43		
		DS100A	U403S-01P	Adult soft 成人软指套/1.1m	C7			
		D-Y-S	U803-01P	Multi-site Y型多功能/1.1m	C5			
MP20/30/40/ PM50/60/70/90 VN4,VM6,VM8 (*M tech)		See M	See M	1.1m	See M	U708M-43		
			*U410-91M	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			*U410S-91M	Adult soft 成人软指套/3.0m	D2			
			*U810-91M	Multi-site Y型多功能/3.0m	D1			
RGB medical								
(*BCI)		3044	U403-06	Adult finger 成人指夹/1.1m	C3	U708-76		
			U403S-06	Adult soft 成人软指套/1.1m	C3			
		3078	U803-06	Adult ear clip 成人耳夹/1.1m	C5			
Schiller ODM/Bruker								
CM8 SM783,784,785, TM910		See RGB medical	See RGB medical	See RGB medical	See RGB medical	U708-17		
			U410-17	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			U410S-17	Adult soft 成人软指套/3.0m	D1			
			U810-17	Multi-site Y型多功能/3.0m	C8			

Schiller ODM/Bruker 下页续接


OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头	
<b>Schiller ODM/Brüker</b>									
TM7, Argus LCM, (*M tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-17			
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6				
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4				
				<b>*U410-17M</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
				<b>*U410S-17M</b>	Adult soft 成人软指套/3.0m	D1			
				<b>*U810-17M</b>	Multi-site Y型多功能/3.0m	C8			
Propaq CS, Encore, Micropaq (New) (M) (*M tech)		See M	See M	1.1m	See M	U708M-15			
		LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8				
		LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1	No need 无需转接线			
			<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8				
Propaq CS, Encore Series, Atlas (* non-Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-01			
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5				
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3				
TM 910-Defigard 500(M) (* non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m/3.0m	See Nellcor	U708-63			
Argus VCM (*Oximax)		DS100A	<b>U403-01X</b>	Adult finger 成人指夹/1.1m	C5	U710-70			
		D-YS	<b>U803-01X</b>	Multi-site Y型多功能/1.1m	C5				
<b>Spacelabs太空</b>									
90419/01, 90651A/07, 90501, Ultraview 90496 (*Novamatrix)			<b>U410-03</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线			
			<b>U110-03</b>	Pediatric finger 儿童指夹/3.0m	D1				
			<b>U410S-03</b>	Adult soft 成人软指套/3.0m	D1				
			<b>U110S-03</b>	Pediatric soft 儿童软指套/3.0m	D1				
			<b>U210S-03</b>	Infant soft 婴儿软指套/3.0m	D1				
			<b>U310-03</b>	Neonate wrap 新生儿捆绑/3.0m	D1				
			<b>U610-03</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8				
			<b>U610B-03</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8				
			<b>U810-03</b>	Multi-site Y型多功能/3.0m	C8				
90351-0/6, 90465, 90466, 90456, PC exress (*Non-Oximax)		DS100A	<b>U403-01</b>	Adult finger 成人指夹/1.1m	C3	U708-04			
		DS100A	<b>U403S-01</b>	Adult soft 成人软指套/1.1m	C5				
		D-YS	<b>U803-01</b>	Multi-site Y型多功能/1.1m	C3				
				<b>U410-04</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
				<b>U110-04</b>	Pediatric finger 儿童指夹/3.0m	D1			
				<b>U410S-04</b>	Adult soft 成人软指套/3.0m	D1			
				<b>U110S-04</b>	Pediatric soft 儿童软指套/3.0m	D1			
				<b>U210S-04</b>	Infant soft 婴儿软指套/3.0m	D1			
				<b>U310-04</b>	Neonate wrap 新生儿捆绑/3.0m	D1			
				<b>U610-04</b>	Adult/Neo wrap 成人/新生儿捆绑/3.0m	C8			
				<b>U610B-04</b>	Pedia/infant wrap 儿童/婴儿捆绑/3.0m	C8			
				<b>U810-04</b>	Multi-site Y型多功能/3.0m	C8			
		<b>U910-04</b>	Adult ear clip 成人耳夹/3.0m	D1					

OEM Manufacturer & Models 原厂设备及型号	Sensor Connector 探头接头端	REF# OEM 原厂型号	REF# UNIMED 尤迈型号	Sensor Descriptions 探头描述	Price Code 价格代码	Ext - Cable 转接线	Cable Plug 机器端接头	Sensor Conn. 探头端接头
<b>Spacelabs太空</b>								
Ultraview (*Non-Oximax)		See Nellcor Non-Oximax	See Nellcor Non-Oximax	1.1m	See Nellcor	U708-74		
			<b>U410-74</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			<b>U410S-74</b>	Adult soft 成人软指套/3.0m	D2			
			<b>U810-74</b>	Multi-site Y型多功能/3.0m	D1			
Ultraview (M Tech)			<b>U410-74M</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			<b>U410S-74M</b>	Adult soft 成人软指套/3.0m	D2			
Ultraview (Oximax Tech)			<b>U410-74P</b>	Adult finger 成人指夹/3.0m	D1	No need 无需转接线		
			<b>U410S-74P</b>	Adult soft 成人软指套/3.0m	D2			
PRO2 (*BCI)		3044	<b>U403-06</b>	Adult finger 成人指夹/1.1m	C3	No need 无需转接线		
			<b>U403S-06</b>	Adult soft 成人软指套/1.1m	C5			
		3043	<b>U803-06</b>	Multi-site Y型多功能/1.1m	C3			
<b>Zoll</b>								
M series AED (*M tech)		1863/ LNCS DCI	<b>*U403-49</b>	Adult finger 成人指夹/1.1m	C4	U708M-15		
		1863/ LNCS DCI	<b>*U403S-49</b>	Adult soft 成人软指套/1.1m	C6			
		2258/ LNCS YI	<b>*U803-49</b>	Multi-site Y型多功能/1.1m	C4			
		LNOP DC-12	<b>*U410-15</b>	Adult finger 成人指夹/3.0m	C8			
		LNOP DC-12	<b>*U410S-15</b>	Adult soft 成人软指套/3.0m	D1			
		<b>*U810-15</b>	Multi-site Y型多功能/3.0m	C8				
<b>Zondan和心重典 (China)</b>								
Apollo N2, Apollo N4 Apollo N5 (M tech optional) DS120C			<b>U410-45</b>	Adult finger 成人指夹/3.0m	C7	No need 无需转接线		
			<b>U410S-45</b>	Adult soft 成人软指套/3.0m	C8			
			<b>U810-45</b>	Multi-site Y型多功能/3.0m	C7			
ZD120B ZD120D ZD120E (*Digital Tech)			<b>U410-66D</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-66D</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-66D</b>	Multi-site Y型多功能/3.0m	C8			
<b>Unicare(China)</b>								
(*Digital tech)			<b>U410-66D</b>	Adult finger 成人指夹/3.0m	C8	No need 无需转接线		
			<b>U410S-66D</b>	Adult soft 成人软指套/3.0m	D1			
			<b>U810-66D</b>	Multi-site Y型多功能/3.0m	C8			

Please be Reminded : The other types of sensors like U100, U100S, U200S, U300, U600, U600B and U900 of the listed respected brands are available upon demands.


注意：其他类型的兼容于各厂家的探头，如U100儿童软指套，U200婴儿软指套，U300新生儿捆绑式，U600包裹式及U900成人耳夹探头都可以选择。





## Compatible with major patient monitor brands 兼容主要品牌

Nellcor (耐普特)      BCI  
 Ohmeda (欧美达)      Datascope  
 Novamatrix (诺医)      Datex (德恩)  
 Criticare                  Nihon Kohden (光电)  
 ...







## ② Disposable SpO<sub>2</sub> Sensor 一次性血氧探头

### Specification 产品规格

Cable material 线材	PVC Jacket, white, 3.2mm diameter PVC 外被, 白色, 3.2mm 直径	Connector 接头	PVC plug & Golden plated pins PVC 插头 & 镀金插针
Wavelength 波长	660nm/905nm/940nm	Photo Detector (size) 接收管 (大小)	4mm/8mm (size) 4mm/8mm (面积)
Cable length 线长	0.9m	Sensor tape (See the Next picture for Reference) 探头绑带 (见下页图)	Blue non-adhesive foam (N500 series) 3M textile adhesive (U500 series) 3M white adhesive foam (F500 series) 蓝色海绵 (N500 系列) 3M 无纺布 (U500 系列) 3M 白色海绵 (F500 系列)
Package 包装	24 pcs per box 24个/盒		

OEM REF# 原产型号	Unimed REF# Textile Adhesive 无纺布粘性	Unimed REF# Foam Adhesive 海绵粘性	Unimed REF# Non-Adhesive 海绵非粘性	Patient Sizes 病人适用范围	Sensor Connector 探头接头端	Price Code (Per pcs) 价格代码 (个)	Adapter Cable required 连接电缆要求
<b>Nellcor Non oximax</b>							
D25	U503-01	F503-01	N543-01	Adult(>30kg) 成人(>30kg)		B1	1. Standard cables 2. U701-01X or U701-01P is required used with Oximax See Spo2 reusable section
D20	U523-01	F523-01	N533-01	Pediatric(15-40kg) 儿童(15-40kg)		B1	
I20	U533-01	F533-01	N533-01	Infant(1-15kg) 婴儿(1-15kg)		B1	
N25	U543-01	F543-01	N543-01	Neonate(<3kg) 新生儿(<3kg)		B1	
<b>Ohmeda*欧美达</b>							
SAS-AF	U503-02	F503-02	N543-02	Adult(>30kg) 成人(>30kg)		B2	U708-72 cable is required 需配转接线U708-72
SAS-AP	U523-02	F523-02	N533-02	Pediatric(15-40kg) 儿童(15-40kg)		B2	
SAS-AP	U533-02	F533-02	N533-02	Infant(1-15kg) 婴儿(1-15kg)		B2	
SAS-AF	U543-02	F543-02	N543-02	Neonate(<3kg) 新生儿(<3kg)		B2	

OEM REF# 原产型号	Unimed REF# Textile Adhesive 无纺布粘性	Unimed REF# Foam Adhesive 泡棉粘性	Unimed REF# Non-Adhesive 泡棉非粘性	Patient Sizes 病人适用范围	Sensor Connector 探头接头端	Price Code (Per pcs) 价格代码 (个)	Adapter Cable required 连接电缆要求
<b>Nonin燕牌</b>							
6000CA/7000A	U503-08	F503-08	N543-08	Adult(>30kg) 成人(>30kg)		B1	Standard cables 标准转接线
6000CP/7000P	U523-08	F523-08	N533-08	Pediatric(15-40kg) 儿童(15-40kg)		B1	
6000CI/7000I	U533-08	F533-08	N533-08	Infant(1-15kg) 婴儿(1-15kg)		B1	See Spo2 reusable section 参照重复性探头部分
6000CN/7000N	U543-08	F543-08	N543-08	Neonate(<3kg) 新生儿(<3kg)		B1	
<b>Criticare</b>							
570SD	U503-05	F503-05	N543-05	Adult(>30kg) 成人(>30kg)		B2	Standard cables 标准转接线
571SD	U523-05	F523-05	N533-05	Pediatric(15-40kg) 儿童(15-40kg)		B2	
572SD	U533-05	F533-05	N533-05	Infant(1-15kg) 婴儿(1-15kg)		B2	See Spo2 reusable section 参照重复性探头部分
573SD	U543-05	F543-05	N543-05	Neonate(<3kg) 新生儿(<3kg)		B2	
<b>Nihon Kohden光电</b>							
TL-251T	U503-16	F503-16	N543-16	Adult(>30kg) 成人(>30kg)		B2	Standard cables 标准转接线
TL-252T	U523-16	F523-16	N533-16	Pediatric(15-40kg) 儿童(15-40kg)		B2	
TL-252T	U533-16	F533-16	N533-16	Infant(1-15kg) 婴儿(1-15kg)		B2	See Spo2 reusable section 参照重复性探头部分
TL-253T	U543-16	F543-16	N543-16	Neonate(<3kg) 新生儿(<3kg)		B2	
<b>BCI</b>							
1300	U503-06	F503-06	N543-06	Adult(>30kg) 成人(>30kg)		B1	Standard cables 标准转接线
1301	U523-06	F523-06	N533-06	Pediatric(15-40kg) 儿童(15-40kg)		B1	
1303	U533-06	F533-06	N533-06	Infant(1-15kg) 婴儿(1-15kg)		B1	See Spo2 reusable section 参照重复性探头部分
1302	U543-06	F543-06	N543-06	Neonate(<3kg) 新生儿(<3kg)		B1	
<b>Datex德恩</b>							
SAS-AF	U503-09	F503-09	N543-09	Adult(>30kg) 成人(>30kg)		B1	Standard cables 标准转接线
SAS-AP	U523-09	F523-09	N533-09	Pediatric(15-40kg) 儿童(15-40kg)		B1	
SAS-AP	U533-09	F533-09	N533-09	Infant(1-15kg) 婴儿(1-15kg)		B1	See Spo2 reusable section 参照重复性探头部分
SAS-AF	U543-09	F543-09	N543-09	Neonate(<3kg) 新生儿(<3kg)		B1	
<b>Datascope</b>							
0998-00-0076-05	U503-07	F503-07	N543-07	Adult(>30kg) 成人(>30kg)		B2	Standard cables 标准转接线
0998-00-0076-04	U523-07	F523-07	N533-07	Pediatric(15-40kg) 儿童(15-40kg)		B2	
0998-00-0074-03	U533-07	F533-07	N533-07	Infant(1-15kg) 婴儿(1-15kg)		B2	See Spo2 reusable section 参照重复性探头部分
0998-00-0074-04	U543-07	F543-07	N543-07	Neonate(<3kg) 新生儿(<3kg)		B2	

OEM REF# 原产型号	Unimed REF# Textile Adhesive 无纺布粘性	Unimed REF# Foam Adhesive 泡棉粘性	Unimed REF# Non-Adhesive 泡棉非粘性	Patient Sizes 病人适用范围	Sensor Connector 探头接头端	Price Code (Per pcs) 价格代码 (个)	Adapter Cable required 连接电缆要求
<b>Novamatrix*诺医</b>							
AS110	U503-03	F503-03	N543-03	Adult(>30kg) 成人(>30kg)		B2	U708-03 cables required 需配转接线U708-03
AS120	U523-03	F523-03	N533-03	Pediatric(15-40kg) 儿童(15-40kg)		B2	
AS120	U533-03	F533-03	N533-03	Infant(1-15kg) 婴儿(1-15kg)		B2	See Spo2 cable section 参照重复性探头部分
AS140	U543-03	F543-03	N543-03	Neonate(<3kg) 新生儿(<3kg)		B2	
<b>MEK</b>							
	U503-26	F503-26	N543-26	Adult(>30kg) 成人(>30kg)		B2	Standard cables 标准转接线
	U523-26	F523-26	N533-26	Pediatric(15-40kg) 儿童(15-40kg)		B2	
	U533-26	F533-26	N533-26	Infant(1-15kg) 婴儿(1-15kg)		B2	See Spo2 reusable section 参照重复性探头部分
	U543-26	F543-26	N543-26	Neonate(<3kg) 新生儿(<3kg)		B2	
<b>GE OxyTip+</b>							
OXY-AP-10	U503-127	F503-127	N543-127	Adult(>30kg) 成人(>30kg)		B4	Standard cables 标准转接线
OXY-AP-10	U523-127	F523-127	N533-127	Pediatric(15-40kg) 儿童(15-40kg)		B4	
OXY-AP-10	U533-127	F533-127	N533-127	Infant(1-15kg) 婴儿(1-15kg)		B4	See Spo2 reusable section 参照重复性探头部分
OXY-AP-10	U543-127	F543-127	N543-127	Neonate(<3kg) 新生儿(<3kg)		B4	
<b>GE TruSignal+</b>							
TS-AF-10	U503-117	F503-117	N543-117	Adult(>30kg) 成人(>30kg)		B4	Standard cables 标准转接线
TS-AF-10	U523-117	F523-117	N533-117	Pediatric(15-40kg) 儿童(15-40kg)		B4	
TS-AF-10	U533-117	F533-117	N533-117	Infant(1-15kg) 婴儿(1-15kg)		B4	See Spo2 reusable section 参照重复性探头部分
TS-AF-10	U543-117	F543-117	N543-117	Neonate(<3kg) 新生儿(<3kg)		B4	



TruSignal+ Disposable Sensor



### ③ Interface Cable

## 血氧转接电缆

#### Specification

#### 产品规格

Cable material 电缆材料	TPU Jacket, Light gray, Blue TPU 外壳, 浅灰色, 蓝色
Connector 插头	Nylon plug & Golden plated pins 尼龙 插头 & 镀金插针
Cable length 电缆长度	1.1m/2.2m/3.0m
Cable diameter 电缆直径	4.0mm/5.0mm
Package 包装	1 pcs per bag 每条一包



REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Artema 阿替玛						
U708-80	59831-0	Nellcor,non-Oximax	2.2m			C3
BCI						
U708-76		BCI	2.2m			B8
U708-06		BCI	2.2m			B8
Biolight 宝莱特						
U710X-51		Nellcor Oximax	3.0m			C3
U708-68		Digital sensor 数字探头	2.2m			C3
Bionet						
U708-58		Nellcor,non-Oximax	2.2m			B8
U708-108		Bionet	2.2m			C3
Criticare						
U708-05	518LD	CSI	2.2m			C3
U708-75	518DD	CSI	2.2m			B8
Datascope						
U708-07	0012-00-0516-02	DATASCOPE	2.2m			C3
Dolphin 海豚						
U708-50		Dolphin digital sensor Dolphin数字探头	2.2m			C3
Drager 德尔格						
U708-28	5720084	Nellcor,non-Oximax	2.2m			D4
EDAN 理邦						
U708X-191	MS1-30131 (M9,M9A,M8,M8A,M8B)	Nellcor-Oximax	2.2m			C3
U708-29	M13-30691 (M9,M9A,M8,M8A,M8B, M3,M3B)	EDAN Oximetry	2.2m			B8

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
GE healthcare						
U708-21	2006644-001 E9004GE	Nellcor,non-Oximax	2.2M			C3
U710X-21	2021406-001	Nellcor,Oximax	3.0M			C7
*U708-21M	M technology	Nellcor Non-Oximax	2.2M			C3
*U708M-21	LNC-10-GE 2027263-002 (REF 2016)	M tech.(Masimo LNCS)	2.2M			C3
*U708M-21P	20169041-001	M tech.(Masimo LNOP)	2.2M			C7
U708-09	OXY-C3	Datex	2.2m			C3
U708-81	Datex>>Oxytip Used with Oxy-C3 cable	Oxytip	2.2m			C3
U708-55	30344358	Hellige(Nellcor)	2.2m			C3
U708-25	4033CAX	Corometrics(Nellcor)	2.2m			C3
U708-72*	All Ohmeda Oximeters	Ohmeda Used with Unimed N500 or U500 Disposable sensors	2.2m			C3
U708-160	OXY-MC3	GE OxyTip+ sensor	2.2m			C7
U708-102 U704-102	OXY-OL3	GE OxyTip+ sensor	2.2m 1.1m			C5
U708-209	OXY-SL3	GE OxyTip+ sensor	2.2m			C5
U708-121 U710-121	OXY-ES3	GE OxyTip+ sensor	2.2m 3.0m			C7
U708-81 U701-81	OXY-SLA OXY-SLC	GE OxyTip+ sensor	2.2m 0.3m			C5
U708-111	TS-G3	GE OxyTip+ sensor	2.2m			C5
U708-112	TS-H3	GE OxyTip+ sensor	2.2m			C5
U708-116	TS-M3	GE OxyTip+ sensor	2.2m			C7
U708-119	TS-N3	GE OxyTip+ sensor	2.2m			C5
Goldway 金科威						
U708-11		Nellcor Non-Oximax	2.2m			C3
U708-37		Nellcor Non-Oximax	2.2m			B8

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Goldway 金科威						
U708-38		Nellcor Non-Oximax	2.2m			B8
U710X-38		Nellcor Oximax	3.0m			C3
Invivo 鹰威						
U708-14		Nellcor,non-Oximax	2.2m			C3
*U708-14M		(Masimo) Nellcor, non-Oximax	2.2m			C3
*U708M-14		Masimo LNCS sensor	2.2m			C3
Kontron 康强						
U708-13		Nellcor,non-Oximax	2.2m			C3
U708-53		Nellcor,non-Oximax	2.2m			B8
M Technology						
*U708-15	M technology	(Masimo) Nellcor, non-Oximax	2.2m			C3
*U708M-15	LNC-4(REF 2017) LNC-10(REF 1814) LNC-14(REF 2013)	Masimo LNCS sensor	2.2m			C3
*U708M-15P	PC 04 / PC 08 / PC 12	Masimo LNOP sensor	2.2m			C3
*U704M-49	LNC-4-Ext (REF 2021) For LNCS Cable	Masimo LNCS sensor	1.1m			B8
*U701M-125	LNCS to LNOP LNCS to PC series (REF 1816)	Masimo LNCS sensor	0.3m			C3
MEK						
U708-56		MEK (Nellcor non-oximax)	2.2m			C3
Mennen 美能						
U708-54		Nellcor,non-Oximax	2.2m			C3
U708-24		Nellcor,non-Oximax	2.2m			D4
*U708-59M	M technology	M tech. (Nellcor, non-Oximax)	2.2m			C3
*U708M-59	M technology	Masimo LNCS sensor	2.2m			C3
U710X-59		Nellcor Oximax	3.0m			C3

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Mindray/迈瑞						
U708-29	0010-20-42594 Accutorr V, DPM 3, DPM 4, DPM 5, Passport V	Nellcor,non-Oximax	2.2M			B8
*U708-29M	M technology For old PM 6000 only	Nellcor,non-Oximax	2.2M			B8
*U708M-29	M technology For old PM 6000 only	Masimo LNCS	2.2M			B8
*U708-30M	M technology (0010-30-42625) DPM 3	Nellcor,non-Oximax	2.2M			B8
*U708M-30	DPM 4 DPM 5	Masimo LNCS	2.2M			B8
U708-22	0515-30-11221	Nellcor,non-Oximax	2.2M			C3
*U708-22M	M technology 0010-30-12452	non-Oximax Nellcor	2.2M			C3
*U708M-22	M technology 0010-30-12452	Masimo LNCS	2.2M			C3
U708-31	CSI Technology	CSI	2.2M			C3
*U708M-48	0010-30-42738 T5/T8 DPM 6, DPM 7	Masimo LNCS	2.2m			C3
U708-48	0010-20-42710 (C) T5/T8 DPM 6, DPM 7 Nellcor Non-Oximax	Nellcor,non-Oximax	2.2m			C3
U710X-48	0010-20-42712 T5/T8, DPM 6, DPM 7 Nellcor Oximax	Nellcor Oximax	3.0m			C4
*U708M-15	(0012-00-1652),1.2m (0012-00-1599),3.1m (0012-00-1653),4.3m Accutorr V, Duo, Trio, Passport 2, PassportV, Spectrum/OR, V12, V21	Masimo LNCS				C3
U708-69	0012-00-1254 Trio, Passport 2, Spectrum	Nellcor,non-Oximax	2.2m			C3
U710-70	0012-00-1464 Accutorr V, Duo, Trio, Passport 2, Passport V, Spectrum/OR, V12, V21	Nellcor Oximax	3.0m			C6
U710X-180	0010-20-42595 DPM4	Nellcor Oximax	3.0m			C3
Nellcor/Covidien 奈普特						
U708-69	SCP-10/MC-10	Nellcor, Non-Oximax	2.2m			C3
U710-85	M-200-13	Nellcor, Non-Oximax	3.0m			D7
U704-71	DEC-4	Nellcor-Oximax	1.1m			B8
U708-71	DEC-8	Nellcor-Oximax	2.2m			B8
U710-70	DOC-10	Nellcor-Oximax	3.0m			C6
U710-70X	DOC-10 + Oximax adapter	Nellcor- Non Oximax	3.0m			C8

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Nellcor/Covidien 奈普特						
U701-01X U708-01X	Oximax adapter	Nellcor, Non-Oximax	0.3m 2.2m			C3
Nihon Kohden 光电						
U710-36	JL-900P BSM/PVM	NIHON KOHDEN	3.0m			C7
U710X-36	JL-650P	Oximax	3.0m			C7
U710M-36	JL-631P	M Tech	3.0m			C7
U710-106	JL-302T OPV-1500,OLV	NIHON KOHDEN	3.0m			C7
Nonin 燕牌						
U708-08		Nonin	2.2m			C3
Novamatrix 诺医						
U708-103	8853	Novamatrix	2.2m			C3
Pace tech/Generra						
U708-37		Nellcor non-Oximax	2.2m			B8
U708-38		Nellcor non-Oximax	2.2m			B8
U708-39		Nellcor non-Oximax	2.2m			C3
Philips/HP 飞利浦						
U708-40	M1940A/B	Philips	2.2m			C3
U708-41	M1941A/B		2.2m			C3
U708-20	M1900A/B	Nellcor non-Oximax	2.2m			C3
*U708M-20P	M1020-61102 CMS M1020A module.	Masimo LNOP	2.2m			C3
U708-43 U710-43	M1943A M1943AL	Nellcor non-Oximax	2.2m 3.0m			C3
U710X-43	M1943NL	Nellcor,Oximax	3.0m			C7
*U708-43M	M Technology	Nellcor non-Oximax	2.2m			C3
*U708M-43	989803148221 LNC MP10 (REF 2281)	Masimo LNCS sensor	2.2m			C3
*U708M-43P	M1020-61100	Masimo LNOP sensor	2.2m			C7



# 4 Accessories 配件

REF# 型号	OEM REF# 原厂型号	Sensor Tech. 传感器技术	Cable length 电缆线长度	Monitor end Conn. 机器端接头	Sensor end Conn. 探头端接头	Price code 价格代码
Schiller席勒						
U708-17		BCI	2.2m			C3
*U708-17M	M tech.	Nellcor non-Oximax	2.2m			C3
*U708M-17	M tech.	Masimo LNCS sensor	2.2m			C3
U708-63 ODAM/Bruker		BCI	2.2m			C3
Siemens/Drager 西门子/德尔格						
U704-23	3368433	Nellcor non-Oximax	1.1m			C6
U708-23	3375834	Nellcor non-Oximax	2.2m			C6
U710X-23		Nellcor Oximax	3.0m			C7
*U704M-23P	MS18680	Masimo* LNOP	1.1m			C6
*U708M-23P	MS17041	Masimo* LNCS	2.2m			C6
*U704M-23	MS24303	Masimo* LNCS	1.1m			C6
*U708M-23	MS17522	Masimo* LNCS	2.2m			C6
U708-33	4533840	Nellcor non-Oximax	2.2m			C3
Spacelabs 太空						
U708-04	700-0002-00 90351-D/G, 90465, 90466, 90467, 90489,M77	Nellcor non-Oximax	2.2m			C3
U708-74	700-0030-00 90496,91496-U	Nellcor non-Oximax	2.2m			C6
U710X-74	700-0792-00 91496-N,91220	Nellcor-Oximax	3.0m			C7
*U708M-74	700-0906-00 91496-M	Masimo LNCS	2.2m			C6
*U708M-74P	700-0789-00 91496-M	Masimo LNCS	2.2m			C6
U701-04NV	175-0646-00 90351-D/G,90465, 90466,90467, 90489, 90651-08)	Novamatrix- Hypertronic Con	0.3m			C3
U701-74NV	700-0029-00 (90496,91496-U)	Novamatrix- Hypertronic Con	0.3m			C6
Welchallyn 伟伦						
U704X-104	001966-E (PIC 40TM Defibrillator PIC 50TM Defibrillator) 001955-U	Nellcor Oximax	2.2m			C7
U704-104	(PIC40 Defibrillator PIC50 Defibrillator)	Nellcor	1.1m			C7

\*Refurbish only



REF# 型号	Product Name 品名	Purpose of Use 用途	Package 包装	Price 价格	Figures 图标
Y04-E041	Ear sensor hanger 成人耳夹挂钩	Used with Unimed ear clip sensor U900 series or Other brands ear sensor(Gray) 用于尤迈900系列耳夹探头或其它品牌探头(灰)	1 pcs	A5	A
UEC-01	Adult ear clip Or as animal small size lingua clip 成人耳夹或作为小动物舌夹	Used with D-YS or Unimed U800 sensor(Black) 用于尤迈800系列或者D-YS探头(黑)	1 pcs	A5	B
UVC-01	Animal lingua clip,Big size 大号动物夹	Used with D-YS or Unimed U800 sensor(Gray) 用于尤迈800系列或者D-YS探头(灰)	1 pcs	A8	C
UVC-02	Animal lingua clip,Middle size 中号动物夹	Used with D-YS or Unimed U800 sensor(Gray) 用于尤迈800系列或者D-YS探头(灰)	1 pcs	A7	D
Y04-E019	Reusable Silicon Y sensor wrap tape 重复用Y型探头硅胶绑带	Used with Nellcor D-YS or Unimed U800 series sensors 用于尤迈800系列或者D-YS探头	1 pcs	A6	E
E023-12	Reusable 3 holes foam tape For All patient sizes 重复用3孔泡绵绑带	Used with Nellcor D-YS or Unimed U800 /U600 series sensors 用于尤迈800系列或者D-YS探头	12 pcs Per bag	B5/bag	F
E024-25	Disposable 3M adhesive textile tape For Adult and neonate 一次性3M无纺布绑带 (成人, 新生儿)	Used with Nellcor D-YS or Unimed U800 series sensors 用于尤迈800系列或者D-YS探头	25 pcs Per bag	B7/bag	G
E025-25	Disposable 3M adhesive textile tape for Pediatric and infant 一次性3M无纺布绑带 (儿童, 婴儿)	Used with Nellcor D-YS or Unimed U800 series sensors 用于尤迈800系列或者D-YS探头	25 pcs Per bag	B7/bag	H
T002-12	Velcro foam wrap tape 维克罗泡绵绑带	To fix the sensor cable onto the patient wrist 用于固定电缆在病人的手腕上	12 pcs per bag	C1.4/bag	I

# Content 内容

<b>1. Reusable NIBP Cuffs 重复用无创血压袖套</b>	<b>87</b>
1.1 Reusable NIBP Bladderless Cuff 重复用无囊血压袖套	88
1.1a Cuff without connector 无接头血压袖套	88
1.1b Single tube Cuff with connector 带接头单管血压袖套	89
1.1c Double tube Cuff with connector 带接头双管血压袖套	90
1.2 Reusable NIBP Cuff with inflation bag 重复用有囊血压袖套	91
1.2a Cuff without connector 无接头血压袖套	91
1.2b Single tube Cuff with connector 带接头单管血压袖套	92
1.2c Double tube Cuff with connector 带接头双管血压袖套	93
<b>2. Disposable NIBP Cuff 一次性无创血压袖套</b>	<b>94</b>
2.1 Disposable neonate NIBP Cuff 一次性新生儿血压袖套	95
2.1a Disposable neonate TPU Cuff 一次性TPU新生儿袖套	95
2.1b Disposable neonate soft fiber Cuff 一次性无纺布新生儿袖套	95
2.2 Disposable NIBP soft fiber Cuff 一次性无纺布血压袖套	95
2.2a Disposable Single tube Cuff 一次性单管血压袖套	95
2.2b Disposable Double tube Cuff 一次性双管血压袖套	96
<b>3. Patient NIBP Adapter Air Hose 无创血压连接管</b>	<b>97</b>
<b>4. NIBP connectors &amp; tubing 无创血压接头&amp;血压管</b>	<b>101</b>
4.1 Universal Y Adapter Y型接头	102
4.2 Connectors for Cuffs and Hoses 袖套及连接管接头	104
4.3 Air Tubing 血压管	105
<b>5. Manual Inflation System 手动充气装置</b>	<b>105</b>
5.1 Reusable cuff with Inflation System 充气装置带重复用袖套	105
5.2 Disposable cuff with Inflation System 充气装置带一次性血压袖套	105
<b>6. Wall manometer accessories 壁挂式血压测量附件</b>	<b>107</b>
6.1 Inflation System with cuff 带袖套血压测量装置	107
6.2 Inflation System without cuff 不带袖套血压测量装置	108
6.3 Parts 部件	108
<b>7. Tourniquet 止血带</b>	<b>109</b>

## NIBP Cuff & Accessories

### 无创血压袖套及配件

#### Certification 产品认证



CE Mark  
CE 标志



FDA 510(k) clearance  
K120364(Disposable cuff)  
K112544(Reusable cuff)



SFDA医疗器械注册证:  
证书号码 No. 1200026

#### Compliance with 符合标准

ANSI/AAMI SP10-2002,  
EN 60601-1:1990,  
ISO10993-1,5,10:2009  
EN980:2008  
EN1041:2008



## 1.1 Reusable NIBP Bladderless Cuff 重复用无囊血压袖套



1.1a NIBP Cuff without connector  
无接头血压袖套

## 1 Reusable NIBP Cuff 重复用血压袖套

### Specification

#### 产品描述

- Reusable for multiple patient use
  - Gentle nylon and TPU materials cuff with 50cm tube
  - Large selection of 9 types cuffs matching the different sizes patients, coded by tubing colors
  - Variety of connectors to fit all major brands monitoring systems
  - Convenient to use and easy to clean
  - Latex free
- 可多病人重复使用
  - 复合尼龙，TPU材质，50厘米气管
  - 9种尺寸满足不同病患者需求，气管颜色表明袖套尺寸大小
  - 配合不同型号接头可使袖套适用于各种不同监护仪系列
  - 使用方便，易于清洗
  - 不含乳胶

Patient Size 适用病人	Limb Circumference 臂周长	REF# (Single tube) 型号 (单管)	REF # (Double tube) 型号 (双管)	Tubing Color 气管颜色	Package 包装	Unit Price 单价
Large adult long 大成人加长	35.5-46cm	U1889NS	U1889ND		1 pcs	B4
Large adult 大成人	35.5-46cm	U1869NS	U1869ND		1 pcs	B3
Adult long 成人加长	27.5-36.5cm	U1886NS	U1886ND		1 pcs	B1
Adult 成人	27.5-36.5cm	U1880NS	U1880ND		1 pcs	A8
Small adult 小成人	20.5-28.5cm	U1885NS	U1885ND		1 pcs	A8
Pediatric 儿童	13.8-21.5cm	U1881NS	U1881ND		1 pcs	A8
Infant 婴儿	9-14.8cm	U1882NS	U1882ND		1 pcs	A8
Neonate 新生儿	5.8-10.8cm	U1883NS	U1883ND		1 pcs	A8
Thigh 大腿	45-56cm	U1884NS	U1884ND		1 pcs	B5
Package 组合1: Adult 成人 × 2 Small adult 小成人 × 1 Adult long 成人加长 × 1 Large adult 大成人 × 1 Large adult long 大成人加长 × 1		U1801NS-06 U1880NS × 2 U1885NS × 1 U1886NS × 1 U1869NS × 1 U1889NS × 1	U1801ND-06 U1880ND × 2 U1885ND × 1 U1886ND × 1 U1869ND × 1 U1889ND × 1	Variety组合	6 pcs	C5
Package 组合2: Adult 成人 × 1 Small adult 小成人 × 1 Pediatric 儿童 × 1 Adult long 成人加长 × 1 Large adult 大成人 × 1 Thigh 大腿 × 1		U1802NS-06 U1880NS × 1 U1885NS × 1 U1881NS × 1 U1886NS × 1 U1869NS × 1 U1884NS × 1	U1802ND-06 U1880ND × 1 U1885ND × 1 U1881ND × 1 U1886ND × 1 U1869ND × 1 U1884ND × 1	Variety组合	6 pcs	C7

1.1b Single-tube Cuff with Connector

带接头单管血压袖套

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Cuff end Connector 连接端	Tubing Color 气管颜色	Package 包装	Unit Price 单价	
Large adult long 大成人加长	35.5-46cm	U1889NS-C12	For Philips(HP) Siemens Drager Datascope Spacelabs		1 pcs	B4	
Large adult 大成人	35.5-46cm	U1869NS-C12				1 pcs	B4
Adult long 成人加长	27.5-36.5cm	U1886NS-C12				1 pcs	B2
Adult 成人	27.5-36.5cm	U1880NS-C12				1 pcs	B1
Small adult 小成人	20.5-28.5cm	U1885NS-C12				1 pcs	B1
Pediatric 儿童	13.8-21.5cm	U1881NS-C12				1 pcs	B1
Infant 婴儿	9-14.8cm	U1882NS-C12				1 pcs	B1
Neonate 新生儿	5.8-10.8cm	U1883NS-C12				1 pcs	B1
Thigh 大腿	45-56cm	U1884NS-C12				1 pcs	B5
Large adult long 大成人加长	35.5-46cm	U1889NS-C03		For Colin Spacelabs Datascope		1 pcs	B4
Large adult 大成人	35.5-46cm	U1869NS-C03					1 pcs
Adult long 成人加长	27.5-36.5cm	U1886NS-C03				1 pcs	B2
Adult 成人	27.5-36.5cm	U1880NS-C03				1 pcs	B1
Small adult 小成人	20.5-28.5cm	U1885NS-C03				1 pcs	B1
Pediatric 儿童	13.8-21.5cm	U1881NS-C03				1 pcs	B1
Infant 婴儿	9-14.8cm	U1882NS-C03				1 pcs	B1
Neonate 新生儿	5.8-10.8cm	U1883NS-C03				1 pcs	B1
Thigh 大腿	45-56cm	U1884NS-C03				1 pcs	B5
Large adult long 大成人加长	35.5-46cm	U1889NS-C18	For Protocol			1 pcs	B4
Large adult 大成人	35.5-46cm	U1869NS-C18					1 pcs
Adult long 成人加长	27.5-36.5cm	U1886NS-C18				1 pcs	B2
Adult 成人	27.5-36.5cm	U1880NS-C18				1 pcs	B1
Small adult 小成人	20.5-28.5cm	U1885NS-C18				1 pcs	B1
Pediatric 儿童	13.8-21.5cm	U1881NS-C18				1 pcs	B1
Infant 婴儿	9-14.8cm	U1882NS-C18				1 pcs	B1
Neonate 新生儿	5.8-10.8cm	U1883NS-C18				1 pcs	B1
Thigh 大腿	45-56cm	U1884NS-C18				1 pcs	B5
Package 组合1: Adult 成人×2 Small adult 小成人×1 Adult long 成人加长×1 Large adult 大成人×1 Large adult long 大成人加长×1		U1801NS-C12-06			Variety组合	6 pcs	C5
		U1801NS-C03-06			Variety组合	6 pcs	C5
		U1801NS-C18-06		Variety组合	6 pcs	C5	
Package 组合2: Adult 成人×1 Small adult 小成人×1 Pediatric 儿童×1 Adult long 成人加长×1 Large adult 大成人×1 Thigh 大腿×1		U1802NS-C12-06		Variety组合	6 pcs	C6	
		U1802NS-C03-06		Variety组合	6 pcs	C6	
		U1802NS-C18-06		Variety组合	6 pcs	C6	

1.1c Double tube Cuff with connector

带接头双管血压袖套

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Cuff end Connector 连接端	Tubing Color 气管颜色	Package 包装	Unit Price 单价	
Large adult long 大成人加长	35.5-46cm	U1889ND-C0808	For GE Corometrics DINAMAP MDE		1 pcs	B4	
Large adult 大成人	35.5-46cm	U1869ND-C0808				1 pcs	B4
Adult long 成人加长	27.5-36.5cm	U1886ND-C0808				1 pcs	B2
Adult 成人	27.5-36.5cm	U1880ND-C0808				1 pcs	B1
Small adult 小成人	20.5-28.5cm	U1885ND-C0808				1 pcs	B1
Pediatric 儿童	13.8-21.5cm	U1881ND-C0808				1 pcs	B1
Infant 婴儿	9-14.8cm	U1882ND-C0808				1 pcs	B1
Neonate 新生儿	5.8-10.8cm	U1883ND-C0808				1 pcs	B1
Thigh 大腿	45-56cm	U1884ND-C0808				1 pcs	B5
Large adult long 大成人加长	35.5-46cm	U1889ND-C1818		For GE Eagle Dash, Tram. Solar, MAC-LAB system		1 pcs	B4
Large adult 大成人	35.5-46cm	U1869ND-C1818					1 pcs
Adult long 成人加长	27.5-36.5cm	U1886ND-C1818				1 pcs	B2
Adult 成人	27.5-36.5cm	U1880ND-C1818				1 pcs	B1
Small adult 小成人	20.5-28.5cm	U1885ND-C1818				1 pcs	B1
Pediatric 儿童	13.8-21.5cm	U1881ND-C1818				1 pcs	B1
Infant 婴儿	9-14.8cm	U1882ND-C1818				1 pcs	B1
Neonate 新生儿	5.8-10.8cm	U1883ND-C1818				1 pcs	B1
Thigh 大腿	45-56cm	U1884ND-C1818				1 pcs	B5
Large adult long 大成人加长	35.5-46cm	U1889ND-C1718	For Datex Ohmeda			1 pcs	B4
Large adult 大成人	35.5-46cm	U1869ND-C1718					1 pcs
Adult long 成人加长	27.5-36.5cm	U1886ND-C1718				1 pcs	B2
Adult 成人	27.5-36.5cm	U1880ND-C1718				1 pcs	B1
Small adult 小成人	20.5-28.5cm	U1885ND-C1718				1 pcs	B1
Pediatric 儿童	13.8-21.5cm	U1881ND-C1718				1 pcs	B1
Infant 婴儿	9-14.8cm	U1882ND-C1718				1 pcs	B1
Neonate 新生儿	5.8-10.8cm	U1883ND-C1718				1 pcs	B1
Thigh 大腿	45-56cm	U1884ND-C1718				1 pcs	B5
Package 组合1: Adult 成人×2 Small adult 小成人×1 Adult long 成人加长×1 Large adult 大成人×1 Large adult long 大成人加长×1		U1801ND-C0808-06			Variety组合	6 pcs	C5
		U1801ND-C1818-06			Variety组合	6 pcs	C5
		U1801ND-C1718-06		Variety组合	6 pcs	C5	
Package 组合2: Adult 成人×1 Small adult 小成人×1 Pediatric 儿童×1 Adult long 成人加长×1 Large adult 大成人×1 Thigh 大腿×1		U1802ND-C0808-06		Variety组合	6 pcs	C6	
		U1802ND-C1818-06		Variety组合	6 pcs	C6	
		U1802ND-C1718-06		Variety组合	6 pcs	C6	

# 1.2 Reusable NIBP Cuff with inflation bag 重复用有囊血压袖套












1.2a NIBP Cuff without connector  
无接头血压袖套

Patient Size 适用病人	Limb Circumference 臂周长	REF# (Single tube) 型号 (单管)	REF # (Double tube) 型号 (双管)	Tubing Color 气管颜色	Package 包装	Unit Price 单价
Large adult long 大成人加长	33-47cm	U1889S	U1889D		1 pcs	B4
Large adult 大成人	33-47cm	U1869S	U1869D		1 pcs	B3
Adult long 成人加长	25-35cm	U1886S	U1886D		1 pcs	B1
Adult 成人	25-35cm	U1880S	U1880D		1 pcs	A8
Small adult 小成人	20-28cm	U1885S	U1885D		1 pcs	A8
Pediatric 儿童	18-26cm	U1881S	U1881D		1 pcs	A8
Infant 婴儿	10-19cm	U1882S	U1882D		1 pcs	A8
Neonate 新生儿	6-11cm	U1883S	U1883D		1 pcs	A8
Thigh 大腿	46-66cm	U1884S	U1884D		1 pcs	B5
Package 组合1: Adult 成人×2 small adult 小成人×1 Adult long 成人加长×1 Large adult 大成人×1 large adult long 大成人加长×1		U1801S-06 U1880S×2 U1885S×1 U1886S×1 U1869S×1 U1889S×1	U1801D-06 U1880D×2 U1885D×1 U1886D×1 U1869D×1 U1889D×1	Variety组合	6 pcs	C5
Package 组合2: Adult 成人×1 small adult 小成人×1 Child 儿童×1 Adult long 成人加长×1 Large adult 大成人×1 Thigh 大腿×1		U1802S-06 U1880S×1 U1885S×1 U1881S×1 U1886S×1 U1869S×1 U1884S×1	U1802D-06 U1880D×1 U1885D×1 U1881D×1 U1886D×1 U1869D×1 U1884D×1	Variety组合	6 pcs	C7

1.2b Single tube Cuff with connector  
带接头单管血压袖套

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Cuff end Connector 连接端	Tubing Color 气管颜色	Package 包装	Unit Price 单价
Large adult long 大成人加长	33-47cm	U1889S-C12	For Philips(HP) Siemens/Drager Datascope Spacelabs		1 pcs	B4
Large adult 大成人	33-47cm	U1869S-C12		1 pcs	B4	
Adult long 成人加长	25-35cm	U1886S-C12		1 pcs	B2	
Adult 成人	25-35cm	U1880S-C12		1 pcs	B1	
Small adult 小成人	20-28cm	U1885S-C12		1 pcs	B1	
Pediatric 儿童	18-26cm	U1881S-C12		1 pcs	B1	
Infant 婴儿	10-19cm	U1882S-C12		1 pcs	B1	
Neonate 新生儿	6-11cm	U1883S-C12		1 pcs	B1	
Thigh 大腿	46-66cm	U1884S-C12		1 pcs	B5	
Large adult long 大成人加长	33-47cm	U1889S-C03		For Colin Spacelabs Datascope		1 pcs
Large adult 大成人	33-47cm	U1869S-C03	1 pcs		B4	
Adult long 成人加长	25-35cm	U1886S-C03	1 pcs		B2	
Adult 成人	25-35cm	U1880S-C03	1 pcs		B1	
Small adult 小成人	20-28cm	U1885S-C03	1 pcs		B1	
Pediatric 儿童	18-26cm	U1881S-C03	1 pcs		B1	
Infant 婴儿	10-19cm	U1882S-C03	1 pcs		B1	
Neonate 新生儿	6-11cm	U1883S-C03	1 pcs		B1	
Thigh 大腿	46-66cm	U1884S-C03	1 pcs		B5	
Large adult long 大成人加长	33-47cm	U1889S-C18	For Protocol			1 pcs
Large adult 大成人	33-47cm	U1869S-C18		1 pcs	B4	
Adult long 成人加长	25-35cm	U1886S-C18		1 pcs	B2	
Adult 成人	25-35cm	U1880S-C18		1 pcs	B1	
Small adult 小成人	20-28cm	U1885S-C18		1 pcs	B1	
Pediatric 儿童	18-26cm	U1881S-C18		1 pcs	B1	
Infant 婴儿	10-19cm	U1882S-C18		1 pcs	B1	
Neonate 新生儿	6-11cm	U1883S-C18		1 pcs	B1	
Thigh 大腿	46-66cm	U1884S-C18		1 pcs	B5	
Package 组合1: Adult 成人×2 Small adult 小成人×1 Adult long 成人加长×1 Large adult 大成人×1 Large adult long 大成人加长×1		U1801S-C12-06			Variety组合	6 pcs
		U1801S-C03-06		Variety组合	6 pcs	C5
		U1801S-C18-06		Variety组合	6 pcs	C5
Package 组合2: Adult 成人×1 Small adult 小成人×1 Pediatric 儿童×1 Adult long 成人加长×1 Large adult 大成人×1 Thigh 大腿×1		U1802S-C12-06		Variety组合	6 pcs	C6
		U1802S-C03-06		Variety组合	6 pcs	C6
		U1802S-C18-06		Variety组合	6 pcs	C6

1.2c Double tube Cuff with particular connector  
带接头双管血压袖套

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Cuff end Connector 连接端	Tubing Color 气管颜色	Package 包装	Unit Price 单价	
Large adult long 大成人加长	33-47cm	U1889D-C0808			1 pcs	B4	
Large adult 大成人	33-47cm	U1869D-C0808			1 pcs	B4	
Adult long 成人加长	25-35cm	U1886D-C0808			1 pcs	B2	
Adult 成人	25-35cm	U1880D-C0808			1 pcs	B1	
Small adult 小成人	20-28cm	U1885D-C0808			1 pcs	B1	
Pediatric 儿童	18-26cm	U1881D-C0808			1 pcs	B1	
Infant 婴儿	10-19cm	U1882D-C0808			1 pcs	B1	
Neonate 新生儿	6-11cm	U1883D-C0808			1 pcs	B1	
Thigh 大腿	46-66cm	U1884D-C0808			1 pcs	B5	
Large adult long 大成人加长	33-47cm	U1889D-C1818				1 pcs	B4
Large adult 大成人	33-47cm	U1869D-C1818				1 pcs	B4
Adult long 成人加长	25-35cm	U1886D-C1818				1 pcs	B2
Adult 成人	25-35cm	U1880D-C1818				1 pcs	B1
Small adult 小成人	20-28cm	U1885D-C1818			1 pcs	B1	
Pediatric 儿童	18-26cm	U1881D-C1818			1 pcs	B1	
Infant 婴儿	10-19cm	U1882D-C1818			1 pcs	B1	
Neonate 新生儿	6-11cm	U1883D-C1818			1 pcs	B1	
Thigh 大腿	46-66cm	U1884D-C1818			1 pcs	B5	
Large adult long 大成人加长	33-47cm	U1889D-C1718				1 pcs	B4
Large adult 大成人	33-47cm	U1869D-C1718				1 pcs	B4
Adult long 成人加长	25-35cm	U1886D-C1718				1 pcs	B2
Adult 成人	25-35cm	U1880D-C1718				1 pcs	B1
Small adult 小成人	20-28cm	U1885D-C1718			1 pcs	B1	
Pediatric 儿童	18-26cm	U1881D-C1718			1 pcs	B1	
Infant 婴儿	10-19cm	U1882D-C1718			1 pcs	B1	
Neonate 新生儿	6-11cm	U1883D-C1718			1 pcs	B1	
Thigh 大腿	46-66cm	U1884D-C1718			1 pcs	B5	
Package 组合1: Adult 成人 × 2 Small adult 小成人 × 1 Adult long 成人加长 × 1 Large adult 大成人 × 1 Large adult long 大成人加长 × 1		U1801D-C0808-06			Variety组合	6 pcs	C5
		U1801D-C1818-06			Variety组合	6 pcs	C5
		U1801D-C1718-06			Variety组合	6 pcs	C5
Package 组合2: Adult 成人 × 1 Small adult 小成人 × 1 Pediatric 儿童 × 1 Adult long 成人加长 × 1 Large adult 大成人 × 1 Thigh 大腿 × 1		U1802D-C0808-06			Variety组合	6 pcs	C6
		U1802D-C1818-06		Variety组合	6 pcs	C6	
		U1802D-C1718-06		Variety组合	6 pcs	C6	




## ② Disposable NIBP Cuff 一次性无创血压袖套

### Specification 产品描述

- Non-woven or TPU material for maximum patient comfort
- Single & double tube cuff with 50 cm length hose (Neonate cuff with 20cm hose)
- Cuffs with particular connectors to fit major brands monitoring systems
- Large selection of 13 sizes cuffs matching the different sizes patients
- Printings in particular color indicate the cuff sizes.
- Latex free
- 舒适的无纺布或TPU材料供选择
- 单/双管血压袖套带50cm长血压气管 (除新生儿20cm外)
- 特定的袖套接头适用不同品牌的监护仪系列
- 13种尺寸大小满足不同大小病患者使用需求
- 特定的颜色表明袖套的尺寸大小
- 不含乳胶





### ③ Patient NIBP Adapter Air Hose

## 无创血压连接管

Original Maker & Models 原厂机型	Patient Size 适用病人	Plug to Monitor 机器端插头	Plug to Cuff 袖套端插头	Hose Length 连接管长度	REF# 产品型号	Original REF# 原装型号	Package 包装	Unit Price 单价
Philips All models except A1/A3	Adult Pediatric 成人/儿童			1.5m 2.5m 3.0m	HS15-16-15 HS-16-15 HS30-16-15	M1599B M1598B	1 pcs	B7
	Infant Neonate 婴儿/新生儿			1.5m 2.5m 3.0m	HSN15-16-20 HSN-16-20 HSN30-16-20	M1596B/C M1597B/C	1 pcs	C1
	Philips A1/A3 Nellcor, Siemens, Datascope, Spacelabs, Kontron, Mindray	Adult Pediatric 成人/儿童			2.5m 3.0m	HS-15-15 HS30-15-15	M3918A(Philips)	1 pcs
	Infant Neonate 婴儿/新生儿			2.5m	HSN-15-20		1 pcs	C1
	GE-Dinamap Pro, Compact, MPS Pro Care series	Adult Pediatric 成人/儿童			2.5m 3.6m 7.2m	HD-25-09 HD36-25-09 HD72-25-09	107365 107366	1 pcs
Infant Neonate 婴儿/新生儿				2.5m 3.6m	HDN-25-27 HDN36-25-27	107368	1 pcs	C4
Adult Pediatric 成人/儿童				2.5m 3.6m	HD-25-17 HD36-25-17	107363	1 pcs	C3
GE-Datex-Ohmeda AS/3, CS/3, S/5 Cardiocap 5/II, Light	Adult Pediatric 成人/儿童			1.8m 2.5m 3.0m 6.0m	HD18-22-17/18 HD-22-17/18 HD30-22-17/18 HD60-22-17/18	895732 877235 879739	1 pcs	C3
	Infant Neonate 婴儿/新生儿			2.5m 3.0m	HDN-22-27 HDN30-22-27	9461-205 877514	1 pcs	C4
	Infant Neonate 婴儿/新生儿			3.6m	HDN36-22-20	9461-215	1 pcs	C3
GE-Marquette Eagle4000, Solar Tram100/200/300	Adult Pediatric 成人/儿童			2.5 3.6m	HD-22-17 HD36-22-17	9461-203	1 pcs	C3
	Adult Pediatric 成人/儿童			2.5m 3.6m	HD-22-09 HD36-22-09	9461-214	1 pcs	C3
	Infant Neonate 婴儿/新生儿			2.5m 3.0m	HDN-22-27 HDN30-22-27	9461-205 877514	1 pcs	C4
GE-Marquette -Dash, Eagle, SOLAR Tram400/450/451 800/850/851, Coro 120	Adult Pediatric 成人/儿童			2.5m 3.6m	HD-24-17 HD36-24-17	414873-001	1 pcs	C3
	Adult Pediatric 成人/儿童			2.5m 3.6m	HD-24-09 HD36-24-09	9461-217	1 pcs	C3
	Infant Neonate 婴儿/新生儿			2.5m	HDN-24-27	414874-001	1 pcs	C4
	Infant Neonate 婴儿/新生儿			3.6m	HDN36-24-S20	9461-216	1 pcs	C3
	Infant Neonate 婴儿/新生儿			3.6m	HDN36-24-20	9461-218	1 pcs	C3



Original Maker & Models 原厂机型	Patient Size 适用病人	Plug to Monitor 机器端插头	Plug to Cuff 袖套端插头	Hose Length 连接管长度	REF# 产品型号	Original REF# 原装型号	Package 包装	Unit Price 单价
GE-Critikon Invivo, Kontron, Nihon Kohden	Adult Pediatric 成人/儿童			2.5m 3.6m 7.2m	HD-09-09 HD36-09-09 HD72-09-09	008841 08842	1 pcs	C3
	Adult Pediatric 成人/儿童			3.6m	HD36-09-17	088845	1 pcs	C3
	Infant Neonate 婴儿/新生儿			2.5m	HDN-09-27	008840	1 pcs	C4
Nihon Kohden BSM/PVB	Adult Pediatric 成人/儿童			1.5m 2.5m 3.5m	HD15-30-18 HD-30-18 HD35-30-18	YN-900P YN-901P	1 pcs	C6
	Infant Neonate 婴儿/新生儿			1.5m 2.5m 3.5m	HDN15-30-20 HDN-30-20 HDN35-30-20	YN-920P YN-921P	1 pcs	C7
	Adult Pediatric 成人/儿童			1.5m 2.5m 3.0m 3.5m	HD15-37-18 HD-37-18 HD30-37-18 HD35-37-18	YN-100T YN-101T YN-101T	1 pcs	C6
Life Scope N OPV-1500 Life Scope N OPV-1500K	Adult Pediatric 成人/儿童			1.5m 2.5m 3.0m 3.5m	HD15-37-18 HD-37-18 HD30-37-18 HD35-37-18	YN-100T YN-101T YN-101T	1 pcs	C6
	Adult Pediatric 成人/儿童			2.5m	HS-04-07		1 pcs	B7
Criticare, Datascope, Spacelabs, Colin, Welch Allyn, Nihon Kohden	Adult Pediatric 成人/儿童			2.5m	HS-04-07		1 pcs	B7
	Infant Neonate 婴儿/新生儿			2.5m	HSN-04-20		1 pcs	C1
Colin	Adult Pediatric 成人/儿童			2.5m 3.5m	HS-31-32 HS35-31-32	A015ZZ	1 pcs	C1
	Infant Neonate 婴儿/新生儿			2.5m 3.5m	HSN-31-20 HSN35-32-20	AL012Z	1 pcs	C1
Drager/Siemens SC6000- 7000/8000/ 9000 Infinity	Adult Pediatric 成人/儿童			2.5m 3.7m	HS-40-15 HS37-40-15	MP00953 (1275 275 EH50u)	1 pcs	C6
	Infant Neonate 婴儿/新生儿			2.4m 2.5m	HSN24-40-20 HSN-40-20	2870298	1 pcs	C7
Physio Control Lifepack 20	Adult Pediatric 成人/儿童			2.5m	HS-12-15		1 pcs	B7
	Infant Neonate 婴儿/新生儿			2.5m	HSN-12-20		1 pcs	C1
MEK MP700, MP1000	Adult Pediatric 成人/儿童			2.5m	HS-09-09		1 pcs	B7
	Infant Neonate 婴儿/新生儿			2.5m	HSN-09-20		1 pcs	C1

Original Maker & Models 原厂机型	Patient Size 适用病人	Plug to Monitor 机器端插头	Plug to Cuff 袖套端插头	Hose Length 连接管长度	REF# 产品型号	Original REF# 原装型号	Package 包装	Unit Price 单价
Welch Allyn Atlas	Adult Pediatric 成人/儿童			2.5m	HS-04-09		1 pcs	B7
	Infant Neonate 婴儿/新生儿			2.5m	HSN-04-20		1 pcs	C1
	Adult Pediatric 成人/儿童			2.5m	HD-41-09		1 pcs	C3
	Infant Neonate 婴儿/新生儿			2.5m	HDN-41-27		1 pcs	C3
Infinium OMNI/OMNI EXPRESS/ OMNI EXPRESS VS OMNI III/OMNI III	Adult Pediatric 成人/儿童			2.1m	HS-36-2T		1 pcs	B7






## 4.1 Universal Y Adapter Y型接头



4.1a 1-tube air hose to 2-tube BP cuff, 10pcs/pack  
单管转双管袖套Y型接头, 10个/包




## 4 NIBP Connectors & tubing 无创血压接头与血压管



REF#型号	BP3T-Y	BP3T-1718	BP3T-1717	BP3T-0909	BP3T-2T2T	BP3T-2020
Image图片						
Price单价	B7	C6	C6	D1	C2	D1
REF#型号	BP12-3T	BP12-1718	BP12-1717	BP12-0909	BP12-2T2T	BP12-2020
Image图片						
Price单价	C2	C8	C8	D3	C6	D3
REF#型号	BP18-3T	BP18-1718	BP18-1717	BP18-0909	BP18-2T2T	BP18-2020
Image图片						
Price单价	C2	C8	C8	D3	C6	D3
REF#型号	BP04-3T	BP04-1718	BP04-1717	BP04-0909	BP04-2T2T	BP04-2020
Image图片						
Price单价	C8	D3	D3	D4	D1	D4
REF#型号	BP08-3T	BP08-1718	BP08-1717	BP08-0909	BP08-2T2T	BP08-2020
Image图片						
Price单价	C2	C8	C8	D3	C6	D3
REF#型号	BP2T-3T	BP2T-1718	BP2T-1717	BP2T-0909	BP2T-2T2T	BP2T-2020
Image图片						
Price单价	C2	C8	C8	D3	C6	D3
REF#型号	BP21-3T	BP21-1718	BP21-1717	BP21-0909	BP21-2T2T	BP21-2020
Image图片						
Price单价	C2	C8	C8	D3	C6	D3

4.1b 2-tube air hose to 1-tube BP cuff, 10 pcs/pack  
 双管转单管Y型接头, 10个/包



REF#型号	BP3T-Y	BP15-3T	BP04-3T	BP18-3T	BP2T-3T	BP20-3T
Image图片						
Price单价	B7	D1	C8	C2	C2	C6
REF#型号	BP3T-1718	BP15-1718	BP04-1718	BP18-1718	BP2T-1718	BP20-1718
Image图片						
Price单价	C6	D4	D3	C8	C8	D1
REF#型号	BP3T-1818	BP15-1818	BP04-1818	BP18-1818	BP2T-1818	BP20-1818
Image图片						
Price单价	C6	D4	D3	C8	C8	D1
REF#型号	BP3T-0808	BP15-0808	BP04-0808	P18-0808	BP2T-0808	BP20-0808
Image图片						
Price单价	C2	D3	D1	C6	C6	C8
REF#型号	BP3T-1212	BP15-1212	BP04-1212	BP18-1212	BP2T-1212	BP20-1212
Image图片						
Price单价	C6	D4	D3	C8	C8	D1
REF#型号	BP3T-0505	BP15-0505	BP04-0505	BP18-0505	BP2T-0505	BP20-0505
Image图片						
Price单价	D1	D4	D3	D3	D3	D3

## 4.2 Connectors for Cuffs and Hoses 袖套及连接管接头



Unit: Per pcs

REF#型号	BP03	BP04	BP05	BP06	BP07	BP08
Fig. 图片						
Match	(BP04)>> BP06	(BP03)>>BP07	>>BP20	(BP07)>>BP03	(BP06)>>BP04	>> BP09
Comp 兼容性	Mennen (Old) Schiller Bruker/Odam Spacelabs(Old)	Same as left	CAS, Drager Kontron(Amphora) Physio control(Old) Schiller Bruker/Odam	BCI Monitor (Old) Datascope(Old)	Same as left	Welch Allyn, GE Critikon(2) Invivo(2), Kontron(2) Mennen(2) Nihon Kohden(2)
Price单价	A1	A5	A3	A1	A5	A1
REF#型号	BP09	BP12	BP12-P	BP15	BP17	BP18
Fig. 图片						
Match	>>BP08	(BP12-P)>>BP15	(BP12)>>BP15	>>BP12	>>BP18	>>BP17
Comp 兼容性	GE-Critikon Invivo, Kontron, Nihon Kohden	Criticare, Datascope Drager/Siemens Physio control Mindray, Philips/HP Spacelabs, Nellcor	Same as BP12	Schiller Air hose connector	Nihon Kohden Datex(2)	Zoll Series M Datex(2) Marquette(2)
Price单价	A3	A2	A1	A6	A2	A2
REF#型号	BP20	BP21	BP2T-M	BP34(BP31)	BP35(BP32)	
Fig. 图片						
Match	>>BP05	>>BS21	Two end thread	>.BP35	Monitor plug	
Comp 兼容性	Regularly as the connector of air hose for neonate disposable BP cuffs	Pace Tech	General purpose	Colin Gambro Tyco( NPB 5500)	>>BP31	
Price单价	A3	A2	A2	A5	A5	

## 4.3 Air Tubing 血压管

Figures图片	REF# 型号	Application 用途	Inner Diameter 内径	Color 颜色	Price / m 价格 / 米
	Y02-H003	Adult/Pediatric single tube 成人/儿童单管	ø3.50mm	Light gray 浅灰色	A3
	Y02-H001	Adult/Pediatric double tube 成人/儿童双管	ø3.50mm	Light gray 浅灰色	A4
	Y02-H004	Neonate/ Infant single tube 新生儿/婴儿单管	ø2.70mm	Light blue 浅蓝色	A4
	Y02-H002	Neonate/ Infant double tube 新生儿/婴儿双管	ø2.70mm	Light blue 浅蓝色	A6



## 5 Manual Inflation System 手动充气装置



### 5.1 Reusable Cuff with Inflation System (2-tube reusable cuff, with bulb&valve)

见图A 充气装置带重复用血压袖套 (双管重复用袖套, 带充气球和阀门)

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Package 包装	Package Price 包装价格
Large adult long 大成人加长	33-47cm	U1889ND-B01	5 pcs/box	D3
Large adult 大成人	33-47cm	U1869ND-B01	5 pcs/box	D2
Adult long 成人加长	25-35cm	U1886ND-B01	5 pcs/box	D1
Adultb 成人	25-35cm	U1880ND-B01	5 pcs/box	C8
Small adult 小成人	20-28cm	U1885ND-B01	5 pcs/box	C8
Pediatric 儿童	18-26cm	U1881ND-B01	5 pcs/box	C8

### 5.2 Disposable cuff with Inflation System (2-tube disposable cuff, with bulb&valve)

见图B 充气装置带一次性血压袖套 (双管一次性血压袖套, 带充气球和阀门)

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Package 包装	Package Price 包装价格
Large adult long 大成人加长	33-47cm	U1710D-B02	5 pcs/box	C8
Large adult 大成人	33-47cm	U1720D-B02	5 pcs/box	C7
Adult long 成人加长	25-35cm	U1730D-B02	5 pcs/box	C6
Adultb 成人	25-35cm	U1740D-B02	5 pcs/box	C4
Small adult 小成人	20-28cm	U1750D-B02	5 pcs/box	C4
Pediatric 儿童	18-26cm	U1760D-B02	5 pcs/box	C4

# 6 Wall Manometer Accessories

## 挂壁式血压计附件



A

B

### 6.1 Inflation system with reusable cuff

(Coiled tubing, Adapter with screw connectors, bulb & Valve, and 2-tube reusable cuff)

#### 带重复用血压袖套充气装置

(弹簧管, 带螺纹接头及充气球&阀门, 重复用双管血压袖套)

Patient Size 适用病人	Limb Circumference 臂周长	REF# 型号	Package 包装	Package Price 包装价格
Adult 成人	25-35cm (with Clip,带环)	U1880HD-CT-B01	5 pcs/box	D4
Large adult long 大成人加长	35.5-46cm	U1889ND-CT-B01	5 pcs/box	D5
Large adult 大成人	35.5-46cm	U1869ND-CT-B01	5 pcs/box	D5
Adult long 成人加长	27.5-36.5cm	U1886ND-CT-B01	5 pcs/box	D4
Adult 成人	27.5-36.5cm	U1880ND-CT-B01	5 pcs/box	D4
Small adult 小成人	20.5-28.5cm	U1885ND-CT-B01	5 pcs/box	D4
Pediatric 儿童	13.8-21.5cm	U1881ND-CT-B01	5 pcs/box	D4

### 6.2 Inflation system without cuff

(Coiled tubing, w or w/o connectors and bulb & Valve)

#### 不带袖套的充气装置

(弹簧管, 带或不带接头及充气球&阀门)



Used with 2-tube Cuff w/o Conn. 用于不带接头的双管袖套	Used with 2-tube Cuff W/BP09 Conn. 用于带BP09接头的双管袖套	Used with 1-tube Cuff w/o Conn. 用于不带接头的单管袖套	Used with 2-tube Cuff W/BP15 Conn. 用于带BP15接头的双管袖套	Package 包装	Package Price 包装价格
REF# SP/S-CD-B01	REF# SP/S-C0909-B01	REF# SP/S-CS-B01	REF# SP/S-C15-B01	10 pcs/box	D5

### 6.3 Parts

#### 部件

Figures 图片	REF# 型号	Description 产品描述	Package 包装	Package Price 包装价格
	U1880HS U1880HD	Reusable adult (25-35cm) 1 or 2-tube BP cuff with clip 带环单或双管成人重复用袖套, 臂周长 (25-35厘米)	1 pcs/box	B6
	Y02-H005	Coiled black single tube, 3.8m, 4.00mm 黑色弹簧单管, 3.8米, 内径4毫米	1 pcs/box	A7
	Y15-S001	Pocket Aneroid gauge, Reusable, gray 掌上无液压力表, 可重复使用, 灰色	1 pcs/box	A7
	Y14-L001	Bulb & Valve, reusable, black 充气球&阀门, 可重复使用, 黑色	5 pcs/box	B4

# 7 Tourniquet

## 止血带

### Specification

#### 产品描述

- Strong gentle nylon materials for thousands circles of inflation and deflation.
- Reusable for multiple patient use
- Convenient and easy to clean
- Come with 50cm PVC tube
- Latex free
- 复合尼龙材料可经受数千次循环充气
- 可多病人重复使用
- 方便且容易清洗
- 带50cm PVC气管
- 不含乳胶



Size Type 尺寸类型		Cuff Sizes(mm) 尺寸 (mm)	REF # (1-tube) 型号 (单管)	REF # (2-tube) 型号 (双管)	Package 包装	Package Price 包装价格
Small	小号	50X400	T1880S	T1880D	1 pcs	A8
Medium	中号	70X500	T1881S	T1881D	1 pcs	B1
Large	大号	80X800	T1882S	T1882D	1 pcs	B2
Variety(each of the above item) 组合 (含以上型号各1个)		See above 见以上尺寸	T1801S-03	T1801D-03	3 pcs	B8

# Original Brands Medical Devices

## 代理/经销产品

### Pronk Technologies (USA) — Patient Multi-Parameter Simulator病人多参数模拟器系列



Oximetry Simulator  
血氧模拟器



Multi-Parameter Simulator  
多参数模拟器



Non-Invasive BP Simulator  
无创血压模拟器



Multi-Parameter Simulator set  
多参数模拟器套件

### Welchallyn 伟伦 (USA)



喉镜



检眼镜/检耳镜



头灯

### Pulse Oximeter 脉搏血氧仪



Tuffsat Pulse Oximeter  
手掌式脉搏血氧仪  
Manufacturer制造商: GE Healthcare



N-65 Oximax Pulse Oximeter  
手掌式脉搏血氧仪  
Manufacturer制造商: Covidien



N-560 Oximax Pulse Oximeter  
台式脉搏血氧仪  
Manufacturer制造商: Covidien

### Original monitoring accessories 原装监护附件



Philips Medical



GE Healthcare



Covidien



Masimo

### Anesthesia & Respiratory Supplies 麻醉呼吸耗材



Anaesthesia & Ventilation Circuit  
麻醉呼吸回路



Resuscitator  
人工呼吸球



Humidifier  
湿化器

### Compatible Fetal Transducer(USA) 兼容性胎监传感器



US/TOCO Transducer  
胎心/宫缩探头  
Compatible with major brands  
兼容主要品牌厂家

### Medical Oxygen Sensor(Germany) 医用氧电池



Oxygen Sensor & Analyzer  
氧电池及氧浓度测试仪  
Compatible with major brands  
兼容主要品牌厂家

• For detail information of the above devices like specifications, prices and availabilities, please contact Unimed sales representatives.  
了解以上产品详情, 请联系尤迈医疗销售代表

# Content

## 内容

1. IBP Adapter Cable for transducer 有创血压电缆 ..... 113
2. IBP Convert Cable 有创电缆转接线 ..... 117
3. IBP Transducer adapter 传感器转接线 ..... 119
4. IBP Disposable Transducer 一次性有创压传感器 ..... 120

# IBP Interface Cable & Transducer

## 有创血压电缆及传感器

### Certification 产品认证




CE Mark  
CE 标志



FDA 510(k) clearance No.: K082959  
Patient Cable & Leadwire / 病人电缆及导联线

### Compliance with 符合标准

ANSI/AAMI EC53-1998(R)2001  
EN 60601-1:1990  
ISO10993-1:2009  
ISO10993-5:2009  
ISO10993-10:2010



# 1 IBP Adapter cable for Transducer

## 有创血压电缆

Compatible with major patient monitor brands  
兼容主要品牌


- Datex-Ohmeda 德恩-欧美达
- Datascope-Mindray 迈瑞
- Marquette 马奎
- Nihon kohden 日本光电
- Philips 飞利浦
- Spacelabs 太空
- Siemens 西门子
- ...

Specification  
产品规格

Cable material 电缆物料	TPU Jacket, Gray TUP料, 灰色
Connector 连接器	Nylon plug & Golden plated pins 尼龙插座和镀金针
Cable length 电缆长度	4.0m/0.3m
Cable Diameter 电缆直径	5.00mm

Available cables for disposable transducer  
用于一次性传感器的电缆

- Argon (AG)  
阿尔贡
- B.Braun (BB)  
贝朗
- B.D. (BD)  
贝达
- Edwards (ED)  
爱德华
- Medex/Abbort (MX)  
美代/雅培
- Utah (UT)  
尤它
- PVB



Reference Number Rule Description  
产品型号说明

BC-HP-UT/2	
BC	IBP cable 有创压电缆
HP	Factory/model (see others above) 厂家/型号 (其他, 看如上所示)
UT	Utah( transducer type) 尤它 (传感器类型),其他型号见如上示
2	the number of channel, default is single channel 通道数; 不写通道数, 则默认为单通道

Available cables for reusable transducer  
用于可重复使用传感器的电缆





Transducer Type 传感器类型	Argon 阿尔贡	B.Braun 贝朗	BD 贝达	Edward 爱德华	Medex / Abbott 美代/雅培	Utah 尤它	Medex Logical 美代	PVB	Cable Length 电缆长度	Cable Connector on Monitor end 仪器端电缆接头
<b>Cable Connector on transducer end 传感器端电缆接头</b>										
<b>Criticare</b>										
Model 2200	BC-CSI-AG C3	BC-CSI-BB C3	BC-CSI-BD C8	BC-CSI-ED C3	BC-CSI-MX C3	BC-CSI-UT C3	BC-CSI-MX1 D1	BC-CSI-PVB C3	4.0m	
<b>Colin</b>										
	BC-CL-AG C3	BC-CL-BB C3	BC-CL-BD C8	BC-CL-ED C3	BC-CL-MX C3	BC-CL-UT C3	BC-CL-MX1 D1	BC-CL-PVB C3	4.0m	
<b>Biolight</b>										
V6 Vital Sign View A2/A3/A5/A6	BC-BLT-AG C3	BC-BLT-BB C3	BC-BLT-BD C8	BC-BLT-ED C3	BC-BLT-MX C3	BC-BLT-UT C3	BC-BLT-MX1 D1	BC-BLT-PVB C3	4.0m	
<b>GE-Datex</b>										
Airshield, Burdick, Kone	BC-DX-AG C3	BC-DX-BB C3	BC-DX-BD C8	BC-DX-ED C3	BC-DX-MX C3	BC-DX-UT C3	BC-DX-MX1 D1	BC-DX-PVB C3	4.0m	
<b>GE-Hellige</b>										
	BC-HL-AG C3	BC-HL-BB C3	BC-HL-BD C8	BC-HL-ED C3	BC-HL-MX C3	BC-HL-UT C3	BC-HL-MX1 D1	BC-HL-PVB C3	4.0m	
<b>GE-Marquette</b>										
Eagle/dash1000/2000 3000/3100/4000Solar 7000/8000/9000/Tram Module 100-600, Cath-Lab, MAC-Lab Dash †/Tram(all versions) E-PSMP/E-PRESTN/PDM	BC-MQ-AG C3	BC-MQ-BB 22339301 C3	BC-MQ-BD 2016995-001 C8	BC-MQ-ED 2021197-001 C3	BC-MQ-MX 2021196-001 C3	BC-MQ-UT 2016997-001 C3	BC-MQ-MX1 D1	BC-MQ-PVB C3	4.0m	
	BC-MQ-AG/2 C8	BC-MQ-BB/2 C8	BC-MQ-BD/2 2016995-003 D3	BC-MQ-ED/2 2021197-003 C8	BC-MQ-MX/2 2021196-003 C8	BC-MQ-UT/2 2016997-003 C8	BC-MQ-MX1/2 D4	BC-MQ-PVB/2 C8		
<b>Mostcare</b>										
	BC-MC-AG C3	BC-MC-BB C3	BC-MC-BD C8	BC-MC-ED C3	BC-MC-MX C3	BC-MC-UT C3	BC-MC-MX1 D1	BC-MC-PVB C3	4.0m	
<b>Mennen</b>										
Envoy- vitalogik 5500	BC-MN-AG C3	BC-MN-BB C3	BC-MN-BD C8	BC-MN-ED C3	BC-MN-MX C3	BC-MN-UT C3	BC-MN-MX1 D1	BC-MN-PVB C3	4.0m	
<b>Mindray/Datascope</b>										
Beneview T5/T6/T8 DPM 6, DPM 7	BC-MR-AG C3	BC-MR-BB C3	BC-MR-BD 040-000053-00 C8	BC-MR-ED 040-000054-00 C3	BC-MR-MX 040-000052-00 C3	BC-MR-UT C3	BC-MR-MX1 D1	BC-MR-PVB C3	4.0m	
<b>Datascope</b> Passport V, V12, V21	BC-DT1-AG C3	BC-DT1-BB C3	BC-DT1-BD 001C-30-70757 C8	BC-DT1-ED 040-000053-00 C3	BC-DT1-MX 040-000054-00 C3	BC-DT1-UT 040-000052-00 C3	BC-DT1-MX1 D1	BC-DT1-PVB C3	4.0m	
<b>Philips/HP</b>										
78352A, M1165A, MerlinViridia 24, CMS 24,M3/M4,M1175A, M1350B,M1205A, M1350A,M1353A, 8030A,8040A	BC-HP-AG C3	BC-HP-BB M1634A C3	BC-HP-BD C8	BC-HP-ED C3	BC-HP-MX C3	BC-HP-UT C3	BC-HP-MX1 D1	BC-HP-PVB C3	4.0m	
<b>Siemens</b>										
Sirecust 400-1281, SC8000	BC-SM1-AG C3	BC-SM1-BB C3	BC-SM1-BD C8	BC-SM1-ED C3	BC-SM1-MX C3	BC-SM1-UT C3	BC-SM1-MX1 D1	BC-SM1-PVB C3	4.0m	
	BC-SM1-AG/2 C8	BC-SM1-BB/2 C8	BC-SM1-BD/2 D3	BC-SM1-ED/2 C8	BC-SM1-MX/2 C8	BC-SM1-UT/2 C8	BC-SM1-MX1/2 D4	BC-SM1-PVB/2 C8		
SC 6000,7000,9000	BC-SM2-AG C3	BC-SM2-BB C3	BC-SM2-BD C8	BC-SM2-ED C3	BC-SM2-MX C3	BC-SM2-UT C3	BC-SM2-MX1 D1	BC-SM2-PVB C3	4.0m	
	BC-SM2-AG/2 C8	BC-SM2-BB/2 C8	BC-SM2-BD/2 D3	BC-SM2-ED/2 C8	BC-SM2-MX/2 C8	BC-SM2-UT/2 C8	BC-SM2-MX1/2 D4	BC-SM2-PVB/2 C8		

Transducer Type 传感器类型	Argon 阿尔贡	B.Braun 贝朗	BD 贝达	Edward 爱德华	Medex / Abbott 美代/雅培	Utah 尤它	Medex Logical 美代	PVB	Cable Length 电缆长度	Cable Connector on Monitor end 仪器端电缆接头
<b>Cable Connector on transducer end 传感器端电缆接头</b>										
<b>Drager</b>										
Cato, Cicero EM, Vitara Parameterbox PB 8800; PM8010; PM8014; PM 8040;PM 8060, Sulla, UM3, UM3.1	BC-DG-AG C3	BC-DG-BB C3	BC-DG-BD C8	BC-DG-ED C3	BC-DG-MX C3	BC-DG-UT C3	BC-DG-MX1 D1	BC-DG-PVB C3	4.0m	
Infinity, Gamma X XI Delta/delta XI Kappa/kappa XI Vista XI	BC-S-AG/2 D5	BC-S-BB/2 D5	BC-S-BD/2 D5	BC-S-ED/2 D5	BC-S-MX/2 D5	BC-S-UT/2 D5	BC-S-MX1/2 D6	BC-S-PVB/2 D5	4.0m	
<b>Datascope</b>										
System 80/82/83/84 90/90T/95/98/97IABP 850/870/P2/P3/2000 2000A/2001/2001A 2002/2002A/2200I 3000/Passport	BC-DT-AG C3	BC-DT-BB C3	BC-DT-BD 0012-00-1245 C8	BC-DT-ED C3	BC-DT-MX C3	BC-DT-UT C3	BC-DT-MX1 D1	BC-DT-PVB C3	4.0m	
<b>Fukuda Denshi</b>										
DYNASCOPE/Expert DS5100E/DS5300W DS5400	BC-FD-AG C3	BC-FD-BB C3	BC-FD-BD C8	BC-FD-ED C3	BC-FD-MX C3	BC-FD-UT C3	BC-FD-MX1 D1	BC-FD-PVB C3	4.0m	
<b>Nihon Kohden</b>										
Lifescop 6/9/12/14 AP-311P -800PA BMS 2100-8500AOEC 3200-8120A/OMP-3221E -n7203/Mercury 9510K	BC-NK1-AG C3	BC-NK1-BB C3	BC-NK1-BD C8	BC-NK1-ED C3	BC-NK1-MX C3	BC-NK1-UT C3	BC-NK1-MX1 D1	BC-NK1-PVB C3	4.0m	
New	BC-NK2-AG C8	BC-NK2-BB C8	BC-NK2-BD JP-900P D3	BC-NK2-ED JP-920P/3.5m JP-921P/1.5m C8	BC-NK2-MX C8	BC-NK2-UT JP-960P C8	BC-NK2-MX1 D1	BC-NK2-PVB C8	4.0m	
<b>Kontron</b>										
Minimon 7131-7136 Fetalmon 2A 2BFetalogic 4000 Fetasonde 5/5A/5B K2000 Ballon Pump KAAT II/Module 10-7286 Supermon	BC-KTN-AG C3	BC-KTN-BB C3	BC-KTN-BD C8	BC-KTN-ED C3	BC-KTN-MX C3	BC-KTN-UT C3	BC-KTN-MX1 D1	BC-KTN-PVB C3	4.0m	
<b>Spacelabs</b>										
Ultraview 90496, 91496	BC-SL-AG C5	BC-SL-BB C5	BC-SL-BD D1	BC-SL-ED 700-0293-00 C5	BC-SL-MX 700-0295-00 C5	BC-SL-UT C5	BC-SL-MX1 D1	BC-SL-PVB C5	4.0m	
	BC-SL-AG/2 D1	BC-SL-BB/2 D1	BC-SL-BD/2 D3	BC-SL-ED/2 700-0292-00 D1	BC-SL-MX/2 700-0294-00 D1	BC-SL-UT/2 D1	BC-SL-MX1/2 D4	BC-SL-PVB/2 D1		
<b>Spacelabs &amp; General</b>										
BCI, IVY,Citikon, Invivo,Mennen,NDE, PB,Nellcor,Tektronix, Viatak, Mindray	BC-6P-AG B8	BC-6P-BB B8	BC-6P-BD C8	BC-6P-ED B8	BC-6P-MX B8	BC-6P-UT B8	BC-6P-MX1 D1	BC-6P-PVB B8	4.0m	
<b>Stockert</b>										
	BC-SK-AG C8	BC-SK-BB C8	BC-SK-BD D3	BC-SK-ED C8	BC-SK-MX C8	BC-6P-UT C8	BC-SK-MX1 D1	BC-SK-PVB C8	4.0m	

# 2 IBP Convert Cable

## 有创电缆转接线

### Reference Number Rule Description

#### 产品型号说明

BC-SL-AD/2	
BC	IBP cable 有创血压电缆
SL	Factory/model (see others on previous page) 厂家/型号 (其他参考前页)
AD	Convert Cable.AD: Converted to AAMI Adapter(MQ: to marquette;SM1 or SM2: to Siemens) 转接线 AD: 连接6针通用有创电缆 (MQ: 连接marquette电缆; SM1或SM2: 连接西门子电缆)
2	2-channel,default is single channel 2-通道数; 不写通道数, 则默认为单通道



Cable end Connector 电缆端接头	Unimed REF# 尤迈型号	OEM PN 原厂型号	Used with 连接电缆	Reference 备注	Cable Length 电缆长度	Monitor End Connector 仪器端接头	Price code 价格代码
GE Marquette							
	BC-MQ-AD/2		6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous for Marquette 参考前页面展示的马葵型号	0.3m		D3
	BC-MQ-MQ/2	2005772-001	Marquette IBP adapter Cable 马葵有创血压转接线	See the previous for Marquette 参考前页面展示的马葵型号	0.3m		D5
Spacelabs							
	BC-SL-AD	700-0791-00	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Spacelabs 参考前页面展示的太空型号	0.3m		D1
	BC-SL-AD/2	700-0028-00	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Spacelabs 参考前页面展示的太空型号	0.3m		D5
Datascope/Mindray							
	BC-MR-AD	0010-20-42795	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Mindray. 参考前页面展示的心迈型号	0.3m		C3
	BC-DT1-AD	040-000096-00	6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Mindray 参考前页面展示的心迈型号	0.3m		C3
	BC-DT1-DT		6 pins general IBP adapter cable 6针通用有创血压转接线	See the previous page for Mindray & Datascope 参考前页面展示的心迈 & 达施康型号	0.3m		C3
Siemens/Drager							
	BC-S-SM1/2	5731281	SM1 conn. 2XIBP Cable SM1接头 2根有创血压电缆	See the previous Page for Siemens 参考前页面展示的西门子型号	0.3m		D6
	BC-S-SM2/2	5592147	SM2 conn. 2XIBP Cable SM2接头 2根有创血压电缆	See the previous Page for Siemens 参考前页面展示的西门子型号	0.3m		D6
	BC-SM2-SM1/2	5588095 Drager Gamma XL and Delta XL	SM1 conn. 2XIBP Cable SM1接头 2根有创血压电缆	See the previous Page for Siemens 参考前页面展示的西门子型号	0.3m		D5
	BC-SM1-SM2		SM2 conn. IBP Cable SM2接头 有创血压电缆	See the previous Page for Siemens 参考前页面展示的西门子型号	0.3m		C3
	BC-SM2-SM1	3368383	SM1 conn. IBP Cable SM1接头 有创血压电缆	See the previous Page for Siemens 参考前页面展示的西门子型号	0.3m		C3

### ③ IBP Transducer Adapter 有创血压传感器转接线




Reference Number Rule Description  
产品型号说明

BC-MX-AG
BC IBP cable有创血压电缆
MX Factory/model, used with Medex adaptor Cable厂家/型号, 与美代有创血压电缆配用。
AG Factory/model, used with Argon disposable transducer厂家/型号, 与阿尔贡一次性传感器配用。 (Factory/model see others above)(具体厂家/型号, 参考前面页的描述)

Used with brands IBP Adapter cable & Major types disposable transducer  
连接各品牌IBP电缆及主要一次性有创压传感器

Transducer Type 传感器类型	Argon 阿尔贡	B.Braun 贝朗	BD 贝达	Edward 爱德华	Medex / Abbott 美代/雅培	Utah 尤它	Medex Logical 美代	Cable Length 电缆长度
With transducer conn. 与传感器接头	NA		NA	NA			NA	
With IBP cable conn. 与电缆接头								
		BC-BB-AG C3	BC-BD-AG C5	BC-ED-AG C3	BC-MX-AG C3	BC-UT-AG C3	BC-ML-AG C3	0.3m
B.Braun								
	BC-AG-BB C3		BC-BD-BB C5	BC-ED-BB C3	BC-MX-BB C3	BC-UT-BB C3	BC-ML-BB C3	0.3m
BD								
	BC-AG-BD C3	BC-BB-BD C3		BC-ED-BD C3	BC-MX-BD C3	BC-UT-BD C3	BC-ML-BD C3	0.3m
Edward								
	BC-AG-ED C3	BC-BB-ED C3	BC-BD-ED C5		BC-MX-ED C3	BC-UT-ED C3	BC-ML-ED C3	0.3m
Medex/Abbott								
	BC-AG-MX C3	BC-BB-MX C3	BC-BD-MX C5	BC-ED-MX C3		BC-UT-MX (97005-C) C3	BC-ML-MX C3	0.3m
Utah								
	BC-AG-UT C3	BC-BB-UT C3	BC-BD-UT C5	BC-ED-UT C3	BC-MX-UT C3		BC-ML-UT C3	0.3m
Medex Logical								
	BC-AG-MX1 C3	BC-BB-MX1 C3	BC-BD-MX1 C8	BC-ED-MX1 C3	BC-MX-MX1 C3	BC-UT-MX1 C3		0.3m



• Please refer to the additional Transducer catalogue  
请参考独立的单页目录册

### ④ Disposable Transducer 一次性有创压传感器



UTAH

MEDEX

B.BRAUN

BD

EDWARD

# Content

## 内容

<b>1. Reusable Temp. Probe 重复用体温探头</b> .....	<b>123</b>
1.1 Reusable Temp. Probe 重复用体温探头 .....	124
1.2 Temp. Adapter Cable for reusable probe 重复用体温探头连接电缆 .....	127
<b>2. Disposable Temp. Probe 一次性体温传感器</b> .....	<b>129</b>
2.1 Disposable Skin, Esophageal/Rectal Temp. Probe 一次性体表, 体腔体温传感器 .....	130
2.2 Disposable Tympanic Temp. Probe 一次性耳鼓体温传感器 .....	131
2.3 Foley Catheter Temp. Probe 一次性导尿管体温传感器 .....	131
2.4 Disposable Esophageal Stethoscope Temp. Probe 一次性听诊体温传感器 .....	132
2.4a Accessories 附件 .....	132
2.5 Temp. Adapter Cable for disposable probe 一次性体温探头转接线 .....	133
<b>3. Smart Temp. Probe 快速体温探头</b> .....	<b>135</b>
3.1 Smart Temp.Probe 快速体温探头 .....	136

# Temperature Probe & Interface Cable

## 体温探头与连接电缆

### Certification 产品认证



### Compliance with 符合标准

ISO 80601-2-56  
BS EN12470-4:2001  
EN 60601-1:1990  
ISO10993-1,5,10:2009

### Compatibility 兼容性

Major brands Temperature monitoring system  
主要品牌体温监护系统



## 1.1 Reusable Temp. Probe 重复用体温探头

### Reference Number Rule Description 产品型号说明

TSM2-15AS	
T	Temperature probe 体温探头
SM	Factory/model (See others left) 厂家/型号见左页
2	Dual channel(Default is single channel) 双通道(默认为单通道)
15	1.5m cable length(Default:=3.0m) 长度为1.5米(默认长度为3.0米)
A	Adult(P= Pediatric) 成人(P=儿童)
S	Skin(G=general purpose) 体表(G=通用)



A. Adult general 成人通用 C. Adult skin 成人体表  
B. Pediatric general 儿童通用 D. Pediatric skin 儿童体表

## ① Reusable Temp. Probe 重复用体温探头

Compatible with major patient  
monitor brands 兼容主要品牌

Air Shielded(AS) Spacelabs(SL)  
Datex-Ohmeda(DX) Siemens(SM)  
Marquette0(MQ) YSI(YSI)  
Philips(HP)











### Features 产品特征




Accurate measurement 精确测量  
Quick response 快速读取数据  
Low noise signal transfer 抗干扰性强  
Fully compatible with OEM technology 完全兼容原厂技术  
12 month warranty 12个月保修期

### Specification 产品规格

Probe type 探头类型	General purpose; Skin 通用、体表
Cable length 电缆长度	3.0m/1.5m;unless requirement applied 3.0米、1.5米(特殊要求除外)
Cable diameter 线材直径	2.60mm(Skin purpose) 4.00mm(General purpose) 2.60毫米(体表探头) 4.00毫米(通用探头)
Resister type 电阻类别	NTC/YSI400/2.252K@25°C NTC/10K@25°C
Cable Material 线材材质	TPU jacket ,Light gray TUP外皮, 灰色
Connector 接头	Golden plated pins 镀金端子

OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	OEM REF# OEM型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
Air Shield						
	TAS-15AS		Adult skin 成人体表	1.5m		C3
	TAS-15PS		Pediatric skin 儿童体表	1.5m		C4
	TAS-15AG		Adult general 成人通用	1.5m		C3
	TAS-15PG		Pediatric general 儿童通用	1.5m		C4
Artema-S&W						
	TSW-AS		Adult skin 成人体表	3.0m		C3
	TSW-PS		Pediatric skin 儿童体表	3.0m		C4
	TSW-AG		Adult general 成人通用	3.0m		C3
	TSW-PG		Pediatric general 儿童通用	3.0m		C4
ATOM						
Infant Warmer: Model V-3200D/V-505 Infant Incubator: Model V-2100G/V-2200 Transport Incubator:V-707V-808	TAT-15PS		Pediatric skin 儿童体表	1.5m		C4
Biolight						
V6 Vital Sign	TBLT-AS		Adult skin 成人体表	3.0m		C3
AnyView	TBLT-PS		Pediatric skin 儿童体表	3.0m		C4
A2/A3/A5/A6/A8	TBLT-AG		Adult general 成人通用	3.0m		C3
	TBLT-PG		Pediatric general 儿童通用	3.0m		C4
Beijing Julongsanyou						
VInfant Radiant Warmer BNT-1000/BNT-2000/BNT-2000B	T20-15PS		Pediatric skin 儿童体表	1.5m		B7
	T20-PS		Pediatric skin 儿童体表	3.0m		B7
Beijing Julongsanyou						

OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	OEM REF# OEM型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
<b>Beijing Julongsanyou</b>						
Infant Incubator B-800/B-1000/B-2000/B-3000	TBNT-15PS		Pediatric skin 儿童体表	1.5m		D1
	TBNT-PS		Pediatric skin 儿童体表	3.0m		D1
<b>David</b>						
HKN-90/HKN-93/HKN-93A	TDV-15PS		Pediatric skin 儿童体表	1.5m		D1
HKN-93B/HKN-9010/HKN-2000	TDV-PS		Pediatric skin 儿童体表	3.0m		D1
<b>GE Marquette</b>						
Eagle, Tram, Solar, Mac-lab Cath lab	TMQ-AS		Adult skin 成人体表	3.0m		C3
	TMQ-PS		Pediatric skin 儿童体表	3.0m		C4
	TMQ-AG		Adult general 成人通用	3.0m		C3
	TMQ-PG		Pediatric general 儿童通用	3.0m		C4
	TMQ2-AS	Dual temp 双探头	Adult skin(2) 成人体表(2)	3.0m		D1
	TMQ2-AS/AG	Dual temp 双探头	Adult skin & general 成人体表 & 通用	3.0m		D1
<b>Philips</b>						
All	THP-AS	21078A	Adult skin 成人体表	3.0m		C3
	THP-PS		Pediatric skin 儿童体表	3.0m		C4
	THP-AG	21075A,21076A	Adult general 成人通用	3.0m		C3
	THP-PG		Pediatric general 儿童通用	3.0m		C4
<b>Mindray</b>						
T5/T8 DPM 6,DPM 7	TMR-AS	0011-30-37393	Adult skin 成人体表	3.0m		C3
	TMR-PS	0011-30-37395	Pediatric skin 儿童体表	3.0m		C4
	TMR-AG	0011-30-37392	Adult general 成人通用	3.0m		C3
	TMR-PG	0011-30-37394	Pediatric general 儿童通用	3.0m		C4
<b>Datascope</b> Passport V, V12, V21	TDT-AS	040-000057-00	Adult skin 成人体表	3.0m		C4
	TDT-PS	040-000058-00	Pediatric skin 儿童体表	3.0m		C5
	TDT-AG	040-000055-00	Adult general 成人通用	3.0m		C4
	TDT-PG	040-000056-00	Pediatric general 儿童通用	3.0m		C5
<b>Siemens/Drager</b>						
	TSM-15AS	4329822(1.5m)	Adult skin 成人体表	1.5m		C4
	TSM-AS	5204669(3.0M)		3.0m		
	TSM-15PS		Pediatric skin 儿童体表	1.5m		C5
	TSM-PS			3.0m		
	TSM-15AG	4329889(1.5m)	Adult general 成人通用	1.5m		C4
	TSM-AG	5204644(3.0m)		3.0m		
	TSM-15PG	4329848(1.5m)	Pediatric general 儿童通用	1.5m		C5
	TSM-PG	5204651(3.0m)		3.0m		
<b>Spacelabs</b>						
	TSL-AS		Adult skin 成人体表	3.0m		C7
	TSL-PS		Pediatric skin 儿童体表	3.0m		C8
	TSL-AG		Adult general 成人通用	3.0m		C7
	TSL-PG		Pediatric general 儿童通用	3.0m		C8
	TSL2-AS		Adult skin(2) 成人体表(2)	3.0m		D3
	TSL2-AS/AG		Adult skin & general 成人体表 & 通用	3.0m		D3
<b>YSI</b>						
YSI 400 Series	T2252-AS	409B	Adult skin 成人体表	3.0m		B6
	T2252-PS	429	Pediatric skin 儿童体表	3.0m		B7
	T2252-AG	401	Adult general 成人通用	3.0m		B6
	T2252-PG	402	Pediatric general 儿童通用	3.0m		B7
YSI 700 Series	T700-AS	709B	Adult skin 成人体表	3.0m		C3
	T700-PS	729	Pediatric skin 儿童体表	3.0m		C4
	T700-AG	701	Adult general 成人通用	3.0m		C3
	T700-PG	702	Pediatric general 儿童通用	3.0m		C4

OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	OEM REF# OEM型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
<b>10K@25°C</b>						
	T10-AS/T10-09AS		Adult skin 成人体表	3.0m/0.9m		B6/B5
	T10-PS		Pediatric skin 儿童体表	3.0m		B7
	T10-AG		Adult general 成人通用			B6
	T10-PG		Pediatric general 儿童通用	3.0m		B7
<b>Ohmeda</b>						
Ohmeda Incubator: Giraffe/Omnibed	TOH-PS		Neonate skin 新生儿体表	3.0m		C4
<b>Air Temperature</b>						
	T100-GP		For air Temp. measurement 空气温度测试	3.0m		C4

Unimed REF# 尤迈型号	Descriptions & Application 描述及应用	PCS /Per box 片/盒	Price code 价格代码
T004-16	Disposable temperature probe cover Used with reusable skin temp. probe to fixed the probe onto the skin and reflect the ambient light to ensure the accurate measurement.  一次性体温探头贴片 该贴片使用在重复性体温探头上面, 将探头端固定在人体表皮从而确保探头读取数据的精确性。	16	A5



The probe cover package not includes the reusable probe. Separated order required.  
体温探头贴片不包含重复用体温探头在内。

# 1.2 Temp. Adapter Cable for reusable probe

## 重复用体温探头连接电缆

### Reference Number Rule Description 产品型号说明

TSM-AD/2	
T	Temperature probe 体温探头
SM	Factory/model (See other on previous page) 厂家/型号见前页
2	Dual channel(Default is single channel) 双通道(默认为单通道)
AD	Adapter, used with reusable YSI probes 转接线, 与YSI重复用探头使用



OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Probe end Connector 探头端插头	Price code 价格代码
<b>Artema-S&amp;W</b>						
	<b>TSW-AD</b>	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
<b>Biolight</b>						
V6 Vital Sign AnyView A2/A3/A5/A6/A8	<b>TBLT-AD</b>	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
<b>GE Marquette</b>						
Eagle, Tram, Solar, Mac-lab Cath lab	<b>TMQ-AD</b>	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
	<b>TMQ400-AD/2</b> 402015-004	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C4
	<b>TMQ700-AD/2</b> 402015-004	Used with reusable YSI700 compatible probes 与YSI700系列重复用体温探头连接使用	0.3m			C4

OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Probe end Connector 探头端插头	Price code 价格代码
<b>Philips Medical</b>						
All	<b>THP-AD</b> 989803163281	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
<b>Mindray</b>						
T5/T8 DPM 6, DPM 7	<b>TMR-AD</b>	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
<b>Datascope</b> Passport V, V12, V21	<b>TDT-AD</b>	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
<b>Siemens/Drager</b>						
	<b>TSM-AD</b> 5198333	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C3
<b>Spacelabs</b>						
	<b>TSL-AD</b>	Used with reusable YSI400 compatible probes 与YSI400系列重复用体温探头连接使用	0.3m			C5
90496, 91496	<b>TSL-AD/2</b> 700-0031-00	Used with reusable YSI400 or YSI700 compatible probes 与YSI400或YSI700系列重复用体温探头连接使用	0.3m			D5
90496/91496 306655-002/306644-001	<b>CO-SL-AD</b> 700-0027-00	CO adapter cable 心排量连接电缆	0.3m			D3


## 2.1 Disposable Skin, Esophageal/Rectal Temp. Probe

一次性体表，体腔体温传感器

Reference Number Rule Description

产品型号说明

TMQ-DAG-20	
T	Temperature probe 体温传感器
MQ	Factory/model (See left) 厂家/机型 (见左页)
D	Disposable Probe 一次性体温探头
A	Adult(P= Pediatric) 成人(P=儿童)
G	General purpose (S=Skin) 通用(S=体表)
20	20 pcs per box (10=10pcs/box) 20个/盒(10=10个/盒)



OEM Technology 原厂技术	REF# 型号	OEM REF# OEM型号	Application 应用	Cable length 线长	Probe adapter 探头转接头端	Adapter cable 连接线	Package 包装	Price 价格代码
400 series	TMQ-D-20	M1024222	Skin 体表	75cm/90cm		See the next page of adapter cable 见下页连接线	20	D5
	TMQ-DAG-20	M1024229	General 成人体腔	75cm/90cm		See the next page of adapter cable 见下页连接线	20	D5
Airshield	TAS30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TAS30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Artema/S&W	TSW30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TSW30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Biolight	TBLT30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TBLT30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Beijing Julongsanyou Infant Radiant Warmer BNT-1000BNT-2000BNT-2000B	T2015-D-10		Skin 体表	150cm		No need 不需要连接线	10	D4
	T2030-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
Beijing Julongsanyou Infant Incubator B-800B-1000B-2000B-3000	TBNT15-D-10		Skin 体表	150cm		No need 不需要连接线	10	D4
	TBNT30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
Drager/Siemens	TDG-D-20		Skin 体表	75cm/90cm		NA	20	D5
	TDG-DAG-20		General 成人体腔	75cm/90cm		NA	20	D5
Mindray T5/T8	TMR30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TMR30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Mindray - Datascope	TDT30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TDT30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Ohmeda,10K	TOH30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TOH30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Philips Medical	THP-D-20	21091A	Skin 体表	75cm/90cm		See the next page of adapter cable 见下页连接线	20	D5
	THP-DAG-20	21090A,M1837A	General 成人体腔	75cm/90cm		See the next page of adapter cable 见下页连接线	20	D5
Philips Medical	TPM30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TPM30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4
Siemens/Drager	TSM16-D-10	MP00995(1.6m) MP00997(1.6m)	Skin 体表	160cm		No need 不需要连接线	10	D4
	TSM30-D-10	MP00996(3.0m) MP00998(3.0m)	Skin 体表	300cm		No need 不需要连接线	10	D4
	TSM16-DAG-10	MP00991(1.6m)	Adult General	160cm		No need 不需要连接线	10	D4
	TSM30-DAG-10	MP00992(3.0m)	Adult General	300cm		No need 不需要连接线	10	D4
	TSM16-DPG-10	MP00993(1.6m)	Pediatric General	160cm		No need 不需要连接线	10	D4
TSM30-DPG-10	MP00994(3.0m)	Pediatric General	300cm		No need 不需要连接线	10	D4	
Spacelabs	TSL30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D5.2.5
	TSL30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D5.2.5
YSI400	TYSI30-D-10		Skin 体表	300cm		No need 不需要连接线	10	D4
	TYSI30-DAG-10		General 成人体腔	300cm		No need 不需要连接线	10	D4

## ② Disposable Temp. Probe

### 一次性体温传感器

Compatible with major patient monitor brands 兼容主要品牌

Air Shield(AS) Mindray-Datascope(DT)  
Drager(DG) Philips(HP)  
Marquette(MQ) YSI400(YSI)  
Mindray(MR) Ohmeda(OH)

Features 产品特征

Single patient use 单个病人使用  
Use in Esophagus/rectum/nasopharynx/axilla, skin, ear, bladder  
可用在食道/直肠/鼻咽/腋窝, 皮肤, 耳朵, 膀胱处  
Accurate measurement 精准测量  
Quick response 快速读取数据  
Low noise signal transfer 抗干扰性强  
Fully compatible with OEM technology 完全兼容原厂技术

Specification  
产品规格

Probe type 探头类型	Esophagus/rectum/nasopharynx/axilla, skin, ear, bladder 食道/直肠/鼻咽/腋窝, 皮肤, 耳朵, 膀胱
Cable length 电缆长度	0.75m/0.9m/3.0m; unless requirement applied 0.75m/0.9m/3.0m; 特殊要求除外
Resistor type 电阻类别	NTC/YSI400/2.252K@25°C
Cable Material 线材材质	PVC jacket ,white PVC外皮, 白色
Connector 接头	Golden plated pins 镀金端子



## 2.2 Disposable Tympanic Temp. Probe

### 一次性耳鼓体温传感器

#### Specification

#### 产品描述

- Monitoring the tympanic membrane temperature in the ear canal, well correlates the core body temperature.
- Dual foam protection with soft memory cell foam to minimize the risk of injury to the tympanic membrane, and seal out the ambient air.
- PVC push tube makes insertion and pulling easily, and protects the lead wire from kink and break.
- Thinner lead wire is more pliable and provides accurate temperature readings.
- 用于监测耳膜温度，被认为最接近核心体温
- 双层记忆泡棉降低对耳膜的损伤，并能有效阻隔外界空气
- 使用PVC 插管能轻易插入而不使线纠结，轻松拉出而不易拉断
- 电缆线小巧柔韧，并能提供精确的体温读数



OEM Technology 原厂技术	REF# 型号	OEM REF# OEM型号	Application 应用	Probe end Connector 探头端插头	Adapter cable 连接线	Package 包装	Price 价格代码
400 series	TMQ-PE		Pediatric Tympanic Temperature probe 400系列儿童耳鼓体温传感器		See the next page of adapter cable 见下页连接线	1	B2
	TMQ-AE		Adult Tympanic Temperature probe 400系列成人耳鼓体温传感器			1	B2

## 2.3 Foley Catheter Temp. Probe

### 一次性导尿管体温传感器

#### Specification

#### 产品描述

- Monitoring the bladder temperature during urinary drainage
- Commonly used to monitor core body temperature during open-heart surgery or critically ill patients
- Provides accurate temperature readings
- Silicone construction and latex-free prevent irritation and sensitive reaction.
- Water-proof design avoids liquid seepage
- 导尿的同时，进行体温监测
- 广泛用于监测心脏切开手术或一些危重病人的核心体温
- 产品读数精确
- 硅胶和无乳胶的使用避免刺激与过敏反应
- 防水设计防止液体渗入



OEM Technology 原厂技术	REF# 型号	OEM REF# OEM型号	Application 应用	Probe end Connector 探头端插头	Adapter cable 连接线	Package 包装	Price 价格代码
400 series	TMQ-CA-14FR		Adult Foley catheters Temperature Probe,14Fr 400系列成人导尿管体温传感器，14Fr		See the next page of adapter cable 见下页连接线	1	B6
	TMQ-CA-16FR		Adult Foley catheters Temperature Probe,16Fr 400系列成人导尿管体温传感器，16Fr			1	B6
	TMQ-CA-18FR		Adult Foley catheters Temperature Probe,18Fr 400系列成人导尿管体温传感器，18Fr			1	B6

## 2.4 Disposable Esophageal Stethoscope Temp. Probe

### 一次性听诊体温传感器

#### Specification

#### 产品描述

- Monitoring the temperature in the esophagus
- Anesthesiologist can hear the heart and lung sounds through the ear piece while monitoring the temperature.
- Acoustic balloon design to enhance sound quality
- Smooth, heat-bonded cuff to tube transition minimizes potential trauma and irritation
- Can use with esophageal stethoscope(Optional)
- 用于食道温度监测
- 监测体温的同时，可传送心脏和肺呼吸的声音到麻醉师的耳朵
- 原声的球囊设计可提高声音质量
- 原声球囊使用圆滑的热粘合技术减少对食道的创伤与刺激
- 可选择性的与食道听诊器一起使用



OEM Technology 原厂技术	REF# 型号	OEM REF# OEM型号	Application 应用	Probe end Connector 探头端插头	Adapter cable 连接线	Package 包装	Price 价格代码
400 series	TMQ-OP-9FR		Pediatric Esophageal Stethoscope Temperature probe, 9Fr, with ear piece 400系列儿童食道体温传感器，9Fr，含食道听诊器		See the next page of adapter cable 见下页连接线	1	B3
	TMQ-OP-12FR		Adult Esophageal Stethoscope Temperature probe, 12Fr, with ear piece 400系列成人食道体温传感器，12Fr，含食道听诊器			1	B3
	TMQ-OP-18FR		Adult Esophageal Stethoscope Temperature probe, 18Fr, with ear piece 400系列成人食道体温传感器，18Fr，含食道听诊器			1	B3
	TMQ-OP-24FR		Adult Esophageal Stethoscope Temperature probe, 24Fr, with ear piece 400系列成人食道体温传感器，24Fr，含食道听诊器			1	B3

## 2.4a Accessories

### 附件

#### Stethoscope Ear Piece 食道听诊器

#### Specification

#### 产品描述

- Use with disposable Esophageal Stethoscope Temp Probe
- Connect ear piece extension tubing to the stethoscope Luer connector for listening to heart and lung sounds.
- Acoustical design to get the good sound quality
- 与一次性听诊体温探头配套使用
- 将食道听诊器连接管连续接到听诊体温传感器的鲁尔接头，听取病人的心脏及肺呼吸的声音
- 原声设计提高声音质量



OEM Technology 原厂技术	REF# 型号	OEM REF# OEM型号	Application 应用	Probe end Connector 探头端插头	Adapter cable 连接线	Package 包装	Price 价格代码
400 series	TSTP-EP		Stethoscope Ear Piece 食道听诊器			1	A3



## 2.5 Temp. Adapter Cable for disposable probe

### 一次性体温探头转接线


#### Reference Number Rule Description

#### 产品型号说明

THP-30-AD	
T	Temperature Adapter Cable 体温转接线
HP	Factory/model (See others on the previous page) 厂家/机型 (见前页)
30	3.0m cable length 3.0米线长
AD	used with YSI 400 series disposable probe(PH: Philips disposable probe) 需与YSI 400系列一次性体温探头一起使用 (PH=飞利浦一次性体温探头)

OEM & Models 原厂品牌 / 机型	Probe end Conn 探头端接头	UNIMED REF# 尤迈型号	Probe end Conn 探头端接头	UNIMED REF# 尤迈型号	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
	For 400 Disposable 适用400系列一次性体温探头		For Philips Disposable 适用飞利浦一次性体温探头				
<b>Artema-S&amp;W</b>							
		TSW-30-AD		TSW-30-PH	3.0m		C3
<b>Biolight</b>							
		TBLT-30-AD		TBLT-30-PH	3.0m		C3
<b>GE Marquette</b>							
Eagle, Tram, Solar, Mac-lab Cath lab		TMQ-30-AD 2021701-001		TMQ-30-PH	3.0m		C3
<b>Philips Medical</b>							
All		THP-30-AD		THP-30-PH (21082A),3.0m (21082B),1.5m	3.0m		C3

OEM & Models 原厂品牌 / 机型	Probe end Conn 探头端接头	UNIMED REF# 尤迈型号	Probe end Conn 探头端接头	UNIMED REF# 尤迈型号	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
	For 400 Disposable 适用400系列一次性体温探头		For Philips Disposable 适用飞利浦一次性体温探头				
<b>Mindray</b>							
T5/T8 DPM 6, DPM 7		TMR-30-AD		TMR-30-PH	3.0m		C3
Passport V, V12, V21		TDT-30-AD 040-000091-00		TDT-30-PH	3.0m		C3
<b>Siemens</b>							
		TSM-30-AD		TSM-30-PH	3.0m		C3
<b>Spacelabs</b>							
		TSL-30-AD		TSL-30-PH	3.0m		C5
<b>YSI 400</b>							
		TYSI-30-AD (165640),1.3m (165641),2.8m (9445-801),3.6m		TYSI-30-PH (989803162591)3.0m (989803162601)1.5m	3.0m		B8





### 3.1 Smart Temp. Probe 快速体温探头



#### Reference Number Rule Description 产品型号说明

TOP-AR-001	
T	Temperature 体温
O	Oral口腔 (Rectal直肠)
P	Probe 探头
AR	Factory Code 厂家代码 (AR=Alaris,GE=GE Healthcare,MR=Mindray)
001	Factory /model 厂家/型号

## ③ Smart Temp. Probe 快速体温探头

Compatible with major patient monitor brands 兼容主要品牌

Major brands temperature monitoring system  
主要品牌体温监护系统

#### Features 产品特征

- Accurate measurement 精确测量
- 15s Quick response 15秒快速读取数据
- Low noise signal transfer 抗干扰性强
- Fully compatible with OEM technology 完全兼容原厂技术
- 12 month warranty 12个月保修期

#### Specification 产品规格

Probe type 探头类型	Oral;Rectal 口腔; 直肠
Cable length 电缆长度	4m unless requirement applied 4米 特殊要求除外
Resister type 电阻类别	NTC NTC
Cable Material 线材材质	TPU jacket,Light gray TPU外皮, 灰色
Package 包装	1pcs per bag 每条一包

OEM & Models 原厂品牌 / 机型	Unimed REF# 尤迈型号	OEM REF# OEM型号	Application 应用	Cable length 线长(米)	Monitor end plug 机器端插头	Price code 价格代码
<b>Alaris</b>						
VITAL CHECK 4400/4410/4500	TOP-AR-001	3880L	Oral 成人口腔	4m		D5
	TRP-AR-001	3882L	Rectal成人直肠	4m		D5
Turbo Temp 2180C/2180CX01EE	TOP-AR-002	2885C	Oral 成人口腔	4m		D5
	TRP-AR-002	2886C	Rectal成人直肠	4m		D5
Turbo Temp 2185B/2185BX01EE	TOP-AR-003	2887B	Oral 成人口腔	4m		D5
	TRP-AR-003	2888B	Rectal成人直肠	4m		D5
<b>GE Healthcare</b>						
DINAMAP Pro200/400 V1	TOP-AR-001	088012	Oral 成人口腔	4m		D5
	TRP-AR-001	088016	Rectal成人直肠	4m		D5
DINAMAP Pro200/400 V2 ProCare200/400	TOP-AR-002	2008774-001	Oral 成人口腔	4m		D5
	TRP-AR-002	2008775-001	Rectal成人直肠	4m		D5
DINAMAP 9300 Series DINAMAP XL	TOP-GE-001	008975,008978	Oral 成人口腔	4m		D5
	TRP-GE-001	008976	Rectal成人直肠	4m		D5
<b>Mindray/迈瑞</b>						
Accutorr Plus	TOP-MR-001	0206-00-0725-01	Oral 成人口腔	4m		D5
	TRP-MR-001	0206-00-0725-02	Rectal成人直肠	4m		D5



# ETCO2 Sensor

## 呼气末二氧化碳传感器

### Certification 产品认证



### Compliance with 符合标准

- ENISO 14971:2009
- EN 980:2010
- EN 1041:2010
- EN 60601-1-2:2007
- ISO 21647:2004
- ENISO 14155-1:2009

Unimed offers both Mainstream ETCO2 Sensor and External Sidestream ETCO2 Sensor which are the supplement of the host device without the internal ETCO2 module, they are compatible with most of the major machines in market. Our sensors with high performance , accurate measurement and long service life, are also easy to operate and plug-in-measure and with good design.

尤迈系列的呼气末二氧化碳传感器（ETCO2 Sensor）分为主流和外加式旁流两大类，是没有 配备内置呼气末二氧化碳模块的机型的一个补充。能兼容市场上大部的主流机型。产品性能稳定，监测数据精确，使用寿命长，同时也具备了即插即测，操作简单，外观精致等优点。

Description 描述	Unimed REF# 尤迈型号	Sensor Picture 产品图片	Compatibility 兼容性	Connector 接头	Price code 价格代码
External Sidestream CO2 Sensor (with one set of standard accessories) 外加式旁流呼气末二氧化碳传感器 (带一套标准配件)	UECO2-001		Phase in (ISA AX+)		F7
	UECO2-002		Philips (M3014A)		F7
	UECO2-003		Respironics (C5)		F7
	UECO2-004		Zoll (E series,M series)		F7
Mainstream ETCO2 sensor (with one set of standard accessories) 主流呼气末二氧化碳传感器 (带一套标准配件)	UMCO2-001		Phase in (IRMA CO2)		F7
	UMCO2-002		Philips(M3014A)		F7
	UMCO2-003		Respironics (C5)		F7
	UMCO2-004		Zoll (E series,M series)		F7
	UMCO2-005		Mindray (T5,T8)		F7

Pls be Reminded : Unimed provides the accessories as below.  
请注意：尤迈有以下CO2附件供选择.





# Multi-Parameter Patient Cable

## 多参数病人电缆

### Certification 产品认证



CE Mark  
CE 标志



510(K) clearance  
No.: K082959



SFDA医疗器械注册证:  
No.:2210215

### Compliance with 符合标准

ANSI/AAMI EC53:1998  
EN 60601-1:1990  
ISO10993-1,5,10:2009

### Compatibility 兼容性

Major brands ECG or multi-parameter monitoring system  
主要品牌的多参数监护仪系列

REF# 型号	Used with the OEM and models 适用的OEM型号	Descriptions 描述	Cable plug 机器端接头	Package 包装	Price code 价格代码
<b>SM-STE5</b>					
(3368391), 2.5m	Siemens/Drager Infinity 6000 / 7000 / 9000 Drager Infinity	Multimed 5 Pod patient cable, 2.5m or 1.5m 5-lead ECG, VS/Siemens Spo2/Temp/Resp. AHA & IEC		1pcs	D8
<b>SM15-STE5</b> (5950192), 1.5m	Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	多参数电缆, 2.5M或1.5M 心电5导, VS/Siemens 血氧/体温/呼吸 美标及欧标		1pcs	D8
<b>SM-STE6</b>					
(5191221)	Siemens/Drager Infinity 6000 7000 9000 Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL	Multimed 6 Pod patient cable, 2.5m 6-lead ECG, VS/Siemens Spo2/Temp/Resp. AHA & IEC 多参数电缆, 2.5M 心电6导, VS/Siemens 血氧/体温/呼吸 美标及欧标		1pcs	D8.5.8
<b>SM-STE12</b>					
(5589663)	Drager Infinity Delta/Delta XL Kappa/Kappa XLT	Multimed 12 Pod patient cable, 3.0m 12-lead ECG, VS/Siemens Spo2/AUX. AHA & IEC 多参数电缆, 3.0m 心电12导, VS/Siemens 血氧/AUX. 美标及欧标		1pcs	NA
<b>SM-MX5</b>					
(5590539)	Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	NeoMed Pod Patient cable, 2.5m ECG/Resp./Spo2 Dual Temp/FiO2 AHA & IEC 新生儿多参数电缆, 2.5m 心电/呼吸/血氧 双通道体温/氧浓度 美标及欧标		1pcs	D8.5.8
<b>DG-STE6</b>					
(MS20093)	Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	Multimed Plus patient cable, 2.5m Spo2/Temp 6-lead ECG/Resp. AHA & IEC 多参数主电缆, 2.5m 血氧/体温 心电/6导呼吸 美标及欧标		1pcs	E1
<b>DG-STE6-OR</b>					
(MS20094)	Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/Delta XL Kappa/Kappa XLT Vista/Vista XL	Multimed Plus OR patient cable, 2.5m Spo2/Temp 6-lead ECG/Resp. AHA & IEC 多参数主电缆(手术室用), 2.5m 血氧/体温 心电6导/呼吸 美标及欧标		1pcs	NA
<b>SM-STE6</b>					
Y05-T091	Unimed compatible multi-parameter patient cable like SM-STE5. 尤迈兼容多参数病人电缆如SM-STE5	The holder of the multi-parameter cable yoke to fix onto the pole or wall on bedside. 多功能参数病人电缆分线器 挂扣, 将电缆固定在杆或床头		1pcs	C2
Y05-T092	Siemens/Drager original multi-parameter patient cable. Siemens/Drager 原装多参数病人电缆	The purpose is the Same as above 功能同上		1pcs	C6



# Veterinary VitalSign Monitoring Cables & Accessories

## 动物监护电缆及配件



Certification  
产品认证



CE Mark  
CE 标志

Compliance with  
符合标准

ISO 80601-2-61  
ANSI/AAMI EC53:1998  
EN 60601-1:1990  
ISO10993-1,5,10:2009



## 1 SpO2 Sensor & Accessories 血氧探头及配件

Sensor Type 探头类型	Sensor REF# 探头型号	Holder REF# 支架型号	Unit Price 单价	
Rectal/Reflect 直肠/反射式		REF# Y04-E050	\$9.90	
Fig.A	U800V series 1.1m/3.0m cable length  Compatible with all major brands like 兼容所有主要品牌, 如: Nellcor、BCI、Nonin...	Fig.B		
Lingual clip 舌夹式		Holder 1 REF# 支架1型号	Holder 2 REF# 支架2型号	Holder 3 REF# 支架3型号
Fig.C	U800 Multi-Site series 1.1m/3.0m cable length  Compatible with all major brands like 兼容所有主要品牌, 如: Nellcor、BCI、Nonin...	REF# UEC-01 Fig.D	REF# UVC-02 Fig.E	REF# UVC-01 Fig.F
		Size: Small 尺寸: 小号	Size: Medium 尺寸: 中号 L: 58mm 长: 58毫米 W: 16mm 宽: 16毫米	Size: large 尺寸: 大号 L: 68mm 长: 68毫米 W: 32mm 宽: 32毫米
		\$2.90	\$4.90	\$5.90

Contact Unimed sales representative for the REF# per the oximeter to be used.  
根据需求的血氧仪选择合适型号的探头, 产品详情请联系业务代表。

## 2 ECG adapter 心电转换头

Universal Vet ECG adapter  
多功能兽用心电转换头

Used with standard ECG leads with following electrodes  
Snap end; Grabber; 3mm needle or 4mm banana end electrode  
配合带以下电极的标准心电导联线使用  
扣式, 钳式, 3毫米针式或4mm香蕉头电极

Ref to fig. 图片编号	REF# 型号	Package 包装	Price code 价格代码
Fig.G	UVC-05-10 UVC-05	10pcs/bag 1 pcs	D1 A7